SL* -GLUCOSE (Enzymatic-colorimetric, End Point)

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
- Sacks, D.B., Carbohydrates. Tietz Fundamentals of Clinical Chemistry, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 427.
- Dods, R.F., Diabetes Mellitus. Clinical Chemistry: Theory, Analysis, Correlation, 4th Ed. Kaplan, L.A, Pasce, A.J., Kazmierczak, S.D., (Mosby, Inc, Eds St Louis USA), (2003), 580.
- Tietz, N.W., Clinical guide to laboratory tests, 4rd Ed. (W.B. Saunders eds. Philadelphia USA), (2006), 444
- Burrin JM., Price CP., Measurement of blood glucose. Ann.Clin.Biochem. (1985), 22, 131.
- Passey RB, Gillum RL, Fuller JB, Urry FM, Giles ML. Evaluation and comparison of 10 glucose methods and the reference method recommended in the proposed product class standard Clin.Chem(1977), 23, 131.
- Kaplan, A.L., Carbohydrate and metabolites. Clinical Chemistry: Theor, Analysis, Correlation, 2ndEd. Kaplan, L.A, Pasce, A.J., (Mosby, Inc. eds St Louis USA), (1989), 850.
- Vassault A., et al., Ann Bio. Clin., (1986), 44, 686.
- Vassault A., et al., Ann Bio. Clin., (1999), 57, 685.
- Berth, M. & Delanghe, J., Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Act Clin Belg., (2004), 59, 263.
- Young, D.S., Effects of preanalytical variable on clinical laboratory tests, 2ndEd. AACC Press, (1997).
- Young, D.S., Effects of drugs on clinical laboratory tests, 4thEd. AACC Press, (1995).
- <https://www.mayocliniclabs.com>
- Peter E., Gus Koerbin., Methods in Clinical chemistry, Kaplan and Pesce's : Clinical Chemistry:Theory, Analysis, Correlation, (2009), 651.

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	Temperature limitation		Catalogue number
	Manufacture address		Expiration date
	Batch code		Date of manufacture
	In vitro diagnostic medical device		Reagent 1
	Consult instruction for use		