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Sample	Original Concentration (IU/L)	50 IU/L added		100 IU/L added	
		Recovery%	Bias%	Recovery%	Bias%
1	36.2	108.00%	4.64%	111.00%	8.08%
2	312	96.00%	-0.55%	92.00%	-1.94%
3	715	107.00%	0.46%	113%	1.6%

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
1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Kasper D.L. et al. (2015). Harrison's principles of internal medicine.19th ed. Mc Graw Hill Education
5. Marcilac I. Troalen F. et al. (1992) Free human chorionic gonadotropin beta subunit in gonadal and nongonadal neoplasms. Canc Res. 1992; 52:3901-3907.
6. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests.6th ed. Elsevier Inc.
7. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. Elsevier Inc

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INTENDED USE:

For the quantitative determination of human chorionic gonadotropin (hCG) concentration in human serum.

CLINICAL SIGNIFICANCE:

Similarly, to LH, FSH and TSH, human chorionic gonadotropin (hCG) is a member of the glycoprotein family and consists of 2 subunits (α - and β -chains) which are associated as the intact hormone. The α -chains in all four of these glycoprotein hormones are virtually identical, whereas the β -chains have greatly differing structures and are responsible for the respective specific hormonal functions.

Human chorionic gonadotropin consists of a number of isohormones with differing molecular size. The biological action of hCG serves to maintain the corpus luteum during pregnancy, and hence, progesterone production, during the first trimester. The serum of pregnant women contains mainly intact hCG.

Measurement of the hCG concentration permits the diagnosis of pregnancy just one week after conception. The determination of hCG in the 1st trimester of pregnancy is of particular importance.

Elevated values here serve as an indication of hydatid form mole or multiple pregnancy. Depressed values indicate threatening or missed abortion, ectopic pregnancy or intra-uterine death. Elevated values in the absence of a pregnancy are indicative of a tumor. The specific monoclonal antibodies used recognize the hCG-hormone. The PSI hCG titration kit should therefore be used in particular in the diagnosis and monitoring of pregnancy.

PRINCIPLE:

The PSI hCG EIA Kit is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. In the 1st incubation phase hCG present in serum samples or calibrators/controls reacts with biotinylated monoclonal antibody against alpha subunit. The strips are then washed and 2nd HRP-conjugated monoclonal antibody against beta-subunit are added to wells. hCG present in calibrators/controls or samples reacts with both antibodies to form a sandwich complex and also adsorbed to the streptavidin coated microstrips via biotinylated anti-alpha subunit mab. The strips are then washed and, buffered substrate/chromogen reagent (hydrogen peroxide and tetramethylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue color will develop and the intensity of the color is directly proportional to the amount of hCG present in the samples.

hCG-titer ELISA KIT		IVD-REF: PS – HCG.T
		Ver. No: 02

REAGENTS PROVIDED:

Reagent	Quantity/48 Test KIT	Quantity/96 Test KIT	Quantity/192 Test KIT	Reconstitution
Microplate (96 breakable wells) coated with Streptavidin	48 wells	96 wells	192 wells	Ready to Use
Standard 0 (Buffer without hCG compatible with human serum matrix containing preservative)	1 x5.0 mL	1 x5.0 mL	1 x5.0 mL	Ready to Use
Calibrators 2-6 (20, 100, 250, 500, 1000 IU/L) in buffer, containing preservative with traceability to reference material "3 rd ISO WHO 75/539"	5 Vial/0.5 ml	5 Vial/0.5 ml	5 Vial/1.0 ml	Ready to Use
Controls N = 2 in buffer compatible with human serum matrix with preservatives	2 Vial/0.5 ml	2 Vial/0.5 ml	2 Vial/1.0 ml	Ready to Use
Assay Buffer	1 Vial/6.0 ml	1 Vial/12 ml	2 Vial/12 ml	Ready to use
Enzyme Conjugate (red color)	1 Vial/6.0 ml	1 Vial/12 ml	2 Vial/12 ml	Ready to use
Concentrated Wash solution	1 Vial/30ml	1 Vial/30 ml	1 Vial/50 ml	Dilute 20X with distilled water
Substrate-Chromogenic solution (Tetra-methylbenzidine and hydrogen peroxide)	1 Vial/6.0ml	1 Vial/12ml	2 Vial/12 ml	Ready to Use
Stop solution (H2So4 0.12M) Additional information about other additive is considered proprietary.	1 Vial/6.0 ml	1 Vial/6.0 ml	1 Vial/12.0 ml	Ready to Use

Materials required but not supplied with kit

The following materials are required but not provided in the kit:

- Distilled water with conductivity < 1 μ S/cm
- Precision micropipettes for delivery of,25, 50, 100 microliters. An 8-channel pipette or respenser pipette with disposable plastic tips for delivery of 50-100 μ L is useful but not essential. Pipettes for dispensing milliliter volumes.
- Vortex mixer
- Automatic washer or any other device which can deliver 350 microliter wash solution for microplates
- Microplate reader with a wavelength of 450 nm and 620-650 nm, and an absorbance range of 0-3.0

STORAGE AND EXPIRATION DATING OF REAGENTS

- If kept at 2-8 °C, before opening, all kit components are stable until the expiry date as indicated on the label on the outside of the kit box.
- Calibrator concentrations are displayed on vial labels and may vary between lots.
- Do not freeze kits.
- Freshly prepared Working Wash solution should be used on the preparation day.
- Each lot of reagents and calibrators has been standardized to produce the correct reaction.

Do not interchange the reagents or calibrators between lots.

hCG-titer ELISA KIT		IVD-REF: PS – HCG.T
		Ver. No: 02

- Alterations in physical appearance of kit reagents may indicate instability or deterioration. Chromogenic solution TMB should be colorless. A blue color indicates that the reagent has been contaminated and should be discarded.

SPECIMEN COLLECTION AND PREPARATION

- Human serum and EDTA plasma are samples of choice for PSI hCG kit.
- Do not use grossly haemolysed or turbid samples.
- If the measurement cannot be done within 24 hours after sample collection, samples should be kept at 2-8 °C up to 48 hours. For longer storage serum or plasma must be stored at -20°C.
- Avoid repeated freeze and thaw cycles. Frozen specimens must be gently but thoroughly mixed after thawing.

PRECAUTIONS AND WARNINGS

Safety

- For Professional in vitro diagnostic Use Only.
- Before starting the assay, read the instructions completely and carefully. Use the valid version of the instruction for use provided with the kit. Be sure that everything is understood.
- Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., USA) publication No (CDC) 88—8_95 on laboratory safety procedures or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Reagents contain sodium azide (NaN₃) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. Follow local guidelines for disposal of all waste material.
- Avoid any skin contact with all reagents, stop solution contains strong acid. In case of contact, wash thoroughly with soap and water.

Caution

Material used in the preparation of human source reagents has been tested and found to be non-reactive for HIV 1 & 2 antibodies, HCV antibody and hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

Preparation of reagents:

Wash Buffer:

According to volume required, pour the desired volume of wash concentrate into a clean container and dilute 20 fold by adding adequate volume of distilled or deionized water to give a buffered Wash Solution. E.g. 10 mL wash concentrate to 190 mL distilled water. Wash solution should be prepared daily and discard unused working wash solution at the end of the day.

HANDELING NOTES

- Do not use the kit or components beyond expiry date.
- Do not mix up vials caps.
- Bring all the reagents to room temperature prior to use.
- Thoroughly mix all reagents, calibrators/controls and samples by gentle agitation or swirling.
- Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
- Use a clean container and high-quality distilled or deionized water to prepare the working

hCG-titer ELISA KIT		IVD-REF: PS – HCG.T
		Ver. No: 02

wash solution.

7. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
8. For the dispensing of Chromogenic solution and Stop solution avoid pipettes with metal parts.
9. High precision pipettes or automated pipetting equipment will improve the precision.
10. Adhere the incubation times carefully. Tolerance limit of incubation time is $\pm 5\%$
11. To avoid drift, the time between pipetting of the first calibrator and the last sample must not be more than 10 minutes (Time delay). It is strongly recommended not to run more than 5 strips when pipetting manually.
12. Perform all steps consecutively without any delay. Any delay in each step has deleterious effects on final results.
13. Draw a calibration curve for each run, do not use stored curve from previous runs.
14. Incubation with Chromogenic solution, must be done in dark.
15. Opened reagents are stable maximally for 8 weeks if they are not contaminated, and stored in resealed original containers and handled as prescribed. Return to 2-8°C immediately after use

Measurement Procedure

1. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminum pouch containing a desiccant and reseal carefully and store at 2-8°C.
2. Pipette 25 μL of each Calibrator (calibrators 0, 20, 100, 250, 500, 1000 IU/L), control and patient samples in duplicate into the appropriate wells.
3. Pipette 100 μL of Assay Buffer into each well. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid.
4. Gently mix for 15 seconds
5. Cover the plate with an adhesive seal and incubate for 20-minutes at room temperature(20-27°C).
6. Wash the plate according to procedure below.
 - a. Carefully dispense 350 μL of working wash solution into each well.
 - b. Aspirate the content of each well.
 - c. Repeat 4 additional wash for a total of 5 wash cycles.
 - d. After the end of last wash cycle invert the plate and tap it carefully against an absorbent paper.

Note: Automatic strip washer: Follow the manufacturer's instructions for maintenance and wash the required number of wash cycles.

7. Pipette 100 μL of Enzyme conjugate into each well.
8. Cover the plate with an adhesive seal and incubate for 20-minutes at room temperature(20-27°C).
9. Wash the strips 5 times as previously described in procedural note 6.
10. Add 100 μL of Chromogenic solution TMB to each well using the same pipetting technique as in item 3. The Chromogenic solution TMB should be added to the wells as quickly as possible and the time between addition to the first and last well should not exceed 2 min.
11. Incubate for 10 min at room temperature preferentially in the dark.
12. Add 50 μL of Stop Solution, mix and read the absorbance at 450 nm (with a reference filter 620, 630, or 650) in a microplate spectrophotometer within 10 min after addition of Stop Solution.

Important Note: always add stop solution in the same order as chromogenic solution TMB.

hCG-titer ELISA KIT		IVD-REF: PS – HCG.T
		Ver. No: 02


CALCULATION OF RESULTS

1. Construct a calibration curve by plotting the average O.D. (450 nm) for each calibrator on the vertical (Y) axis against the concentration on the horizontal (X) axis, and draw a best fit curve through the points on the graph.
2. To determine concentration of hCG in each sample, first locate the O.D. value on the Y-axis and extend a horizontal line to the calibration curve. At the point of intersection, draw a vertical line to the X-axis and read the corresponding concentration.
3. If a microplate spectrophotometer with built-in data calculation program is used, refer to the manual for the spectrophotometer and create a program using the concentration stated on the label of each of the hCG calibrators. For automatic calculation of hCG results it is recommended to use point-to-point method.

Typical Data

The following data are for illustration only and should never be used instead of the real-time calibration curve.

Row	OD	U/L
1	0.025	0.0
2	0.110	20
3	0.380	100
4	0.800	250
5	1.450	500
6	2.600	1000

Quality control

hCG low and high control serum should be used for validation of the assay series. Ranges of expected results are indicated on the vial labels. If assigned range are not obtained, a complete check of reagents and reader performance should be made and the assay shall be repeated. Each laboratory can also prepare its own serum pools at different levels, which can be used as internal controls in order to ensure the precision of the assay.

Expected values

It is strongly recommended that each laboratory should determine its own normal and abnormal values. In a study conducted with apparently normal healthy adults, using the PSI hCG ELISA the following values are observed:

Pre-menopausal Non-pregnant women	N=158	<10 IU/L
Suspicious	-	10-25 IU/L
Pregnancy	N=120	>25 IU/L

hCG-titer ELISA KIT		IVD-REF: PS – HCG.T
		Ver. No: 02

Pregnant women in different pregnancy weeks

Weeks of pregnancy	2.5 th percentile(IU/L)	97.5 th percentile(IU/L)
4 weeks	5.0	100
5 weeks	200	3000
6 weeks	10000	80000
7-14 weeks	90000	500000
15-26 weeks	5000	80000
27-40 weeks	3000	15000

PERFORMANCE CRITERIA

1. Lower limit of measurement:

The Limit of Blank and Limit of Detection were determined according to CLSI (Clinical and Laboratory Standards Institute) EP17-A2 guideline.

In brief Limit of Blank was determined as the 95th percentile value from n = 60 measurements on standard 0 over several runs. The Limit of blank equals to the concentration that 95% of repeated measurements on zero standard give value below it. In the next step we determined standard deviation of repeated measurements on 3 samples with hCG contents less than 2.0 IU/L determined from a separate method study (ECLIA Method). By adding the LOB calculated at 1st step to 1.645 times the greatest standards deviation calculated in the 2nd step ($LOD = LOB + 1.645SD_L$) LOD calculated as 1.0 IU/L.

2. Precision:

Precision was determined using PSI hCG lit reagents and three pooled human sera according to EP05-A3 Guideline.

Three runs per day (morning, midday and afternoon) by 3 user in 10 working days ($10 \times 3 \times 3$). Results analyzed by fully nested ANOVA. Following results were obtained:

Sample Description	Mean (IU/L)	Repeatability		Within-Laboratory Precision	
		SD	%CV	SD	%CV
Patient Pool	78.2	5.9	7.54	4.2	5.37
Patient Pool	353.8	20.8	5.88	8.1	2.29
Patient Pool	850.6	39.1	4.6	38.3	4.5

3. Specificity:

Specificity evaluation was performed according to guideline CLSI EP07-A2 to assess the cross-reactivity of the assay with glycoproteins hormones LH, FSH and TSH. Samples containing cross-reactants was prepared by spiking natural human sera with material containing high levels of cross-reactants. Modification of serum matrix was no more than 10%. The % cross-reactivity was calculated for each sample using the equation below and normalized to the hCG content:

$$\%cross - reactivity = \frac{\text{mean conc. of spiked sample} - \text{mean conc. of unspiked sample}}{\text{spiked concentration}} \times 100\%$$

The results are summarized in the following table:

hCG-titer ELISA KIT		IVD-REF: PS – HCG.T
		Ver. No: 02

Cross-reactants	concentration	% Cross-reaction
FSH	500 m IU/ml	0.05
LH	500 m IU/ml	0.37
TSH	500 µ IU/ml	0.18

Interference evaluation from common sample abnormalities, hemolysis, Lipemia and icterus were also performed according to aforementioned document and %interference calculated using following formula:

$$\%interference = \frac{\text{measured value} - \text{true value}}{\text{true value}} \times 100\%$$

No interference were seen up to following listed concentration:

Interferent	source	Final concentration of interferent in test serum
Hemoglobin	Fresh haemolysed blood	50 mg/mL
Conjugated Bilirubin	Commercially crystalline Unconjugated Bilirubin	20 mg/dL
Triglycerides	Commercially Intralipid 10%	1000 mg/dL

4. Linearity

Evaluation of Linearity was performed according to CLSI (Clinical and Laboratory Standards Institute) EP 06-A guideline

In brief a sample with concentration of 930 IU/L hCG as a high sample was diluted with another sample with 25 IU/L in different proportions and their hCG, measured in duplicate on each dilution. Measured values compared with expected values which estimated by dilution formula, then % bias and %recovery calculated for each concentration and following results are obtained:

No	Ratio	Expected (IU/L)	Rep 1 (IU/L)	Rep 2 (IU/L)	recovery%	%Bias
1	1	930	930	930	NA	NA
2	0.9	839.5	845	839	100.30%	0.30%
3	0.8	749	751	727	98.66%	-1.34%
4	0.7	658.5	631	667	98.56%	-1.44%
5	0.6	568	515	543	93.13%	-6.87%
6	0.5	477.5	452	457	95.18%	-4.82%
7	0.4	387	375	365	95.61%	-4.39%
8	0.2	206	197	186	92.96%	-7.04%
9	0.1	115.5	125	129	109.96%	9.96%
10	0	25	25	25	NA	NA

NA: Not applicable

5. Trueness:

5.1. Recovery Evaluation (Proportional error)

Recovery was assessed by adding different amounts of hCG to 3 samples with different concentrations. Then

hCG-titer ELISA KIT		IVD-REF: PS – HCG.T
		Ver. No: 02

% recovery and % bias was calculated for each sample using the equations below:

$$\%Recovery = \frac{Concentration\ Recovered}{Concentration\ Added} \times 100$$

$$\%Bias = \frac{\bar{X}\ measured - expected}{Expected} \times 100$$

The results are summarized in the following table:

Sample	Original concentration	50 IU/L added		100 IU/L added	
		Recovery %	Bias%	Recovery %	Bias%
1	45.8 IU/L	109.00%	4.70%	109.00%	6.17%
2	345 IU/L	95.00%	-0.63%	95.00%	-1.12%
3	725 IU/L	112.00%	0.77%	98.00%	-1.21%


5.2. Method comparison:

A comparison of the PSI hCG ELISA kit (y axis) using patient samples (n=88 range:7-941 IU/L) with ECLIA-Cobas E411-Roche diagnostic (HCG stat) (X axis) gave the following results by ordinary linear regression:

$$PSI\ hCG_{ELISA} = 0.9803\ hCG_{ECLIA} - 1.4885$$

$$r = 0.9988$$

$$r^2 = 0.9967$$

6. High dose hook effect

In PSI hCG EIA kit high dose hook effect has not been seen up to 1000000 IU/L


hCG-titer ELISA KIT		IVD-REF: PS – HCG.T
		Ver. No: 02

References:

1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Desirable Specifications for Total Error, Imprecision, and Bias,derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
5. Kasper D.L. et al. (2015). Harrison’s principles of internal medicine.19th ed.(page. 124e-5) Mc Graw Hill Education
6. Marcilac I. Troalen F. et al.(1992) Free human chorionic gonadotropin beta subunit in gonadal and nongonadal neoplasms. Canc Res. 1992;52:3901-3907.
7. Pagana K.D. et al. (2018). Mosby’s Manual of diagnostic and Laboratory tests.6th ed. Elsevier Inc.
8. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. Elsevier Inc.
9. Royal college of pathologists of Australasia. <https://www.datainnovations.com/allowable-total-error-table>

hCG-titer ELISA KIT		IVD-REF: PS – HCG.T
		Ver. No: 02

Procedure in-brief

Read at
450/630 nm