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1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Desirable Specifications for Total Error, Imprecision, and Bias, derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
5. Hoff rand A.W. et al. (2011). Postgraduate Haematology.6th ed (pp 28-29). Wiley Blackwell. John wiley & Sons, Ltd, publication
6. Kasper D.L. et al. (2015). Harrison's principles of internal medicine.19th ed. (pp. 1817). Mc Graw Hill Education
7. McPherson R.A. et al. (2017). HENRY'S clinical diagnosis and management by laboratory Methods. 23rd ed. (pp. 562). Elsevier Inc.
8. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests. (pp. 211-213). 6th ed. Elsevier Inc.
9. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. (pp. 756-9). Elsevier Inc.

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INTENDED USE

For the quantitative determination of Ferritin concentration in human serum or plasma.

CLINICAL SIGNIFICANCE

Ferritin is a soluble protein that along with insoluble hemosiderin are iron-storage proteins in human body. It consists of an approximately spherical apoprotein shell enclosing a core of ferric hydroxyphosphate which can accommodate up to 4000 iron atoms. Human Ferritin is made up from 24 subunits with molecular weight 2000 Dalton of two immunologically distinct types: H and L. The internal cavity of the Ferritin molecule communicates with the exterior via six channels, through which ferrous iron may enter to interact with a ferroxidase center on the Ferritin H subunit or leave after reduction, e.g. by dihydroflavins or ascorbic acid. Basic form of Ferritin is responsible for long-term iron-storage in liver, spleen and bone marrow. Acidic isoform is found in myocardium, placenta, tumor tissue and to a lesser extent in the depot organs. small amount of ferritin normally present in serum contains little iron and consists almost exclusively of L subunits. 1 ng/mL serum Ferritin corresponds to 8 mg total iron stores. And therefore serum Ferritin level measurement is the best, and most convenient laboratory test to estimate iron stores and diagnosis of iron deficiency and other Iron metabolisms disorders.

It has substituted the invasive and semiquantitative histochemical examination of bone marrow aspirate or biopsy as the gold standard for diagnosis of iron deficiency anemia.

Serum Ferritin measurement is a good indicator of iron storage in the body, but it does not reflect iron actually is available for erythropoiesis. Decrease in Ferritin below 15 ng/mL always indicates iron deficiency anemia and often occurs as a result of blood loss, decreased iron uptake, transferrin deficiency or increased iron demand (e.g. pregnancy). It is worthy to mention that, Ferritin is an acute phase reactant and can be increased in infections, acute and chronic inflammation, and malignant tumors; despite Iron deficiency. It also increases in viral and alcoholic hepatitis, chronic renal insufficiency without any relation to iron stores.

Diagnosis should be done based on entire clinical finding of individual patients.

PRINCIPLES OF THE METHOD

The PSI ferritin ELISA Kit is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators, controls and patient samples are incubated in streptavidin coated microplate together with a biotinylated Anti- ferritin monoclonal antibody and an HRP-conjugated monoclonal antibody, both recognize ferritin from different sites. Ferritin present in calibrators/controls or samples is reacted with both antibodies to form a sandwich complex and also adsorbed to the streptavidin coated microstrips via the biotinylated Anti-ferritin Mab. The strips are then washed and, buffered substrate/chromogen reagent (hydrogen peroxide and tetramethylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue color will develop and the intensity of the color is directly proportional to the amount of FERRITIN present in the samples.

Ferritin ELISA KIT		IVD-REF: PS - Ferr
		Ver. No: 02

REAGENTS PROVIDED

Reagent	Quantity/48 Test KIT	Quantity/96 Test KIT	Quantity/192 Test KIT	Reconstitution
Microplate (96 breakable wells) coated with streptavidin	48 wells	96 wells	192 wells	Ready to Use
Calibrators 1-7 (0, 10, 50, 100, 250, 500 and 1000 ng/mL) in buffer, containing preservative with traceability to reference material "WHO IRP 80/602"	7 Vial/0.5 ml	7 Vial/0.5 ml	7 Vial/1.0 ml	Ready to Use
Control serum in buffer compatible with human serum matrix with preservatives	2 Vial/0.5 ml	2 Vial/0.5 ml	2 Vial/1.0 ml	Ready to Use
Enzyme Conjugate (red color)	1 Vial/12.0 ml	2 Vial/12.0 ml	4Vial/12 ml	Ready to use
Concentrated Wash solution	1 Vial/30ml	1 Vial/30 ml	1 Vial/50 ml	Dilute 20X with distilled water
Substrate-Chromogenic solution (Tetramethylbenzidine and hydrogen peroxide)	1 Vial/6.0ml	1 Vial/12ml	2 Vial/12ml	Ready to Use
Stop solution (H2So4 0.12M) Additional information about other additive is considered proprietary.	1 Vial/6.0 ml	1 Vial/6.0 ml	1 Vial/12.0 ml	Ready to Use

Materials required but not supplied with kit

The following materials are required but not provided in the kit:

- Distilled water with conductivity < 1 μ S/cm
- Precision micropipettes for delivery of, 20, 50, 100, 200 microliters. An 8-channel pipette or resenser pipette with disposable plastic tips for delivery of 50-200 μ L is useful but not essential. Pipettes for dispensing milliliter volumes.
- Vortex mixer
- Automatic washer or any other device which can deliver 350 microliter wash solution for microplates
- Microplate reader with a wavelength of 450 nm and 620-650 nm, and an absorbance range of 0-3.0

STORAGE AND EXPIRATION DATING OF REAGENTS

- If kept at 2-8 °C, before opening or reconstitution, all kit components are stable until the expiry date as indicated on the label on the outside of the kit box.
- Calibrator concentrations are displayed on vial labels and may vary between lots.
- Do not freeze kits.
- Freshly prepared Working Wash solution should be used on the preparation day.
- Each lot of reagents and calibrators has been standardized to produce the correct reaction.

Do not interchange the reagents or calibrators between lots.

- Alterations in physical appearance of kit reagents may indicate instability or deterioration. Chromogenic solution TMB should be colorless. A blue color indicates that the reagent has been contaminated and should be discarded.

Ferritin ELISA KIT		IVD-REF: PS - Ferr
		Ver. No: 02

SPECIMEN COLLECTION AND PREPARATION

- Human serum and EDTA plasma are samples of choice for PSI Ferritin kit.
- Do not use grossly haemolysed or turbid samples.
- If the measurement cannot be done within 24 hours after sample collection, samples should be kept at 2-8°C up to 48 hours. For longer storage serum or plasma must be stored at -20°C.
- Avoid repeated freeze and thaw cycles. Frozen specimens must be gently but thoroughly mixed after thawing.

PRECAUTIONS AND WARNINGS

Safety

- For Professional in vitro diagnostic Use Only.
- All the reagents within the kit are calibrated for the direct determination of ferritin in human serum. The kit is not calibrated for the determination of ferritin in saliva, Non-EDTA plasma or other specimens of human or animal origin.
- Before starting the assay, read the working manual completely and carefully. Use the valid version of the instruction for use provided with the kit. Be sure that everything is understood.
- Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., USA) publication No (CDC) 88—8_95 on laboratory safety procedures or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Reagents contain sodium azide (NaN₃) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. Follow local guidelines for disposal of all waste material.
- Avoid any skin contact with all reagents, stop solution contains strong acid. In case of contact, wash thoroughly with soap and water.

Caution

Material used in the preparation of human source reagents has been tested and found to be non-reactive for HIV 1 & 2 antibodies, HCV antibody and hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

Preparation of reagents

Wash Buffer:

According to volume required, pour the desired volume of wash concentrate into a clean container and dilute 20 fold by adding adequate volume of distilled or deionized water to give a buffered Wash Solution. E.g. 10 mL wash concentrate to 190 mL distilled water. Wash solution should be prepared daily and discard unused working wash solution at the end of the day.

HANDELING NOTES

1. Do not use the kit or components beyond expiry date.
2. Do not mix up vials caps.
3. Bring all the reagents to room temperature prior to use.
4. Thoroughly mix all reagents, calibrators/controls and samples by gentle agitation or swirling.
5. Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
6. Use a clean container and high-quality distilled or deionized water to prepare the working wash solution.
7. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of

Ferritin ELISA KIT		IVD-REF: PS - Ferr
		Ver. No: 02

each reagent and sample.

8. For the dispensing of Chromogenic solution and Stop solution avoid pipettes with metal parts.
9. High precision pipettes or automated pipetting equipment will improve the precision.
10. Adhere the incubation times carefully. Tolerance limit of incubation time is $\pm 5\%$
11. Although PSI ferritin EIA kit uses streptavidin coated plate technology and reaction does not proceed before adding conjugate solution, but to avoid drift, arising from low volume sample evaporation or incubation time difference between 1st and last well it is strongly recommended not to run more than 6 strips when pipetting manually (Time delay).
12. Perform all steps consecutively without any delay. Any delay in each step has deleterious effects on final results.
13. Draw a calibration curve for each run, do not use stored curve from previous runs.
14. Incubation with Chromogenic solution, must be done in dark.
15. Only standard 0 may be used to dilute any high serum samples. The use of any other reagent may lead to false results.
16. Opened reagents are stable maximally for 8 weeks if they are not contaminated, and stored in resealed original containers and handled as prescribed. Return to 2-8°C immediately after use

Measurement Procedure

1. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminum pouch containing a desiccant and resealed carefully and store at 2-8°C.
2. Pipette 20 μL of each Calibrator (calibrators 0, 10, 50, 100, 250, 500, 1000 ng/mL), controls and patient samples in duplicate into the appropriate wells.
3. Pipette 200 μL of enzyme conjugate into each well. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid.
4. Gently mix for 15 seconds
5. Cover the plate with an adhesive seal and incubate for 45 minutes at room temperature (20-27°C).
6. Wash the plate according to procedure below.
 - a. Carefully dispense 350 μL of working wash solution into each well.
 - b. Aspirate the content of each well.
 - c. Repeat 4 additional wash for a total of 5 wash cycles.
 - d. After the end of last wash cycle invert the plate and tap it carefully against an absorbent paper.

Note: In automatic strip washer, follow the manufacturer's instructions for maintenance and wash the required number of wash cycles.

7. Add 100 μL of Chromogenic solution TMB to each well using the same pipetting technique as in procedural note 3. The Chromogenic solution TMB should be added to the wells as quickly as possible.
8. Incubate for 15 min at room temperature preferentially in the dark.
9. Add 50 μL of Stop Solution, mix and read the absorbance at 450 nm (with a reference filter 620, 630, or 650) in a microplate spectrophotometer within 10 min after addition of Stop Solution.

Important Note: always add stop solution in the same order as chromogenic solution TMB.

CALCULATION OF RESULTS

1. Construct a calibration curve by plotting the average O.D. (450 nm) for each calibrator on the vertical (Y) axis against the concentration on the horizontal (X) axis, and draw a best fit curve through the

Ferritin ELISA KIT

IVD-REF: PS - Ferr

Ver. No: 02


points on the graph.

- To determine concentration of Ferritin in each sample, first locate the O.D. value on the Y-axis and extend a horizontal line to the calibration curve. At the point of intersection, draw a vertical line to the X-axis and read the corresponding concentration.
- If a microplate spectrophotometer with built-in data calculation program is used, refer to the manual for the spectrophotometer and create a program using the concentration stated on the label of each of the Ferritin calibrators. For automatic calculation of Ferritin results it is recommended to use point-to-point method.

Typical Data

The following data are for illustration only and should never be used instead of the real-time calibration curve.

Row	OD	ng/mL
1	0.03	0
2	0.08	10
3	0.35	50
4	0.65	100
5	1.25	250
6	1.80	500
7	2.53	1000

**Quality control**

Ferritin control serums should be used for validation of the assay series. Ranges of expected results are indicated on the vial label. If assigned range are not obtained, a complete check of reagents and reader performance should be made and the assay shall be repeated. Each laboratory can also prepare its own serum pools at different levels, which can be used as internal controls in order to ensure the precision of the assay.

Expected values

A study with PSI ferritin EIA kit was performed using samples from some apparently healthy persons (ages 20-55) and also referring to some published scientific literature following reference interval corresponding to 2.5th and 97.5th percentiles of results was extrapolated. Each laboratory should first investigate the transferability of the recommended expected value to its own patient population and if it doesn't match, determine its own reference interval.:

Age groups	N	Reference interval (ng/mL)
Men: 20-55 y	125	20-250
Female: 20-55 y	142	10-120
Newborn	-	25-200
1 mo	-	20-600
2-5 mo	-	50-200
6 mo- 15 years	-	7-140

Ferritin ELISA KIT		IVD-REF: PS - Ferr
		Ver. No: 02

PERFORMANCE CRITERIA

1. Lower limit of measurement:

The Limit of Blank and Limit of Detection were determined according to CLSI (Clinical and Laboratory Standards Institute) EP17-A2 guideline.

In brief Limit of Blank was determined as the 95th percentile value from $n = 60$ measurements on standard 0 over several runs. The Limit of blank equals to the concentration that 95% of repeated measurements on zero standard give value below it. In the next step we determined standard deviation of repeated measurements on 3 samples with Ferritin contents less than 2.0 ng/mL determined from a separate method study (ECLIA Method). By adding the LOB calculated at 1st step to 1.645 times the greatest standards deviation calculated in the 2nd step ($LOD = LOB + 1.645SD_L$) LOD calculated as 2.0 ng/mL.

2. Precision:

Precision was determined using PSI Ferritin kit reagents and three pooled human sera according to EP05-A3 Guideline.

Three runs per day (morning, midday and afternoon) by 3 user in 10 working days ($10 \times 3 \times 3$). Results analyzed by fully nested ANOVA. Following results were obtained:

Sample Description	Mean (ng/mL)	Repeatability		Within-Laboratory Precision	
		SD	%CV	SD	%CV
Patient Pool	8.50	0.720	8.47	1.090	12.82
Patient Pool	78.50	4.62	5.89	5.34	6.80
Patient Pool	589.00	4.17	0.71	28.56	4.85

3. Specificity:

Specificity evaluation was performed according to guideline CLSI EP07-A2 to assess the cross-reactivity of the assay with human serum albumin, Alpha-fetoprotein and human transferrin. Samples containing cross-reactants was prepared by spiking natural human sera with material containing high levels of cross-reactants. Modification of serum matrix was no more than 10%. The % cross-reactivity was calculated for each sample using the equation below and normalized to the Ferritin content:

$$\% \text{cross-reactivity} = \frac{\text{mean conc. of spiked sample} - \text{mean conc. of unspiked sample}}{\text{spiked concentration}} \times 100\%$$

The results are summarized in the following table:

Cross-reactants	concentration	% Cross-reaction
Human serum albumin	500 µg/ml	0.08
Alpha- FP	500 µg/ml	0.13
Transferrin	500 µg/ml	0.06

Also interference evaluation from common sample abnormalities, hemolysis, Lipemia and icterus were performed according to aforementioned document and percentage interference was calculated for each sample using the equation below and normalized to the sample ferritin content.

Ferritin ELISA KIT		IVD-REF: PS - Ferr
		Ver. No: 02

$$\%interference = \frac{\text{measured value} - \text{true value}}{\text{true value}} \times 100\%$$

No interference were seen up to following listed concentration:

Interferent	source	Final concentration of interferent in test serum
Hemoglobin	Fresh haemolysed blood	50 mg/mL
Conjugated Bilirubin	Commercially crystalline Unconjugated Bilirubin	20 mg/dL
Triglycerides	Commercially Intralipid 10%	1000 mg/dL

4. Linearity

Evaluation of Linearity was performed according to CLSI (Clinical and Laboratory Standards Institute) EP 06-A guideline

In brief a sample with concentration of 813 ng/mL Ferritin as a high sample was diluted with another sample with 12 ng/mL in different proportions and their ferritin, measured in duplicate on each dilution. Measured values compared with expected values which estimated by dilution formula, then % bias and %recovery calculated for each concentration and following results are obtained:

No	Ratio	Expected (ng/mL)	Rep 1 (ng/mL)	Rep 2 (ng/mL)	recovery%	%Bias
1	1	813	813	813	NA	NA
2	0.9	726.6	731	739	101.16%	1.16%
3	0.8	647.2	664	656	101.98%	1.98%
4	0.7	567.8	534	541	94.66%	-5.34%
5	0.6	488.4	465	452	93.88%	-6.12%
6	0.5	409	416	435	104.03%	4.03%
7	0.4	329.6	318	337	99.36%	-0.64%
8	0.2	170.8	165	156	93.97%	-6.03%
9	0.1	91.4	96	102	108.32%	8.32%
10	0	12	12	12	NA	NA

NA: Not applicable

5. Trueness:

5.1. Recovery Evaluation (Proportional error)

Recovery was assessed by adding 100 ng/mL of Ferritin to 5 samples with different ferritin contents. Then % recovery and % bias was calculated for each sample using the equations below:

$$\%Recovery = \frac{\text{Concentration Recovered}}{\text{Concentration Added}} \times 100$$

$$\%Bias = \frac{\bar{X} \text{ measured} - \text{expected}}{\text{Expected}} \times 100$$

The results are summarized in the following table:

Ferritin ELISA KIT

IVD-REF: PS - Ferr

Ver. No: 02

Sample	Original Concentration (ng/mL)	Added ferritin (ng/mL)	Expected value (ng/mL)	Observed Value (ng/mL)	Recovery%	Bias%
1	21.15	100	121.15	136.75	115.6	12.88%
2	60.8		160.8	168.75	107.95	4.94%
3	129.15		229	228.7	99.55	-0.20%
4	225.1		325	337	112	3.61%
5	728		828	811	83	-2.05%


5.2. Method comparison:

A comparison of the PSI Ferritin ELISA kit (y axis) using patient samples (**n=88 range: 2.1-670 ng/mL**) with **ECLIA-Cobas E411-Roche diagnostic** (X axis) gave the following results by ordinary linear regression:

$$PSI \ EIA = 0.9235 \ CLIA + 3.3128$$

$$r = 0.9923$$

$$r^2 = 0.9847$$

**6. High dose hook effect**

In PSI Ferritin EIA kit high dose hook effect has not been seen up to 5000 ng/mL

Ferritin ELISA KIT		IVD-REF: PS - Ferr
		Ver. No: 02

References:

1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Desirable Specifications for Total Error, Imprecision, and Bias,derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
5. Hoffbrand A.W. et al. (2016). Postgraduate Haematology.7th ed (pp 28-29). Wiley Blackwell. John wiley & Sons, Ltd, publication
6. Kasper D.L. et al. (2015). Harrison's principles of internal medicine.19th ed.(pp. 627-628) Mc Graw Hill Education
7. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests.6th ed. Elsevier Inc.
8. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. Elsevier Inc.


Ferritin ELISA KIT

IVD-REF: PS - Ferr



Ver. No: 02

Procedure in-brief

Ferritin ELISA KIT		IVD-REF: PS - Ferr
		Ver. No: 02

SYMBOLS:

	BATCH NUMBER
	USE BY
	MANUFACTURER
	CONTAINS SUFFICIENT FOR <n> TESTS
	<i>IN VITRO</i> DIAGNOSTIC MEDICAL DEVICE
	TEMPERATURE LIMITATION
	CATALOGUE NUMBER
	CONSULT INSTRUCTIONS FOR USE