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$$\% \text{cross - reactivity} = \frac{\text{mean conc. of spiked sample} - \text{mean conc. of unspiked sample}}{\text{spiked conc.}} \times 100$$

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No	Ratio	Expected (ng/mL)	Rep 1 (ng/mL)	Rep 2 (ng/mL)	recovery%	%Bias
1	1	93	93	93	-	-
2	0.9	83.95	84	82.5	99.17%	0.83%
3	0.8	74.9	73.1	72.5	97.20%	2.80%
4	0.7	65.85	64.3	61.3	95.37%	4.63%
5	0.6	56.8	52.4	53.4	93.13%	6.87%
6	0.5	47.75	45.5	46.2	96.02%	3.98%
7	0.4	38.7	37.5	36.4	95.48%	4.52%
8	0.2	20.6	19.7	21.4	99.76%	0.24%
9	0.1	11.55	10.7	10.5	92.64%	7.36%
10	0	2.5	2.5	2.5	-	-

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Sample	Original concentration	10 ng/mL added		20 ng/mL added		40 ng/mL added	
		Recovery%	Bias%	Recovery%	Bias%	Recovery %	Bias%
1	9.25 ng/mL	112	6.23%	109	6.15%	108	7.95%
2	24.56 ng/mL	110	2.89%	95	-2.24%	107	5.88%
3	57.24 ng/mL	98	-0.30%	112	3.11%	93	-4.15%

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
1. Barth J.H et al.(2018). Observational studies on macroprolactin in a routine clinical laboratory. Clin Chem lab Med.26(56) 1-4.
2. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
3. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
4. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
5. Desirable Specifications for Total Error, Imprecision, and Bias, derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
6. McPherson R.A. et al.(2017).HENRY'S Clinical diagnosis and management by Laboratory Methods. 23rd ed. (pp. 363-364). Elsevier Inc.
7. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests.6th ed. Elsevier Inc.
8. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. Elsevier Inc.

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INTENDED USE

For the quantitative determination of Prolactin concentration in human serum or plasma.

CLINICAL SIGNIFICANCE

Prolactin is synthesized in anterior hypophysis and secreted in a circadian fashion. This hormone is made up of 198 amino acids with approximate molecular weight of 22-23 KDa. It is found in human serum in three different forms. Biologically active nonglycosylated monomeric form(Little) which is the predominant form of human serum prolactin and in a healthy woman consist of 80% of serum prolactin, followed by biologically inactive dimeric form (big form) and tetrameric(Big-big) form which has low biologic activity. Mammary gland is target organ for prolactin and its differentiation and development induced by this hormone. High prolactin concentration has inhibitory effect on steroids production in ovaries and hypophyseal gonadotropin secretion. During pregnancy, serum prolactin level increases considerably due to elevation of estrogens and progesterone concentrations.

Hyperprolactinemia (in men and women) is the main cause of fertility disorders. Stimulatory effect of prolactin on mammary gland leads to post-partum lactation. Prolactin measurement is helpful in diagnosis of oligomenorrhea (infrequent menstrual periods), amenorrhea, galctorrhea, loss of libido, impotence, primary hypothyroidism, polycystic ovary syndrome and azoospermia and hypogonadism in men. Extremely high levels of prolactin most often indicate pituitary adenoma, while other para-neoplastic tumors such as lung cancer can also act as an ectopic centre for prolactin secretion and may lead to hyperprolactinemia.

Because forms other than monomeric prolactin all react with immunoassays, they can result in falsely high PRL results in patients in whom a pathological elevation of PRL is not supported by CT or MRI and clinical findings. Various methods have been developed to eliminate this confusion, including measurement following polyethylene glycol extraction and centrifugal ultrafiltration.

PRINCIPLES OF THE METHOD

The PSI Prolactin ELISA kit is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators, controls and patient samples are incubated in anti-prolactin coated microplate together with a HRP-conjugated anti- Prolactin monoclonal antibody which recognize Prolactin from different sites. Prolactin present in calibrators/controls or samples is reacted with both capture and detection antibodies to form a sandwich complex. The strips are then washed and, buffered substrate/chromogen reagent (hydrogen peroxide and tetramethylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue color will develop and the intensity of the color is directly proportional to the amount of Prolactin present in the samples.

Prolactin ELISA KIT		IVD-REF: PS - PRL
		Ver. No: 02

REAGENTS PROVIDED

Reagent	Quantity/48 Test KIT	Quantity/96 Test KIT	Quantity/192 Test KIT	Reconstitution
Microplate (96 breakable wells) coated with Anti-prolactin Mab	48 wells	96 wells	192 wells	Ready to Use
Calibrators 1-6 (0, 5,10, 25, 50 and 100, ng/mL) in buffer, containing preservative with traceability to reference material "WHO ISO 84/500"	6 Vial/0.5 ml	6Vial/0.5 ml	6 Vial/1.0 ml	Ready to Use
Control serum in buffer compatible with human serum matrix with preservatives	1 Vial/0.5 ml	1 Vial/0.5 ml	1 Vial/1.0 ml	Ready to Use
Enzyme Conjugate (red color)	1 Vial/6.0 ml	1 Vial/12.0 ml	2Vial/12 ml	Ready to use
Concentrated Wash solution	1 Vial/30ml	1 Vial/30 ml	1 Vial/50 ml	Dilute 20X with distilled water
Substrate-Chromogenic solution (Tetramethylbenzidine and hydrogen peroxide)	1 Vial/6.0ml	1 Vial/12ml	2 Vial/12 ml	Ready to Use
Stop solution (H2So4 0.12M) Additional information about other additive is considered proprietary.	1 Vial/6.0 ml	1 Vial/6.0 ml	1 Vial/12.0 ml	Ready to Use

Materials required but not supplied with kit

The following materials are required but not provided in the kit:

- Distilled water with conductivity <math>< 1\mu\text{S}/\text{cm}^3</math>
- Precision micropipettes for delivery of 20, 50, 100 microliters. An 8-channel pipette or responder pipette with disposable plastic tips for delivery of 50-100 μL is useful but not essential.
- Pipettes for dispensing milliliter volumes.
- Vortex mixer
- Automatic washer or any other device which can deliver 350 microliter wash solution for microplates
- Microplate reader with a filter set at 450 nm and an upper OD limit of 3.0 or greater* (see measurement procedure step 9).

STORAGE AND EXPIRATION DATING OF REAGENTS

- If kept at 2-8 °C, before opening or reconstitution, all kit components are stable until the expiry date as indicated on the label on the outside of the kit box.
- Freshly prepared Working Wash solution should be used on the preparation day.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration of reagents.
- The substrate solution (TMB) is sensitive to light and should remain colorless if properly stored. Instability or contamination may be indicated by the development of a blue color, in which case it should not be used.

SPECIMEN COLLECTION AND PREPARATION

- Human serum and EDTA plasma are samples of choice for PSI PRL kit.
- If the measurement cannot be done within 24 hours after sample collection, samples should be kept at 2-8 °C up to 48 hours. For longer storage, serum or plasma must be stored at -20°C.

Prolactin ELISA KIT		IVD-REF: PS - PRL
		Ver. No: 02

Avoid repeated freeze and thaw cycles. Frozen specimens must be gently but thoroughly mixed after thawing

PRECAUTIONS AND WARNINGS

- For Professional in vitro diagnostic Use Only
- All the reagents within the kit are calibrated for the direct determination of Prolactin in human serum. The kit is not calibrated for the determination of prolactin in urine, Non-EDTA plasma or other specimens of human or animal origin.
- Before starting the assay, read this instruction completely and carefully. Use the valid version of the instruction provided with the kit. Be sure that everything is understood
- Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., USA) publication No (CDC) 88—8_95 on laboratory safety procedures or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Reagents contain sodium azide (NaN₃) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. Follow local guidelines for disposal of all waste material.
- Avoid any skin contact with all reagents; stop solution contains strong acid. In case of contact, wash thoroughly with soap and water.

Caution

- Material used in the preparation of human source reagent has been tested and found to be non-reactive for HIV 1 & 2 antibodies, HCV antibody and hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

Preparation of reagents

Wash Buffer:

According to volume required, pour the desired volume of wash concentrate into a clean container and dilute 20 fold by adding adequate volume of distilled or deionized water to give a buffered Wash Solution. E.g. 10 mL wash concentrate to 190 mL distilled water. Wash solution should be prepared daily and discard unused working wash solution at the end of the day.

HANDELING NOTES:

- Do not use the kit or components beyond expiry date.
- Do not mix up vials caps.
- Do not mix reagents from different kit lots.
- All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use.
- A calibrator curve must be established for every run.
- Calibrators, controls and specimen samples should be assayed in duplicate. Vertical alignment is recommended.
- Use a clean container and high quality distilled or deionized water to prepare the Working Wash solution.
- In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
- For the dispensing of Chromogenic and Stop solution, avoid pipettes with metal parts.

Prolactin ELISA KIT		IVD-REF: PS - PRL
		Ver. No: 02

- High precision pipettes or automated pipetting equipment will improve the precision.
- Adhere the incubation times carefully. Tolerance limit of incubation time is $\pm 5\%$
- To avoid drift, arising from low volume sample evaporation or incubation time difference between first and last well, it is strongly recommended not to run more than six strips when pipetting manually (Time delay).
 - Perform all steps consecutively without any delay. Any delay in each step has deleterious effects on measurement results.
 - Draw a calibration curve for each run, do not use stored curve from previous runs.
 - Incubation with Chromogenic solution must be done in dark.
 - Opened reagents are stable maximally for 8 weeks if they are not contaminated, and stored in resealed original containers and handled as prescribed. Return to 2-8°C immediately after use
 - Kit reagents must be regarded as hazardous waste and disposed of according to national regulations

Measurement Procedure

1. Remove the required number of well strips. Reseal the bag and return any unused strips to the refrigerator.
2. Pipette 20 μL of each Calibrator (calibrators 0, 10, 50, 100, 250, 500, 1000 ng/mL), controls and patient samples in duplicate into the appropriate wells.
3. Pipette 100 μL of enzyme conjugate into each well. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid.
4. Gently mix for 30 seconds
5. Cover the plate with an adhesive seal and incubate for 45 minutes at room temperature (20-27°C).
6. Wash the plate according to procedure below.
 - a. Carefully dispense 350 μL of working wash solution into each well.
 - b. Aspirate the content of each well.
 - c. Repeat 4 additional wash for a total of 5 wash cycles.
 - d. After the end of last wash cycle invert the plate and tap it carefully against an absorbent paper.

Note: Automatic strip washer: Follow the manufacturer's instructions for maintenance and wash the required number of wash cycles.

7. Add 100 μL of Chromogenic solution TMB to each well using the same pipetting technique as in procedural note 3. The Chromogenic solution TMB should be added to the wells as quickly as possible.
8. Incubate for 15 min at room temperature preferentially in the dark.
9. Add 50 μL of Stop Solution, mix and read the absorbance at 450 nm (with a reference filter 620, 630, or 650) in a microplate spectrophotometer within 10 min after addition of Stop Solution.

Important Note: always add stop solution in the same order as chromogenic solution TMB.

CALCULATION OF RESULTS

1. Construct a calibration curve by plotting the average O.D. (450 nm) for each calibrator on the vertical (Y) axis against the concentration on the horizontal (X) axis, and draw a best fit curve through the points on the graph.
2. To determine concentration of Prolactin in each sample, first locate the O.D. value on the Y-axis and extend a horizontal line to the calibration curve. At the point of intersection, draw a vertical line to the X-axis and read the corresponding concentration.
3. If a microplate spectrophotometer with built-in data calculation program is used, refer to the manual


Prolactin ELISA KIT		IVD-REF: PS - PRL
		Ver. No: 02

for the spectrophotometer and create a program using the concentration stated on the label of each of the Prolactin calibrators. For automatic calculation of Prolactin results it is recommended to use point-to-point method.

Typical Data

The following data are for illustration only and should never be used instead of the real-time calibration curve.

OD	ng/mL
0.029	0
0.177	5
0.350	10
0.840	25
1.550	50
2.700	100

Quality control

Prolactin control serum should be used for validation of the assay series. Ranges of expected results are indicated on the vial label. If assigned range are not obtained, a complete check of reagents and reader performance should be made and the assay shall be repeated. Each laboratory can also prepare its own serum pools at different levels, which can be used as internal controls in order to ensure the precision of the assay.

Reference interval

A study with PSI Prolactin EIA kit was performed using samples from some apparently healthy persons (ages 20-55), following reference interval corresponding to 2.5th and 97.5th percentiles of results was extrapolated. Each laboratory should first investigate the transferability of the recommended expected value to its own patient population and if it doesn't match, determine its own reference interval.:

Age sex group	N	Reference interval
Men	100	1.8-21.4 ng/ml
Women-pre-menopausal	180	4.6-29.5 ng/ml
Women-post-menopausal	130	1.5-18.5 ng/ml

PERFORMANCE CRITERIA

1. Lower limit of measurement:

The Limit of Blank and Limit of Detection were determined according to CLSI (Clinical and Laboratory Standards Institute) EP17-A2 guideline.

In brief Limit of Blank was determined as the 95th percentile value from n = 60 measurements on standard 0 over several runs. The Limit of blank equals to the concentration that 95% of repeated measurements on zero

Prolactin ELISA KIT		IVD-REF: PS - PRL
		Ver. No: 02

standard give value below it. In the next step we determined standard deviation of repeated measurements on 3 samples with Prolactin contents less than 1.5 ng/mL determined from a separate method study (ECLIA Method). By adding the LOB calculated at 1st step to 1.645 times the greatest standards deviation calculated in the 2nd step($LOD = LOB + 1.645SD_L$) LOD calculated as 1.0 ng/mL.

2. Precision:

Precision was determined using PSI Prolactin Kit reagents and two pooled human sera and a commercial control material according to EP05-A3 Guideline.

Three runs per day (morning, midday and afternoon) by 3 user in 10 working days($10 \times 3 \times 3$). Results analyzed by fully nested ANOVA. Following results were obtained:

Sample Description	Mean (ng/mL)	Repeatability		Within-Laboratory Precision	
		SD	%CV	SD	%CV
Patient Pool	8.521	0.69	8.14	0.77	8.98
Patient Pool	24.645	1.25	5.05	1.54	6.25
L3-Biorad	59.641	3.21	5.39	4.25	7.12

3. Specificity:

Specificity evaluation was performed according to guideline CLSI EP07-A2 to assess the cross-reactivity of the assay with following cross-reactants. Samples containing cross-reactants was prepared by spiking natural human sera with material containing high levels of cross-reactants. Modification of serum matrix was no more than 10%. The % cross-reactivity was calculated for each sample using the equation below and normalized to the Prolactin content:

$$\%cross - reactivity = \frac{\text{mean conc. of spiked sample} - \text{mean conc. of unspiked sample}}{\text{spiked concentration}} \times 100\%$$

The results are summarized in the following table:

Cross-reactants	Concentration	% cross reactivity
hTSH	500 μ IU/ml	0.13
hLH	500 mIU/ml	0.17
hFSH	500 mIU/ml	0.04

Interference evaluation from common sample abnormalities, hemolysis, Lipemia and icterus were also performed according to aforementioned document and percentage interference calculated using following formula:

$$\%interference = \frac{\text{measured value} - \text{true value}}{\text{true value}} \times 100\%$$

No interference were seen up to following listed concentration:

Interferent	source	Final concentration of interferent in test serum
Hemoglobin	Fresh haemolysed blood	50 mg/mL
Conjugated Bilirubin	Commercially crystalline Unconjugated Bilirubin	20 mg/dL

Prolactin ELISA KIT		IVD-REF: PS - PRL
		Ver. No: 02

Triglycerides	Commercially Intralipid 10%	1000 mg/dL
---------------	-----------------------------	-------------------

4. Linearity

Evaluation of Linearity was performed according to CLSI (Clinical and Laboratory Standards Institute) EP 06-A guideline

In brief a sample with concentration of 93 ng/mL Prolactin as a high sample was diluted with another sample with 2.5 ng/mL in different proportions and their Prolactin, measured in duplicate on each dilution.

Measured values compared with expected values which estimated by dilution formula, then % bias and %recovery calculated for each concentration and following results are obtained:

No	Ratio	Expected (ng/mL)	Rep 1 (ng/mL)	Rep 2 (ng/mL)	recovery%	%Bias
1	1	93	92	94	NA	NA
2	0.9	83.95	84	82.5	99.17%	-0.83%
3	0.8	74.9	73.1	72.5	97.20%	-2.80%
4	0.7	65.85	64.3	61.3	95.37%	-4.63%
5	0.6	56.8	52.4	53.4	93.13%	-6.87%
6	0.5	47.75	45.5	46.2	96.02%	-3.98%
7	0.4	38.7	37.5	36.4	95.48%	-4.52%
8	0.2	20.6	19.7	21.4	99.76%	-0.24%
9	0.1	11.55	10.7	1.0.5	92.64%	-7.36%
10	0	2.5	2.1	2.9	NA	NA

NA: Not applicable

5. Trueness:

6. 5.1. Recovery Evaluation (Proportional error)

Recovery was assessed by adding three different amount of prolactin to three sample with different prolactin contents. Then % recovery and % bias was calculated for each sample using the equations below:

$$\%Recovery = \frac{\text{Concentration Recovered}}{\text{Concentration Added}} \times 100$$

$$\%Bias = \frac{\bar{X} \text{ measured} - \text{expected}}{\text{Expected}} \times 100$$

The results are summarized in the following table:

Sample	Original concentration	10 ng/mL Ladded		20 ng/mL Ladded		40 ng/mL Ladded	
		Recovery%	Bias%	Recovery%	Bias%	Recovery%	Bias%
1	9.25 ng/mL	112	6.23%	109	6.15%	108	7.95%
2	24.56 ng/mL	110	2.89%	95	-2.24%	107	5.88%
3	57.24 ng/mL	98	-0.30%	112	3.11%	93	-4.15%

Prolactin ELISA KIT		IVD-REF: PS - PRL
		Ver. No: 02


5.2. Method comparison:

A comparison of the PSI Prolactin ELISA kit (y axis) using patient samples (n=56) with **ECLIA-Cobas E411-Roche diagnostic** (X axis) gave the following results by ordinary linear regression:

$$PSI = 0.9493ECLIA - 0.1558$$

$$r = 0.9753$$

$$r^2 = 0.9513$$

7. High dose hook effect

In PSI Prolactin EIA kit high dose hook effect has not been seen up to 4000 ng/mL

Prolactin ELISA KIT		IVD-REF: PS - PRL
		Ver. No: 02

References:

1. Barth J.H et al.(2018). Observational studies on macroprolactin in a routine clinical laboratory. Clin Chem lab Med.26(56) 1-4.
2. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
3. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
4. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
5. Desirable Specifications for Total Error, Imprecision, and Bias,derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
6. Kasper D.L. et al. (2015). Harrison's principles of internal medicine.19th ed. Mc Graw Hill Education
7. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests.6th ed. Elsevier Inc.
8. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. Elsevier Inc.
9. Royal college of pathologists of Australasia. <https://www.datainnovations.com/allowable-total-error-table>

Prolactin ELISA KIT

IVD-REF: PS - PRL

Ver. No: 02

Procedure in-brief