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1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline. 2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline. 2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Helen Mason. et al. (2014). The Physiology and clinical utility of anti-Müllerian hormone in women. Human Reproduction update. 0(0) pp 1-16
5. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics. 6th ed. Elsevier Inc. 1644
6. Scott M.G. et al. (1989). Hormonal evaluation of female infertility and reproductive disorders. Clin. Chem. 1989. 35:620-3

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INTENDED USE

For the quantitative determination of **Anti-Müllerian Hormone (AMH) or Müllerian-inhibiting hormone(MIH)** concentration in human serum or plasma.

CLINICAL SIGNIFICANCE

AMH is a glycoprotein belonging to a family of growth and differentiation factors called transforming growth factor beta and has been primarily known for its role in male sexual differentiation. All members of this superfamily are involved in the regulation of tissues growth and differentiation. During male fetal differentiation, it is produced by Sertolli cells in testis and is responsible for regression of Müllerian ducts during male fetal sex differentiation. Secretion of AMH by Sertolli cells starts during embryogenesis and continues throughout the life.

AMH is produced as a precursor protein, consisting of 70 KDa monomers. Proteolytic cleavage yields a 55 KDa N-terminal proregion and a 12.5 KDa C-terminal mature region. The pro- and mature homodimers remain non-covalently associated, resulting in a 140 KDa complex in circulation. The mature region of AMH holds the biological activity of the protein, but in contrast to other TGF- β family members, requires the N-terminal proregion to obtain its full activity. It has been suggested that the proregion is involved in protein stability and folding. AMH Type II receptor can only bind to biologically active form of AMH.

In females, AMH plays an important role in the ovarian folliculogenesis. Follicle development in the ovaries comprises two distinct stages: initial recruitment, by which primordial follicles start to mature, and cyclic recruitment, which leads to harmonized growth of small antral follicles, one of them is then selected to ovulate. Cyclic recruitment is under the control of FSH.

AMH expression in granulosa cells starts in primary follicles and in preantral and small antral follicles remains high until a follicle diameter reaches a diameter of around 8 mm. A sharp decline in AMH expression which is seen around 8 mm corresponds to selection of follicle for dominance. This pattern of AMH expression indicates inhibitory role of AMH at two distinct stages of folliculogenesis. First, AMH blocks the transition of follicles from primordial into maturation stages and so plays an important role in regulating the reservoir of primordial follicles. Second, AMH decreases follicle sensitivity to FSH and therefore has a role in the process of follicular selection.

Serum levels of AMH are rarely detectable at birth in females, reach their highest levels after puberty, decrease progressively thereafter with age, and become undetectable at menopause. Serum AMH levels are relatively stable through menstrual cycle with substantial fluctuations being observed in younger women. Serum AMH levels decrease significantly during the use of combined contraceptives.

Serum AMH determination has a great variety of clinical applications. It is mainly used for estimation of ovarian reserve indicating the number of antral and pre-antral follicles, (an index called antral follicle count or AFC). It provides useful information regarding ovarian follicular status and responsiveness to controlled ovarian hyperstimulation. AMH measurement is also used for monitoring of iatrogenic complications such as effects of radiotherapy, chemotherapy and surgery on ovarian reserve. Diagnosis of disorders of sex development in children and monitoring of granulosa cell tumors to detect remaining or recurrent disease are further clinical applications of AMH measurement. AMH increases in polycystic ovary syndrome (PCOS) which affects 5-10% of women. AMH is a good marker for predicting time to menopause.

AMH Gen II ELISA KIT		IVD-REF: PS - AMH
		Ver. No: 02

PRINCIPLES OF THE METHOD

The PSI AMH second generation ELISA kit is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. In the first step calibrators, controls and patient samples are incubated in a monoclonal anti-AMH coated plate. Strips are then washed to remove unbound matrix. At the next step a biotinylated detection antibody which recognizes AMH from different region is added to each well. After incubation and washing HRP-conjugated streptavidin is added to each well.

During these phases, a complex consisting of solid phase capture Ab-AMH-detection biotinylated antibody-HRP conjugated streptavidin is formed. The strips are then washed and, buffered substrate/chromogen reagent (hydrogen peroxide and tetramethylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction, a blue color will develop and the intensity of the color is directly proportional to the amount of AMH present in the samples.

REAGENTS PROVIDED

Reagent	Quantity/48 Test KIT	Quantity/96 Test KIT	Reconstitution
Microplate (96 breakable wells) coated with monoclonal anti-AMH	48 wells	96 wells	Ready to Use
Calibrators 1-7 (0, 0.2, 1.0, 2.5, 5.0, 10 and 20 ng/mL) in buffer, containing preservative	7 Vial/1.0 ml	7 Vial/1.5 ml	Ready to Use
Control serum in buffer compatible with human serum matrix with preservatives	2 Vial/1.0 ml	2 Vial/1.5 ml	Ready to Use
Biotin-Reagent(Blue color)	1 Vial/6.0 ml	1 Vial/12.0 ml	Ready to Use
Conjugate-Reagent(Red color)	1 Vial/6.0 ml	1 Vial/12.0 ml	Ready to use
Concentrated Wash solution	1 Vial/30ml	1 Vial/30 ml	Dilute 1:20 with distilled water
Substrate-Chromogenic solution (Tetramethylbenzidine and hydrogen peroxide)	1 Vial/6.0ml	1 Vial/12ml	Ready to Use
Stop solution (H ₂ SO ₄ 0.12M) Additional information about the other additive is considered proprietary.	1 Vial/6.0 ml	1 Vial/6.0 ml	Ready to Use

Materials required but not supplied with kit

The following materials are required but not provided in the kit:

- Distilled water with conductivity < 1µS/cm
- Precise micropipettes for delivery of 50, 100 microliters. An 8-channel pipette or responder pipette with disposable plastic tips for delivery of 50-200 µL is useful but not essential.
- Pipettes for dispensing milliliter volumes.
- Vortex mixer
- Automatic washer or any other device which can deliver 350 microliter wash solution for microplates

AMH Gen II ELISA KIT		IVD-REF: PS - AMH
		Ver. No: 02

- Microplate reader with a wavelength of 450 nm and 620-650 nm, and an absorbance range of 0-3.0

STORAGE AND EXPIRATION DATING OF REAGENTS

- If kept at 2-8 °C, before opening or reconstitution, all kit components are stable until the expiry date as indicated on their label.
- Calibrators concentration is displayed on the vial labels and may vary between lots.
- Do not freeze kits.
- Freshly prepared Working Wash solution should be used on the preparation day.
- Each lot of reagents and calibrators has been standardized to produce the correct reaction. Do not interchange the reagents or calibrators between lots.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration. Chromogenic solution TMB should be colorless. Blue color indicates that the reagent has been contaminated and should be discarded.

SPECIMEN COLLECTION AND PREPARATION

- Human serum and EDTA plasma are samples of choice for PSI AMH kit.
- Do not use grossly haemolysed or turbid samples.
- If the measurement cannot be done within 24 hours after sample collection, samples should be kept at 2-8°C up to 48 hours. For longer storage serum or plasma must be stored at -20°C.
- Avoid repeated freeze and thaw cycles. Frozen specimens must be gently but thoroughly mixed after thawing.

PRECAUTIONS AND WARNINGS

Safety

- For professional in vitro diagnostic use only.
- All the reagents within the kit are calibrated for the direct determination of AMH in human serum. The kit is not calibrated for the determination of AMH in Non-EDTA plasma or other specimens of human or animal origin.
- Before starting the assay, read the working manual completely and carefully. Use the valid version of the instruction provided with the kit. Be sure that everything is understood.
- Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., USA) publication No (CDC) 88—8_95 on laboratory safety procedures or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Reagents contain sodium azide (NaN₃) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. Follow local guidelines for disposal of all waste material.
- Avoid any skin contact with all reagents. Stop solution contains strong acid in case of contact, wash thoroughly with soap and water.

AMH Gen II ELISA KIT		IVD-REF: PS - AMH
		Ver. No: 02

Caution

Material used in the preparation of human source reagents has been tested and found to be non-reactive for HIV 1 & 2 antibodies, HCV antibody and hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they are potentially infectious.

Preparation of reagents

Wash Buffer:

According to the required volume, pour the desired volume of wash concentrate into a clean container and dilute 20 fold by adding adequate volume of distilled or deionized water to gain buffered Wash Solution. E.g. 10 mL wash concentrate to 190 mL distilled water. Wash solution should be prepared daily and discard unused working wash solution at the end of the day.

HANDELING NOTES

1. Do not use the kit or components beyond the expiry date.
2. Do not mix up vial caps.
3. Bring all the reagents to room temperature prior to use.
4. Thoroughly mix all reagents, calibrators/controls and samples by gentle agitation or swirling.
5. Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
6. Use a clean container and high-quality distilled or deionized water to prepare the working wash solution.
7. In order to avoid cross-contamination, use a clean disposable pipette tip for addition of each reagent and sample.
8. For the dispensing of Chromogenic solution and Stop solution, avoid pipettes with metal parts.
9. High precision pipettes or automated pipetting equipment will improve the precision.
10. Adhere the incubation times carefully. Tolerance limit of incubation time is $\pm 5\%$
11. To avoid drift, arising from incubation time difference between first and last well, it is strongly recommended not to run more than six strips when pipetting manually (Time delay).
12. Perform all steps consecutively without any delay. Any delay in each step has deleterious effects on final results.
13. Draw a calibration curve for each run, do not use stored curve from previous runs.
14. Incubation with Chromogenic solution, must be done in the dark.
15. Only standard 0 may be used to dilute any high serum samples. The use of any other reagent may lead to false results.
16. Opened reagents are stable maximally for 8 weeks if they are not contaminated, and stored in resealed original containers and handled as prescribed. Return to 2-8°C immediately after use.

Measurement Procedure

1. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminum pouch containing a desiccant, reseal carefully and store at 2-8°C.
2. Pipette 100 μL of each Calibrator (calibrators 0, 0.2, 1.0, 2.5, 5.0, 10 & 20 ng/mL), controls and patient samples in duplicate into the appropriate wells.

AMH Gen II ELISA KIT		IVD-REF: PS - AMH
		Ver. No: 02

3. Cover the plate with an adhesive seal and incubate for 60 minutes at room temperature (20-27°C).
4. Wash the plate according to procedure below.
 - i. Carefully dispense 350 µL of working wash solution into each well.
 - ii. Aspirate the content of each well.
 - iii. Repeat 4 additional wash for a total of 5 wash cycles.
 - iv. After the last wash cycle invert the plate and tap it carefully against an absorbent paper.

Note: Automatic strip washer: Follow the manufacturer's instructions for maintenance and wash the required number of wash cycles.

5. Add 100 µL of Biotin-regent to each well and incubate for 30 minutes at room temperature (20-27°C).
6. Wash the strips as described in procedural note 4.
7. Add 100 µL of conjugate-regent to each well and incubate for 15 minutes at room temperature (20-27°C).
8. Wash the strips as described in procedural note 4.
9. Add 100 µL of Chromogenic solution TMB to each well. The Chromogenic solution TMB should be added to the wells as quickly as possible.
10. Incubate for 15 min at room temperature preferentially in the dark.
11. Add 50 µL of Stop Solution, mix and read the absorbance at 450 nm (with a reference filter 620, 630, or 650 nm) in a microplate spectrophotometer within 10 min after addition of Stop Solution.

Important Note: always add stop solution in the same order as chromogenic solution TMB.

CALCULATION OF RESULTS


1. Construct a calibration curve by plotting the average O.D. (450 nm) for each calibrator on the vertical (Y) axis against the concentration on the horizontal (X) axis, and draw the best fit curve through the points on the graph.
2. To determine concentration of AMH in each sample, first locate the O.D. value on the Y-axis and extend a horizontal line to the calibration curve. At the point of intersection, draw a vertical line to the X-axis and read the corresponding concentration.
3. If a microplate spectrophotometer with built-in data calculation program is used, refer to the manual for the spectrophotometer and create a program using the concentration stated on the label of each calibrator. For automatic calculation of AMH results, it is recommended to use point-to-point method.

Typical Data

The following data are for illustration only and should never be used instead of the real-time calibration curve.

AMH Gen II ELISA KIT		IVD-REF: PS - AMH
		Ver. No: 02

Row	OD	ng/mL
1	10.02	0.0
2	20.10	0.2
3	30.26	1.0
4	40.52	2.5
5	50.94	5.0
6	61.65	10.0
7	72.78	20.0

Quality control

AMH control serum should be used for validation of the assay series. Range of expected results are indicated on the vial label. If assigned range is not obtained, a complete check of reagents and reader performance should be made and the assay shall be repeated. Each laboratory can also prepare its own serum pools at different levels, which can be used as internal controls in order to ensure the precision of the assay.

Expected values

A study with PSI AMH EIA kit was performed using samples from some apparently healthy people from different sex and age groups following reference interval corresponding to 2.5th and 97.5th percentiles of results was extrapolated. Each laboratory should first investigate the transferability of the recommended expected value to its own patient population and if it doesn't match, determine its own reference interval:

Age group	Median(ng/mL)	2.5-97.5 th percentile (ng/mL)
Males(>12 years)	5.7	0.7-19
Boys(24 mo-12 y)	60	7.4-243
Females(13-45 years)	2.4	1.0-4.0
Girls(24 mo-12 y)	1.3	<8.9
Girls(<24 mo)	-	<4.7
Females(diminished ovarian reserve)	-	<1.0
Females(Post-Menopausal)	-	ND

PERFORMANCE CRITERIA

1. Lower limit of measurement:

The Limit of Blank and Limit of Detection were determined according to CLSI (Clinical and Laboratory Standards Institute) EP17-A2 guideline.

In brief, Limit of Blank was determined as the 95th percentile value from n = 60 measurements on standard 0 over several runs. The Limit of blank equals to the concentration that 95% of the repeated measurements on zero standard give value below it. In the next step, we determined standard deviation of the repeated measurements on 3 samples with AMH contents less than 0.2 ng/mL determined from a separate method study (ECLIA Method). By adding the LOB calculated at 1st step to 1.645 times the greatest standards deviation calculated in the 2nd step ($LOD = LOB + 1.645SD_L$) LOD calculated as 0.03 ng/mL.

AMH Gen II ELISA KIT		IVD-REF: PS - AMH
		Ver. No: 02

2. Precision:

Precision was determined using PSI AMH kit reagents and three pooled human sera according to EP05-A3 Guideline.

Three runs per day (morning, midday and afternoon) by 3 user in 10 working days (10 × 3 × 3). Results analyzed by fully nested ANOVA. Following results were obtained:

Sample Description	Mean (ng/mL)	Repeatability		Within-Laboratory Precision	
		SD	%CV	SD	%CV
Patient Pool	0.87	0.048	5.22	0.035	4.02
Patient Pool	1.70	0.06	3.53	0.074	4.35
Patient Pool	6.37	0.34	5.38	0.404	6.34

3. Specificity:

Specificity evaluation was performed according to guideline CLSI EP07-A2 to assess the cross-reactivity of the assay with common cross-reactants usually encountered in human serum. Samples containing cross-reactants was prepared by spiking natural human sera with material containing high levels of cross-reactants. The modification of serum matrix was no more than 10%. The % cross-reactivity was calculated for each sample using the equation below and normalized to the AMH content:

$$\% \text{cross-reactivity} = \frac{\text{mean conc. of spiked sample} - \text{mean conc. of unspiked sample}}{\text{spiked concentration}} \times 100\%$$

The results are summarized in the following table:

Cross reactants	Concentration	Interference%
hLH	100 mIU/mL	0.04
hFSH	100 mIU/mL	0.06
TGF β1	100 ng/mL	0.06
Activin A	100 ng/mL	0.07
Inhibin A	20 ng/mL	0.04

PSI AMH Gen II ELISA measurement method is unaffected by icterus (Bilirubin up to 40 mg/dL) Lipemia (Intralipid up to 1000 mg/dL) and hemolysis (Hb up to 50 mg/mL)

4. Linearity

Evaluation of Linearity was performed according to CLSI (Clinical and Laboratory Standards Institute) EP 06-A guideline

In brief, a sample with concentration of 18.5 ng/mL AMH as a high sample was diluted with another sample with 0.07 ng/mL in different proportions and their AMH was measured in duplicate on each dilution. Measured values were compared with expected values which were estimated by dilution formula, then % bias and %recovery were calculated for each concentration. The following results are obtained:

AMH Gen II ELISA KIT		IVD-REF: PS - AMH
		Ver. No: 02

No	Ratio	Expected (ng/mL)	Rep 1 (ng/mL)	Rep 2 (ng/mL)	recovery%	%Bias
1	1	18.50	18.2	18.8	NA	NA
2	0.9	16.66	17.3	17.5	104%	4.5%
3	0.8	14.81	16.4	16.5	111%	11.0%
4	0.7	12.97	14	13.5	106%	6.0%
5	0.6	11.13	12.4	11.9	109%	9.2%
6	0.5	9.28	9.5	8.7	98%	-2.0%
7	0.4	7.44	6.9	7.5	97%	-3.2%
8	0.3	5.60	5.5	5.2	96%	-4.4%
9	0.2	3.75	3.1	3.9	93%	-6.8%
10	0.1	1.91	1.8	1.75	93%	-7.1%
11	0	0.07	0.047	0.087	NA	NA

NA: Not applicable

5. Trueness:

5.1. Recovery Evaluation (Proportional error)

Recovery was assessed by adding 10 ng/mL (using standard 20 ng/mL as source of analyte) of AMH to six samples with different AMH contents. Then % recovery and % bias was calculated for each sample using the equations below:

$$\%Recovery = \frac{Concentration\ Recovered}{Concentration\ Added} \times 100$$

$$\%Bias = \frac{\bar{X}\ measured - expected}{Expected} \times 100$$

The results are summarized in the following table:

Sam. No	Base Conc. ng/mL	Added Analyte	Rep 1 ng/mL	Rep 2 ng/mL	Bias	Recovery
1	0.045	10 ng/mL	10.18	9.35	-2.79%	97.20%
2	1.15		12.31	10.54	2.47%	102.75%
3	2.7		11.68	11.87	-7.28%	90.75%
4	4.5		13.74	13.64	-5.59%	91.90%
5	5.6		15.45	16.45	2.24%	103.50%
6	7.2		17.53	18.47	4.65%	108.00%

AMH Gen II ELISA KIT		IVD-REF: PS - AMH
		Ver. No: 02


5.2. Method comparison:

A comparison of the PSI AMH Gen II ELISA kit (y axis) using patient samples (**n=89 range:0.07-18.5ng/mL**) with EIA BECKMAN COULTER AMH Gen II (X axis) gave the following results by ordinary linear regression:

$$PSI\ EIA = 0.9113\ BC + 0.1149$$

$$r = 0.9892$$

$$r^2 = 0.9785$$

6. High dose hook effect

In PSI AMH Gen II ELISA kit, high dose hook effect has not been seen up to 2000 ng/mL.

References:

1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Helen Mason. et al. (2014). The Physiology and clinical utility of anti-Müllerian hormone in women. Human Reproduction update. 0(0) pp 1-16
5. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. Elsevier Inc.pp 1644
6. Scott M.G. et al. (1989). Hormonal evaluation of female infertility and reproductive disorders. Clin. Chem. 1989. 35:620-30

AMH Gen II ELISA KIT		IVD-REF: PS - AMH
		Ver. No: 02

Procedure in-brief