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1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Desirable Specifications for Total Error, Imprecision, and Bias, derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
5. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests.6th ed. Elsevier Inc.
6. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. Elsevier Inc.
7. Royal college of pathologists of Australasia. <https://www.datainnovations.com/allowable-total-error-table>
8. Scott M.G. et al. (1989). Hormonal evaluation of female infertility and reproductive disorders. Clin. Chem. 1989. 35:620-30

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INTENDED USE:

The Quantitative Determination of Follicle Stimulating hormone (FSH) Concentration in Human Serum.

CLINICAL SIGNIFICANCE:

Follicle-Stimulating Hormone (FSH) and Luteinizing Hormone (LH) are intimately involved in the control of the growth and reproductive activities of the gonadal tissues, which synthesize and secrete male and female sex hormones. The levels of circulating FSH and LH are controlled by these sex hormones through a negative feedback relationship. FSH is a glycoprotein secreted by the gonadotropic cells of the anterior pituitary. Gonadotropin-releasing hormone (GnRH), produced in the hypothalamus, controls the release of FSH from the anterior pituitary. Like other glycoproteins, such as LH, TSH, and HCG, FSH consists of subunits designated as alpha and beta. Hormones of this type have the same alpha subunits; therefore, the biological and immunological properties of each hormone are dependent on the unique beta subunit. In the female, FSH stimulates the growth and maturation of ovarian follicles by acting directly on the receptors located on the granulosa cells; follicular steroidogenesis is promoted and LH production is stimulated. The LH produced then binds to the theca cells and stimulates steroidogenesis. Increased intraovarian estradiol production occurs as follicular maturation advances, thereupon stimulating increased FSH receptor activity and FSH follicular binding. FSH, LH, and estradiol are therefore intimately related in supporting ovarian recruitment and maturation in women.

FSH levels are elevated after menopause, castration, and in premature ovarian failure. The levels of FSH may be normalized through the administration of estrogens, which demonstrate a negative feedback mechanism. Abnormal relationships between FSH and LH, and between FSH and estrogen have been linked to anorexia nervosa and polycystic ovarian disease. Although there are significant exceptions, ovarian failure is indicated when random FSH concentrations exceed 40 mIU/ml.

The growth of the seminiferous tubules and maintenance of spermatogenesis in men are regulated by FSH. However, androgens, unlike estrogens, do not lower FSH levels, therefore demonstrating a feedback relationship only with serum LH. For reasons not fully understood, azoospermic and oligospermic males usually have elevated FSH levels. Tumors of the testes generally depress serum FSH concentrations, but levels of LH are elevated, as determined by radioimmunoassay. It has been postulated that the apparent LH increase may be caused by cross-reactivity with hCG-like substances secreted by tumors of the testes.

High levels of FSH in men may be found in primary testicular failure and Klinefelter syndrome. Elevated concentrations are also present in cases of starvation, renal failure, hyperthyroidism.

PRINCIPLE:

The PSI FSH ELISA kit is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators, controls and patient samples are incubated in streptavidin coated microplate together with mixture of a biotinylated Anti- FSH monoclonal antibody and an HRP-conjugated monoclonal antibody, both recognize FSH from different sites. FSH present in calibrators/controls or samples is reacted with both antibodies to form a sandwich complex and simultaneously captured by streptavidin coated microwells via biotinylated Anti-FSH Mab. The strips are then washed and, buffered substrate/chromogen reagent (hydrogen peroxide and tetramethylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During

FSH ELISA KIT		IVD-REF: PS - FSH
		Ver. No: 02

the enzyme reaction, a blue color will develop and the intensity of the color is directly proportional to the amount of FSH present in the samples.

REAGENTS PROVIDED

Reagent	Quantity/48 Test KIT	Quantity/96 Test KIT	Quantity/192 Test KIT	Reconstitution
Microplate (96 breakable wells) coated with streptavidin	48 wells	96 wells	192 wells	Ready to Use
Calibrators 1-6 (0, 5, 10, 25,50, 100 mIU/mL) in buffer, containing preservative with traceability to reference material “ NIBSC/WHO, 3rd IS, HMG ”	6 Vial/0.5 ml	6Vial/0.5 ml	6 Vial/1.0 ml	Ready to Use
Control serum in buffer compatible with human serum matrix with preservatives	1 Vial/0.5 ml	1 Vial/0.5 ml	1 Vial/1.0 ml	Ready to Use
Enzyme Conjugate (red color)	1 Vial/6.0 ml	1 Vial/12.0 ml	2Vial/12 ml	Ready to use
Concentrated Wash solution	1 Vial/30ml	1 Vial/30 ml	1 Vial/50 ml	Dilute 20X with distilled water
Substrate-Chromogenic solution (Tetramethylbenzidine and hydrogen peroxide)	1 Vial/6.0ml	1 Vial/12ml	2 Vial/12ml	Ready to Use
Stop solution (H2So4 0.12M) Additional information about other additive is considered proprietary.	1 Vial/6.0 ml	1 Vial/6.0 ml	1 Vial/12.0 ml	Ready to Use

Materials required but not supplied with kit

The following materials are required but not provided in the kit:

- Distilled water with conductivity <math>< 1\mu\text{S}/\text{cm}</math>
- Precision micropipettes for delivery of,20, 50, 100, microliters. An 8-channel pipette or resenser pipette with disposable plastic tips for delivery of 50-100 μL is useful but not essential. Pipettes for dispensing milliliter volumes.
- Vortex mixer
- Automatic washer or any other device which can deliver 350 microliter wash solution for microplates
- Microplate reader with a wavelength of 450 nm and 620-650 nm, and an absorbance range of 0-3.0

STORAGE AND EXPIRATION DATING OF REAGENTS

- If kept at 2-8 °C, before opening or reconstitution, all kit components are stable until the expiry date as indicated on the label on the outside of the kit box.
- Calibrator concentrations are displayed on vial labels and may vary between lots.
- Do not freeze kits.
- Freshly prepared Working Wash solution should be used on the preparation day.
- Each lot of reagents and calibrators has been standardized to produce the correct reaction.

Do not interchange the reagents or calibrators between lots.

- Alterations in physical appearance of kit reagents may indicate instability or deterioration. Chromogenic solution TMB should be colorless. A blue color indicates that the reagent has been contaminated and should be discarded.

FSH ELISA KIT		IVD-REF: PS - FSH
		Ver. No: 02

SPECIMEN COLLECTION AND PREPARATION

- Human serum and EDTA plasma are samples of choice for PSI FSH kit.
- Do not use grossly haemolysed or turbid samples.
- If the measurement cannot be done within 24 hours after sample collection, samples should be kept at 2-8°C up to 48 hours. For longer storage serum or plasma must be stored at -20°C.
- Avoid repeated freeze and thaw cycles. Frozen specimens must be gently but thoroughly mixed after thawing.

PRECAUTIONS AND WARNINGS

Safety

- For Professional in vitro diagnostic Use Only.
- All the reagents within the kit are calibrated for the direct determination of FSH in human serum. The kit is not calibrated for the determination of FSH in urine, or Non-EDTA plasma or other specimens of human or animal origin.
- Before starting the assay, read the working manual completely and carefully. Use the valid version of the instruction for use provided with the kit. Be sure that everything is understood.
- Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., USA) publication No (CDC) 88—8_95 on laboratory safety procedures or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Reagents contain sodium azide (NaN₃) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. Follow local guidelines for disposal of all waste material.
- Avoid any skin contact with all reagents, stop solution contains strong acid. In case of contact, wash thoroughly with soap and water.

Caution

Material used in the preparation of human source reagents has been tested and found to be non-reactive for HIV 1 & 2 antibodies, HCV antibody and hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be done as if they were potentially infectious.

Preparation of reagents

Wash Buffer:

According to volume required, pour the desired volume of wash concentrate into a clean container and dilute 20 fold by adding adequate volume of distilled or deionized water to give a buffered Wash Solution. E.g. 10 mL wash concentrate to 190 mL distilled water. Wash solution should be prepared daily and discard unused working wash solution at the end of the day.

HANDELING NOTES

1. Do not use the kit or components beyond expiry date.
2. Do not mix up vials caps.
3. Bring all the reagents to room temperature prior to use.
4. Thoroughly mix all reagents, calibrators/controls and samples by gentle agitation or swirling.
5. Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
6. Use a clean container and high-quality distilled or deionized water to prepare the working wash solution.

FSH ELISA KIT		IVD-REF: PS - FSH
		Ver. No: 02

7. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
8. For the dispensing of Chromogenic solution, and Stop solution avoid pipettes with metal parts.
9. High precision pipettes or automated pipetting equipment will improve the precision.
10. Adhere the incubation times carefully. Tolerance limit of incubation time is $\pm 5\%$
11. Although PSI FSH EIA kit uses streptavidin coated plate technology and reaction does not proceed before adding conjugate solution, but to avoid drift, arising from low volume sample evaporation or incubation time difference between first and last well it is strongly recommended not to run more than six strips when pipetting manually (Time delay).
12. Perform all steps consecutively without any delay. Any delay in each step has deleterious effects on results.
13. Draw a calibration curve for each run, do not use stored curve from previous runs.
14. Incubation with Chromogenic solution, must be done in dark.
15. Only standard zero may be used to dilute any high serum samples. The use of any other reagent may lead to false results.
16. Opened reagents are stable maximally for 8 weeks if they are not contaminated, and stored in resealed original containers and handled as prescribed. Return to 2-8°C immediately after use

Preparation of the Sample:

Serum or plasma samples can be prepared according to common procedures used routinely in clinical laboratory practice. Samples can be stored at 2-8 °C if the assay is carried out within 24 hours, otherwise aliquots should be prepared and stored deep frozen (-20 °C). Frozen samples should be thawed and thoroughly mixed before assaying. Repeated freezing and thawing should be avoided.

PRECAUTION:

1. In one test-run do not combine strips, conjugate and standards from kits which have different lot numbers.
2. Solutions containing TMB and/or peroxide should not come into contact with metals or metal-ions, since this may give rise to unwanted color formation.
3. Allow that samples and all reagents come to ambient temperature before use.

Measurement Procedure

1. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminum pouch containing a desiccant and reseal carefully and store at 2-8°C.
2. Pipette 20 μL of each Calibrator (calibrators 0, 5, 10, 25, 50, 100, mIU/mL), controls and patient samples in duplicate into the appropriate wells.
3. Pipette 100 μL of enzyme conjugate into each well. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid.
4. Gently mix for 15 seconds
5. Cover the plate with an adhesive seal and incubate for 45 minutes at room temperature (20-27°C).
6. Wash the plate according to procedure below.
 - a. Carefully dispense 350 μL of working wash solution into each well.

- b. Aspirate the content of each well.
- c. Repeat 4 additional wash for a total of 5 wash cycles.
- d. After the end of last wash cycle invert the plate and tap it carefully against an absorbent paper.

Note: Automatic strip washer: Follow the manufacturer’s instructions for maintenance and wash the required number of wash cycles.

- 7. Add 100 µL of Chromogenic solution TMB to each well using the same pipetting technique as in procedural note 3. The Chromogenic solution TMB should be added to the wells as quickly as possible.
- 8. Incubate for 15 min at room temperature preferentially in the dark.
- 9. Add 50 µL of Stop Solution, mix and read the absorbance at 450 nm (with a reference filter 620, 630, or 650) in a microplate spectrophotometer within 10 min after addition of Stop Solution.

Important Note: always add stop solution in the same order as chromogenic solution TMB.


CALCULATION OF RESULTS

- 1. Construct a calibration curve by plotting the average O.D. (450 nm) for each calibrator on the vertical (Y) axis against the concentration on the horizontal (X) axis, and draw a best fit curve through the points on the graph.
- 2. To determine concentration of FSH in each sample, first locate the O.D. value on the Y-axis and extend a horizontal line to the calibration curve. At the point of intersection, draw a vertical line to the X-axis and read the corresponding concentration.
- 3. If a microplate spectrophotometer with built-in data calculation program is used, refer to the manual for the spectrophotometer and create a program using the concentration stated on the label of each of the FSH calibrators. For automatic calculation of FSH results it is recommended to use point-to-point method.

Typical Data

The following data are for illustration only and should never be used instead of the real-time calibration curve.

Row	OD	mIU/mL
1	0.028	0.0
2	0.165	5.0
3	0.302	10
4	0.674	25
5	1.312	50
6	2.387	100

Quality control

FSH control serums should be used for validation of the assay series. Ranges of expected results are indicated on the vial label. If assigned range are not obtained, a complete check of reagents and reader performance should be made and the assay shall be repeated. Each laboratory can also prepare its own serum pools at different levels, which can be used as internal controls in order to ensure the precision of the assay.

FSH ELISA KIT		IVD-REF: PS - FSH
		Ver. No: 02

Expected values

In a study conducted with some apparently healthy adults using the PSI FSH ELISA reagents, the following values are observed. Each laboratory should first investigate the transferability of the recommended expected value to its own patient population and if it does not match, determine its own reference interval.

sex	Population	2.5 th - 97.5 th Percentile [mIU/mL]
Female	Follicular phase	2.3-9.2
	Mid-cycle	7-18.6
	Luteal Phase	1.5-6.0
	Post-Menopausal	20-100
Men		1.0-12

PERFORMANCE CRITERIA

Lower limit of measurement:

The Limit of Blank and Limit of Detection were determined according to CLSI (Clinical and Laboratory Standards Institute) EP17-A2 guideline.

In brief Limit of Blank was determined as the 95th percentile value from $n = 60$ measurements on standard 0 over several runs. The Limit of blank equals to the concentration that 95% of repeated measurements on zero standard give value below it. In the next step, we determined standard deviation of repeated measurements on 3 samples with FSH contents less than 0.5 mIU/mL determined from a separate method study (ECLIA Method). By adding the LOB calculated at 1st step to 1.645 times the greatest standards deviation calculated in the 2nd step ($LOD = LOB + 1.645SD_L$) LOD calculated as 1.0 mIU/mL.

Precision:

Precision was determined using PSI FSH kit reagents and two pooled human sera and one Commercial control material according to EP05-A3 Guideline.

Three runs per day (morning, midday and afternoon) by 3 user in 10 working days ($10 \times 3 \times 3$). Results analyzed by fully nested ANOVA. Following results were obtained:

Sample Description	Mean (mIU/mL)	Repeatability		Within-Laboratory Precision	
		SD	%CV	SD	%CV
Patient Pool	6.01	0.28	4.66	0.37	6.16
L2-Biorad	19.62	0.91	4.64	0.94	4.79
Patient Pool	40.6	1.76	4.33	2.42	5.96

Specificity:

Specificity evaluation was performed according to guideline CLSI EP07-A2 to assess the cross-reactivity of the assay with other hypophyseal glycoprotein hormones. Samples containing cross-reactants was prepared by spiking natural human sera with material containing high levels of cross-reactants. Modification of serum matrix was no more than 10%. The percentage cross-reactivity was calculated for each sample using the equation below and normalized to the FSH content:

FSH ELISA KIT		IVD-REF: PS - FSH
		Ver. No: 02

$$\%cross - reactivity = \frac{\text{mean conc. of spiked sample} - \text{mean conc. of unspiked sample}}{\text{spiked concentration}} \times 100\%$$

The results are summarized in the following table:

interferent	Concentration	%interference
TSH	500 µIU/mL	0.05
LH	500 mIU/mL	0.16
hCG	2500 IU/L	0.08

Also interference evaluation from common sample abnormalities, hemolysis, Lipemia and icterus were performed according to aforementioned document and percentage interference was calculated for each sample using the equation below and normalized to the sample FSH content.

$$\%interference = \frac{\text{measured value} - \text{true value}}{\text{true value}} \times 100\%$$

No interference were seen up to following listed concentration:

Interferent	source	Final concentration of interferent in test serum
Hemoglobin	Fresh haemolysed blood	50 mg/mL
Conjugated Bilirubin	Commercially crystalline Unconjugated Bilirubin	20 mg/dL
Triglycerides	Commercially Intralipid 10%	1000 mg/dL

Linearity

Evaluation of Linearity was performed according to CLSI (Clinical and Laboratory Standards Institute) EP 06-A guideline

In brief a sample with concentration of 83 mIU/mL FSH as a high sample was diluted with another sample with 2 mIU/mL in different proportions and their FSH, measured in duplicate on each dilution. Measured values compared with expected values which estimated by dilution formula, then % bias and %recovery calculated for each concentration and following results are obtained:

No	Ratio	Expected (mIU/mL)	Rep 1 (mIU/mL)	Rep 2 (mIU/mL)	recovery%	%Bias
1	1	83	83	83	NA	NA
2	0.9	74.9	73.2	70.4	95.86%	-4.14%
3	0.8	66.8	64.5	68	99.18%	-0.82%
4	0.7	58.7	57.2	61.25	100.89%	0.89%
5	0.6	50.6	48.7	47.9	95.45%	-4.55%
6	0.5	42.5	40.7	39	93.76%	-6.24%
7	0.4	34.4	32.7	31.5	93.31%	-6.69%
8	0.3	26.3	26.5	27.4	102.47%	2.47%
9	0.2	18.2	17.7	18.6	99.73%	-0.27%
10	0.1	10.1	9.4	11.3	102.48%	2.48%
11	0	2.0	2	2	NA	NA

NA: Not applicable

FSH ELISA KIT		IVD-REF: PS - FSH
		Ver. No: 02

Trueness:

Recovery Evaluation (Proportional error)

Recovery was assessed by adding 10 mIU/mL of FSH to different samples with different FSH contents. Then % recovery and % bias was calculated for each sample using the equations below:

$$\%Recovery = \frac{Concentration\ Recovered}{Concentration\ Added} \times 100$$

$$\%Bias = \frac{\bar{X}\ measured - expected}{Expected} \times 100$$

The results are summarized in the following table:

Sam. No	Base Conc. mIU/mL	Added Analyte	Rep 1 mIU/mL	Rep 2 mIU/mL	Bias%	Recovery%
1	47.35	10 mIU/mL	56.7	59.9	1.66%	109.50%
2	34.1		45.2	43.7	0.79%	103.50%
3	12.25		22.3	23.3	2.47%	105.50%
4	1.865		11.8	12.7	3.24%	103.85%
5	59.4		69	70.7	0.65%	104.50%
6	6.9		16.7	17.7	1.78%	103.00%


Method comparison:

A comparison of the PSI FSH ELISA kit (y axis) using patient samples (n=84 range:2.5-88 mIU/mL) with ECLIA-Cobas E411-Roche diagnostic)X axis) gave the following results by ordinary linear regression:

PSI ELISA = 0.9877 ECLIA – 0.8522

r = 0.9982

r² = 0.9965

High dose hook effect


In PSI FSH EIA kit high dose hook effect has not been seen up to 2000 mIU/mL

FSH ELISA KIT		IVD-REF: PS - FSH
		Ver. No: 02

References:

1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Desirable Specifications for Total Error, Imprecision, and Bias, derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
5. Kasper D.L. et al. (2015). Harrison's principles of internal medicine.19th ed. Mc Graw Hill Education
6. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests.6th ed. Elsevier Inc.
7. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. Elsevier Inc.
8. Royal college of pathologists of Australasia. <https://www.datainnovations.com/allowable-total-error-table>
9. Scott M.G. et al. (1989). Hormonal evaluation of female infertility and reproductive disorders. Clin. Chem. 1989. 35:620-30

Procedure in-brief


FSH ELISA KIT

IVD-REF: PS - FSH

Ver. No: 02

SYMBOLS:

	BATCH NUMBER
	USE BY
	MANUFACTURER
	CONTAINS SUFFICIENT FOR <n> TESTS
	<i>IN VITRO</i> DIAGNOSTIC MEDICAL DEVICE
	TEMPERATURE LIMITATION
	CATALOGUE NUMBER
	CONSULT INSTRUCTIONS FOR USE