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$$\%Bias = \frac{\bar{X}\ measured - expected}{Expected} \times 100$$

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Sample	Original Concentration (ng/mL)	1.86 ng/mL added	
		Recovery%	Bias
1	1.48	116.41%	9.03%
2	2.54	99.64%	-0.15%
3	3.52	111.42%	3.87%
4	4.22	115.05%	4.51%

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
1. Catalona W.J. et al. (1995). Evaluation of percentage of free serum prostate-specific antigen to improve specificity of prostate cancer screening. JAMA. 274(15):1214-20.
2. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline. 2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
3. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
4. CLSI. Interference Testing in clinical chemistry; approved guideline. 2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
5. CLSI. User Protocol for Evaluation of qualitative test performance; Approved guideline 2<sup>nd</sup> ed. CLSI document EP 12-A2. Garrett P.E. Clinical and Laboratory Standard Institute, 2008.
6. Desirable Specifications for Total Error, Imprecision, and Bias, derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
7. Gatalon, W.J. et al. (1991). Measurement of prostate specific antigen in serum as a screening test for prostate cancer. The England Journal of Medicine Anthropological Research, 324(17), 1156-61
8. Kasper D.L. et al. (2015). Harrison's principles of internal medicine. 19th ed. (pp. 580-84). Mc Graw Hill Education
9. McPherson R.A. et al (2017). Henry's Clinical Diagnosis and management by Laboratory Methods. 23rd ed. (pp. 1442-43) Elsevier Inc.
10. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests. 6th ed. (pp. 378-81) Elsevier Inc.
11. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics. 6th ed. (pp. 472-75) Elsevier Inc.

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### INTENDED USE:

Pishgaman Sanjesh Isatis, Free PSA ELISA KIT is intended for the quantitative determination of Free form of Prostatic Specific Antigen in human serum and plasma.

### INTRODUCTION:

Prostate-specific antigen is a serine protease from kallikrein family which is synthesized uniquely in the epithelial cells of the prostate gland. Its synthesis is under the control of androgen receptor. Because of its high degree of tissue specificity, it is perhaps the most widely used tumor marker discovered thus far. Most synthesized PSA are transported into semen and is responsible for dissolving the seminal coagulum.

A small amount of PSA get into the bloodstream where it is rapidly inactivated by complexing with protease inhibitors such as alpha-1-antichymotrypsin (ACT),  $\alpha$ -2-macroglobulin (to which it binds covalently), and  $\alpha$ -protease inhibitor (API). These forms are collectively known as complexed PSA (cPSA). About 90 percent of cPSA exists in the form of complex with ACT and is measured in total PSA measurement kits, but complex with  $\alpha$ -2-macroglobulin cannot be measured by conventional total PSA kits. Approximately 10-30% of total PSA exist in the free form without proteolytic activity.

Various surveys have shown prostate cancer prevalence between 36-51% among men aged 70-79 years old. PSA is produced by both nonmalignant and malignant epithelial cells and, as such, is prostate-specific, but not prostate cancer-specific. It may also increase in benign prostatic hypertrophy and prostatitis. Several measures has been proposed to increase diagnostic accuracy of total PSA measurement such as age-based reference interval, free PSA to total PSA ratio (fPSA/tPSA), specific measurement of PSA-ACT complex, PSA velocity (change in PSA level over time), PSA density (PSA level relationship to prostate size). Measurement of fPSA/tPSA ratio is likely to be most informative when total PSA is between 3 or 4 and 10 ng/mL (the so-called “diagnostic gray zone”). The ratio of fPSA/tPSA is lower in prostate cancer patients than in men with BPH and those without prostatic disease. In patients receiving therapy, particularly hormone withdrawal therapy, the fPSA/tPSA ratio cannot be utilized to differentiate prostate hyperplasia from cancer of prostate. An equimolar tPSA determination such as Pishgaman Sanjesh PSA ELISA KIT is the prerequisite for reliable ratios.

### PRINCIPLES OF THE METHOD:

The PSI fPSA ELISA Kit is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators, controls and patient samples are incubated together with a biotinylated Anti-fPSA monoclonal antibody and an HRP-conjugated monoclonal antibody; both recognize fPSA from different sites, in streptavidin coated microstrips. Free PSA present in calibrators/controls or samples is reacted with both antibodies to form a sandwich complex and adsorbed to the streptavidin-coated microstrips via the biotinylated Anti-fPSA Mab. The strip are then washed and, buffered substrate/chromogen reagent (hydrogen peroxide and tetramethylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction, a blue color will develop and the intensity of the color is directly proportional to the amount of fPSA present in the samples.

<b>Free PSA ELISA KIT</b>		IVD-REF: PS – Free PSA
		Ver. No: 02

#### REAGENTS PROVIDED:

Reagent	Quantity/48 Test KIT	Quantity/96 Test KIT	Reconstitution
Microplate (96 breakable wells) coated with Streptavidin	48 wells	96 wells	Ready to Use
Calibrators 1-6 (0, 0.25, 1.0, 2.5, 5.0 and 10 ng/mL) in buffer, containing preservative	6 Vials/1.0 ml	6 Vials/1.0 ml	Ready to Use
Controls N = 2 in buffer compatible with human serum matrix with preservatives	2Vials/1.0 ml	2Vials/1.0 ml	Ready to Use
Enzyme Conjugate (red color)	1 Vial/6.0 ml	1 Vial/12 ml	Ready to use
Concentrated Wash solution	1 Vial/30ml	1 Vial/30 ml	Dilute 20X with distilled water
Substrate-Chromogenic solution (Tetramethylbenzidine and hydrogen peroxide)	1 Vial/6.0ml	1 Vial/12ml	Ready to Use
Stop solution (HCl 2N) Additional information about other additive is considered proprietary.	1 Vial/3.0 ml	1 Vial/6.0 ml	Ready to Use

#### Materials required but not supplied with kit

The following materials are required but not provided in the kit:

- Distilled water with conductivity < 1 $\mu$ S
- Precision micropipettes for delivery of 50, 100 microliters. An 8-channel pipette or resenser pipette with disposable plastic tips for delivery of 50-100  $\mu$ L is useful but not essential.
- Pipettes for dispensing milliliter volumes.
- Vortex mixer
- Automatic washer or any other device which can deliver 350 microliter wash solution for microplates
- Microplate reader with a filter set at 450 nm and an upper OD limit of 3.0 or greater\* (see measurement procedure step 9).

#### STORAGE AND EXPIRATION DATING OF REAGENTS

- If kept at 2-8°C, before opening or reconstitution, all kit components are stable until the expiry date as indicated on the label on the outside of the kit box.
- Freshly prepared Working Wash solution should be used on the preparation day.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration of reagents. The substrate solution (TMB) is sensitive to light and should remain colorless if properly stored. Instability or contamination may be indicated by the development of a blue color, in which case it should not be used.

#### SPECIMEN COLLECTION AND PREPARATION

- Human serum and EDTA plasma are samples of choice for PSI fPSA kit.

<b>Free PSA ELISA KIT</b>		IVD-REF: PS – Free PSA
		Ver. No: 02

- If the measurement cannot be done within 3 hours after sample collection, samples should be kept at 2-8 °C up to 24 hours. For longer storage, serum or plasma must be stored at -20°C.
- Avoid repeated freeze and thaw cycles. Frozen specimens must be gently but thoroughly mixed after thawing.

#### PRECAUTIONS AND WARNINGS

- For Professional in vitro diagnostic Use Only
- All the reagents within the kit are calibrated for the direct determination of fPSA in human serum. The kit is not calibrated for the determination of fPSA in urine, semen, Non-EDTA plasma or other biologic specimens of human or animal origin.
- Before starting the assay, read this instruction completely and carefully. Use the valid version of the instruction provided with the kit. Be sure that everything is understood
- Please refer to the Guidelines for safe Work practice in Human and animal medical diagnostic laboratories recommendations of a CDC-convened, Biosafety Blue Ribbon Panel. MMWR Supplement / Vol. 61, January 6, 2012” or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Reagents contain sodium azide (NaN<sub>3</sub>) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. Follow local guidelines for disposal of all waste material.
- Avoid any skin contact with all reagents; stop solution contains strong acid. In case of contact, wash thoroughly with soap and water.

#### Caution

- Material used in the preparation of human source reagent has been tested and found to be non-reactive for HIV 1 & 2 antibodies, HCV antibody and hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

#### Preparation of reagents:

##### Wash Buffer:

According to volume required, pour the desired volume of concentrated wash solution into a clean container and dilute 1:20 by adding adequate volume of distilled or deionized water to give a buffered Wash Solution. E.g. 10 mL wash concentrate to 190 mL distilled water. Wash solution should be prepared daily and discard unused working wash solution at the end of the day.

#### HANDELING NOTES:

- Do not use the kit or components beyond expiry date.
- Do not mix up vials caps.
- Do not mix reagents from different kit lots.
- All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use.
- A calibrator curve must be established for every run.
- Calibrators, controls and specimen samples should be assayed in duplicate. Vertical alignment is recommended.

<b>Free PSA ELISA KIT</b>		IVD-REF: PS – Free PSA
		Ver. No: 02

- Use a clean container and high quality distilled or deionized water to prepare the working wash solution.
- In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
- For the dispensing of Chromogenic and Stop solution, avoid pipettes with metal parts.
- High precision pipettes or automated pipetting equipment will improve the precision.
- Adhere the incubation times carefully. Tolerance limit of incubation time is  $\pm 5\%$
- Although, fPSI PSA EIA kit uses streptavidin coated plate technology and reaction does not proceed before adding conjugate solution, but to avoid drift, arising from low volume sample evaporation or incubation time difference between first and last well, it is strongly recommended not to run more than six strips when pipetting manually (Time delay).
- Perform all steps consecutively without any delay. Any delay in each step has deleterious effects on measurement results.
- Draw a calibration curve for each run, do not use stored curve from previous runs.
- Incubation with Chromogenic solution must be done in dark.
- Opened reagents are stable maximally for 8 weeks if they are not contaminated, and stored in resealed original containers and handled as prescribed. Return to 2-8°C immediately after use
- Kit reagents must be regarded as hazardous waste and disposed of according to national regulations

## Measurement Procedure:

1. Remove the required number of well strips. Reseal the bag and return any unused strips to the refrigerator.
2. Pipette 50  $\mu\text{L}$  of each Calibrator (calibrators 0, 0.25, 1.0, 2.5, 5.0 and 10 ng/mL), controls and patient samples in duplicate into the appropriate wells.
3. Pipette 100  $\mu\text{L}$  of enzyme conjugate solution into each well. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid.
4. Gently mix for 15 seconds
5. Cover the plate with an adhesive seal and incubate for 45-minutes at room temperature (20-27°C).
6. Wash the plate according to procedure below.
  - i. Carefully dispense 350  $\mu\text{L}$  of working wash solution into each well.
  - ii. Aspirate the content of each well.
  - iii. Repeat 4 additional wash for a total of 5 wash cycles.
  - iv. After the end of last wash cycle invert the plate and tap it carefully against an absorbent paper.

Note: Automatic strip washer: Follow the manufacturer's instructions for maintenance and wash the required number of wash cycles after second incubation step

7. Add 100  $\mu\text{L}$  of Chromogenic solution TMB to each well using the same pipetting technique as in item 3. The Chromogenic solution TMB should be added to the wells as quickly as possible and the time between addition to the first and last well should not exceed 5 min.

<b>Free PSA ELISA KIT</b>		IVD-REF: PS – Free PSA
		Ver. No: 02

- Incubate for 15 min at room temperature preferentially in dark. Avoid exposure to direct sunlight.
- Add 50  $\mu$ L of Stop Solution, mix and read the absorbance at 450 nm in a microplate spectrophotometer within 5 min after addition of Stop Solution.

***Important Note:*** always add stop solution in the same order as chromogenic solution TMB.


## CALCULATION OF RESULTS

- Construct a calibration curve by plotting the average O.D. (450 nm) for each calibrator on the vertical (Y) axis against the concentration on the horizontal (X) axis, and draw a best fit curve through the points on the graph.
- To determine concentration of fPSA in each sample, first locate the O.D. value on the Y-axis and extend a horizontal line to the calibration curve. At the point of intersection, draw a vertical line to the X-axis and read the corresponding concentration.
- If a microplate spectrophotometer with built-in data calculation program is used, refer to the manual for the spectrophotometer and create a program using the concentration stated on the label of each of the PSA calibrators. For automatic calculation of fPSA results, it is recommended to use point-to-point method.

### Typical Data:

The following data are for illustration only and should never be used instead of the real-time calibration curve.

calibrator	OD	ng/mL
1	<b>0.03</b>	<b>0</b>
2	<b>0.11</b>	<b>0.25</b>
3	<b>0.35</b>	<b>1.0</b>
4	<b>0.72</b>	<b>2.5</b>
5	<b>1.41</b>	<b>5.0</b>
6	<b>2.83</b>	<b>10</b>

### Quality control:

Free PSA control serums should be used for validation of the assay series. Ranges of expected results are indicated on the vial label. If assigned range are not obtained, a complete check of reagents and reader performance should be made and the assay shall be repeated. Each laboratory can also prepare its on serum pools at different levels, which can be used as internal controls in order to ensure the precision of the assay.

### Reference values

Free PSA measurement must be used with an equimolar test such as Pishgaman Sanjesh PSA ELISA (PSPS-01) for Total PSA in order to generate reliable fPSA/tPSA ratio. Serum specimens

<b>Free PSA ELISA KIT</b>		IVD-REF: PS – Free PSA
		Ver. No: 02

from 48 men diagnosed as benign prostatic hyperplasia (BPH) and 68 men diagnosed as prostate cancer (PCa) by DRE, were analysed by PSI fPSA and Total PSA ELISA kits and the following results were obtained:

Diagnosis	fPSA/tPSA		
	Median	Range	Mean
<b>BPH</b>	<b>0.19</b>	<b>0.04-0.4</b>	<b>0.2</b>
<b>PCa</b>	<b>0.09</b>	<b>0.02-0.5</b>	<b>0.11</b>

The optimum cut-off in different clinical disciplines depends upon the clinical intended use of results, and desired diagnostic sensitivity and specificity. Diagnostic sensitivity and specificity move inversely, as sensitivity increases, specificity decreases and vice versa. Low specificity correlates with high rate of false negatives and Low sensitivity with high rate of false positives for PCa. Sensitivities (%correctly diagnosed PCa) and specificities (%correctly diagnosed BPH) for different fPSA/tPSA ratio cut-offs are shown below:

**Note: The laboratory must provide information about diagnostic sensitivity and specificity and false results rate for customers as requested.**

fPSA/tPSA	PCa(TP/FP)	BPH(TN/FN)	Diagnostic Sensitivity(CI)	Diagnostic Specificity(CI)
<b>0.09</b>	<b>32 /3</b>	<b>45/36</b>	<b>47.06 % ( 34.83% to 59.55%)</b>	<b>93.75% (82.80% to 98.69%)</b>
<b>0.13</b>	<b>50/9</b>	<b>39/18</b>	<b>73.53% (61.43% to 83.50%)</b>	<b>81.25% (67.37% to 91.05%)</b>
<b>0.23</b>	<b>66/33</b>	<b>15/2</b>	<b>97.06 % ( 89.78% to 99.64%)</b>	<b>31.25 % ( 18.66% to 46.25%)</b>

## PERFORMANCE CRITERIA:

### 1. Detection capabilities:

The Limit of Blank and Limit of Detection were determined according to CLSI (Clinical and Laboratory Standards Institute) EP17-A2 guideline.

In brief Limit of Blank was determined as the 95th percentile value from n = 60 measurements on standard 0 over several runs. The Limit of Blank equals to the concentration that 95% of repeated measurements on zero standard give value below it. In the next step, we determined standard deviation of repeated measurements on 3 samples with fPSA contents less than 0.1 ng/mL determined from a separate method study (ECLIA Method). By adding, the LOB calculated at 1<sup>st</sup> step to 1.645 times the greatest standards deviation calculated in the 2<sup>nd</sup> step ( $LOD = LOB + 1.645SD_L$ ) LOD calculated as 0.05 ng/mL.

### 2. Precision:

Precision was determined using PSI fPSA reagents and three human pooled sera in different levels of measuring range (0.05-10 ng/mL) according to EP05-A3 Guideline.

Three runs per day (morning, midday and afternoon) by 3 user in 10 working days(10 × 3 × 3) performed on these samples Results analyzed by fully nested ANOVA using a spreadsheet. The results are summarized in table below:

<b>Free PSA ELISA KIT</b>		IVD-REF: PS – Free PSA
		Ver. No: 02

Sample Description	Mean (ng/mL)	Repeatability		Within-Laboratory Precision	
		SD	%CV	SD	%CV
Patient Pool	<b>0.33</b>	<b>0.017</b>	<b>5.15</b>	<b>0.036</b>	<b>10.91</b>
Patient Pool	<b>3.04</b>	<b>0.125</b>	<b>4.11</b>	<b>0.428</b>	<b>14.08</b>
Patient Pool	<b>7.20</b>	<b>0.481</b>	<b>6.67</b>	<b>0.581</b>	<b>8.07</b>

### 3. specificity

Interference evaluation from common sample abnormalities, hemolysis, Lipemia and icterus were determined according to CLSI (Clinical and Laboratory Standards Institute) EP07-A2 guideline and %interference calculated using following formula:

$$\%interference = \frac{\text{measured value} - \text{true value}}{\text{true value}} \times 100\%$$

No interference were seen up to following listed concentration:

Compound	Concentration with no Significant interference
Lipemia (Intralipid)	1000 mg/dL
Bilirubin(unconjugated)	20 mg/dL
hemoglobin	500 mg/dL

### 4. Linearity

Evaluation of Linearity was performed according to CLSI (Clinical and Laboratory Standards Institute) EP 06-A guideline. In brief a sample with high fPSA (9.65 ng/mL) was diluted with a low fPSA sample (0.34 ng/mL) in different proportions and their fPSA, measured in duplicate on each dilution. Measured values compared with expected values which estimated by dilution formula, then % bias and %recovery calculated for each concentration and following results are obtained:

No	Ratio	Expected (ng/mL)	Rep 1 (ng/mL)	Rep 2 (ng/mL)	recovery%	%Bias
1	1	9.65	9.33	9.97	NA	NA
2	0.9	8.72	8.39	8.53	97%	-2.9%
3	0.8	7.79	7.48	7.97	99%	-0.8%
4	0.7	6.86	6.88	6.97	101%	1.0%
5	0.6	5.93	6.20	7.06	112%	11.9%
6	0.5	4.99	5.24	5.71	110%	9.6%
7	0.4	4.06	4.59	4.46	111%	11.5%
8	0.3	3.13	3.32	3.33	106%	6.1%
9	0.2	2.20	2.34	2.42	108%	8.3%
10	0.1	1.27	1.29	1.32	103%	2.8%
11	0	0.34	0.34	0.34	NA	NA

<b>Free PSA ELISA KIT</b>		IVD-REF: PS – Free PSA
		Ver. No: 02

**5. Trueness:**

**5.1. Recovery Evaluation (Proportional error)**

Recovery was assessed by adding a given amount (1.86 ng/mL) of fPSA to 4 samples with different fPSA contents. Then % recovery and % bias was calculated for each sample using the equations below:

$$Recovery = 100 \times \frac{\text{concentration recovered}}{\text{concentration added}}$$

$$\%Bias = \frac{\bar{X} \text{ measured} - \text{expected}}{\text{Expected}} \times 100$$

The results are summarized in the following table:

Sample	Original Concentration (ng/mL)	1.86 ng/mL added	
		Recovery %	Bias
1	1.48	116.41%	9.03%
2	2.54	99.64%	-0.15%
3	3.52	111.42%	3.87%
4	4.22	115.05%	4.51%


**5.2. Trueness (Method comparison)**

A comparison of equimolar methods (PSI PSA /Free PSA ELISA Kits) and (CanAg PSA ref no: 340-10/Free PSA ref no: 350-10) were performed using 47 patient samples and fPSA/tPSA ratio obtained by each pair of methods compared by ordinary linear regression. Following results were obtained:

$$PSI = 1.0681CanAg + 0.0106$$

$$r = 0.991457$$

$$r^2 = 0.982987$$

<b>Free PSA ELISA KIT</b>		IVD-REF: PS – Free PSA
		Ver. No: 02

## 6. High dose hook effect

In PSI Free PSA EIA kit high dose hook effect has not been seen up to 150 ng/mL

## REFERENCES:

1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. CLSI. User Protocol for Evaluation of qualitative test performance; Approved guideline 2<sup>nd</sup> ed. CLSI document EP 12-A2. Garrett P.E. Clinical and Laboratory Standard Institute, 2008.
5. Desirable Specifications for Total Error, Imprecision, and Bias, derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
6. Gatalon, W.J. et al. (1991). Measurement of prostate specific antigen in serum as a screening test for prostate cancer. The England Journal of Medicine Anthropological Research, 324(17), 1156-61
7. Kasper D.L. et al. (2015). Harrison’s principles of internal medicine.19th ed.(pp. 580-84). Mc Graw Hill Education
8. McPherson R.A. et al (2017). Henry’s Clinical Diagnosis and management by Laboratory Methods.23rd ed. (pp.1442-43) Elsevier Inc.
9. Pagana K.D. et al. (2018). Mosby’s Manual of diagnostic and Laboratory tests.6th ed. (pp. 378-81) Elsevier Inc.
10. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. (pp.472-75) Elsevier Inc.

<b>Free PSA ELISA KIT</b>		IVD-REF: PS – Free PSA
		Ver. No: 02

**Procedure in-brief**