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$$\%Bias = \frac{\bar{X}\ measured - expected}{Expected} \times 100$$

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Sample	Original Concentration (mIU/mL)	5 mIU/mL added		10 mIU/mL added	
		Recovery %	Bias	Recovery %	Bias
1	5.0	97.6	-1.52%	96.0	-1.90%
2	17.4	104.0	0.89%	101.0	0.36%
3	45.1	106.0	0.60%	108.5	1.54%

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
1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline. 2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline. 2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Desirable Specifications for Total Error, Imprecision, and Bias, derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
5. McPherson R.A. et al. (2017). HENRY'S Clinical diagnosis and management by Laboratory Methods. 23rd ed. (pp. 406-408). Elsevier Inc.
6. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests. 6th ed. (pp. 311-313). Elsevier Inc.
7. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics. 6th ed. (pp. 633-634). Elsevier Inc.

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INTENDED USE:

The Quantitative Determination of Luteinizing hormone (LH) in Human Serum.

CLINICAL SIGNIFICANCE:

Luteinizing hormone (LH) is a glycoprotein hormone secreted in both men and women from the anterior pituitary gland in response to luteinizing hormone-releasing hormone (LH-RH or Gn-RH), which is released by the hypothalamus. LH, also called interstitial cell-stimulating hormone (ICSH) in men, is a glycoprotein with a molecular weight of approximately 30.000 Daltons. LH Like FSH, hCG and TSH is a heterodimer consist of two non-covalently associated polypeptide chains α and β . Alpha chain is identical in all of these hormone, but β chain is unique in each of them and is responsible for their specific physiological activities and immunological characteristics.

LH and FSH released are secreted in pulsatile and pass through bloodstream to the ovaries, where they stimulate the growth and maturation of follicle and hence the biosynthesis of estrogens and progesterone. The highest level of LH occurs at mid-cycle peak that induce ovulation and formation of corpus luteum, which mainly produces progesterone. Primary function of LH in men is stimulation of the interstitial cells (Leydig cells) to produce testosterone.

LH measurement is used to evaluate dysfunction of the hypothalamus-pituitary-gonads system.

It is always requested in conjunction with FSH if any of the congenital disease with chromosomal abnormalities (e. g. Klinefelter syndrome), polycystic ovaries (PCO) indicated. It is also used to find the cause of amenorrhea, menopausal syndrome, and suspected Leydig cell deficiency.

PRINCIPLE:

The PSI LH EIA Kit is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators, controls and patient samples are incubated together with a biotinylated Anti-LH monoclonal antibody and an HRP-conjugated monoclonal antibody; both recognize LH from different sites, in streptavidin coated microstrips. LH present in calibrators/controls or samples is reacted with both antibodies to form a sandwich complex and also adsorbed to the streptavidin coated microstrips via the biotinylated Anti-LH Mab. The strips are then washed and, buffered substrate/chromogen reagent (hydrogen peroxide and tetramethylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction, a blue color will develop and the intensity of the color is directly proportional to the amount of LH present in the samples.

Reagent and material supplied in the kit:

Reagent	Quantity/96 Test KIT	Reconstitution
Microplate (96 breakable wells) coated with streptavidin	96 wells	Ready to Use
Calibrators 1-6 (0,5,10, 25,50and 100 mIU/ml LH) in buffer, containing preservative with traceability to reference material "WHO 2nd International Standard 552/80 "	6 Vial/0.5ml	Ready to Use
Control N = 1 in buffer compatible with human serum matrix with preservatives	1 Vial/0.5 ml	Ready to Use
Enzyme Conjugate (red color)	1 Vial/12 ml	Ready to use
Concentrated Wash solution	1 Vial/30 ml	Dilute 20X with distilled water
Substrate-Chromogenic solution (Tetra-methylbenzedrine and hydrogen peroxide)	1 Vial/12ml	Ready to Use
Stop solution (H2So4 0.12M) Additional information about other additive is considered proprietary.	1 Vial/6.0 ml	Ready to Use

LH ELISA KIT		IVD-REF: PS - LH
		Ver. No: 02

Materials required but not supplied with kit

The following materials are required but not provided in the kit:

- Distilled water with conductivity < 1µS/cm
- Precise micropipettes for delivery of 20, 50, 100, microliters. An 8-channel pipette or respenser pipette with disposable plastic tips for delivery of 50-200 µL is useful but not essential.
- Pipettes for dispensing milliliter volumes.
- Vortex mixer
- Automatic washer or any other device which can deliver 350 microliter wash solution for microplates
- Microplate reader with a wavelength of 450 nm and 620-650 nm, and an absorbance range of 0-3.0

STORAGE AND EXPIRATION DATING OF REAGENTS

- If kept at 2-8 °C, before opening or reconstitution, all kit components are stable until the expiry date as indicated on their label.
- Calibrators' concentration is displayed on the vial labels and may vary between lots.
- Do not freeze kit components.
- Freshly prepared Working Wash solution should be used on the preparation day.
- Each lot of reagents and calibrators has been standardized to produce the correct reaction. Do not interchange the reagents or calibrators between lots.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration. Chromogenic solution TMB should be colorless. Blue color indicates that the reagent has been contaminated and should be discarded.

SPECIMEN COLLECTION AND PREPARATION

- Human serum and EDTA plasma are samples of choice for PSI LH kit.
- Do not use grossly haemolysed or turbid samples.
- If the measurement cannot be done within 24 hours after sample collection, samples should be kept at 2-8°C up to 48 hours. For longer storage, serum or plasma must be stored at -20°C.
- Avoid repeated freeze and thaw cycles. Frozen specimens must be gently but thoroughly mixed after thawing.

PRECAUTIONS AND WARNINGS

Safety

- For Professional in vitro diagnostic use only.
- All the reagents within the kit are calibrated for the direct determination of LH in human serum. The kit is not calibrated for the determination of LH in urine, Non-EDTA plasma or other specimens of human or animal origin.

LH ELISA KIT		IVD-REF: PS - LH
		Ver. No: 02

- Before starting the assay, read the working manual completely and carefully. Use the valid version of the instruction provided with the kit. Be sure that everything is understood.
- Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., USA) publication No (CDC) 88—8_95 on laboratory safety procedures or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Reagents contain sodium azide (NaN₃) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. Follow local guidelines for disposal of all waste material.
- Avoid any skin contact with all reagents. Stop solution contains strong acid in case of contact, wash thoroughly with soap and water.

Caution

Material used in the preparation of human source reagents has been tested and found to be non-reactive for HIV 1 & 2 antibodies, HCV antibody and hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they are potentially infectious.

Preparation of reagents

Wash Buffer:

According to the required volume, pour the desired volume of wash concentrate into a clean container and dilute 20 fold by adding adequate volume of distilled or deionized water to gain buffered Wash Solution. E.g., 10 mL wash concentrate to 190 mL distilled water. Wash solution should be prepared daily and discard unused working wash solution at the end of the day.

HANDELING NOTES

- Do not use the kit or components beyond the expiry date.
- Do not mix up vial caps.
- Bring all the reagents to room temperature prior to use.
- Thoroughly mix all reagents, calibrators/controls and samples by gentle agitation or swirling.
- Assay calibrators, controls and samples in duplicate. Vertical alignment is recommended.
- Use a clean container and high quality distilled or deionized water to prepare the working wash solution.
- In order to avoid cross-contamination, use a clean disposable pipette tip for addition of each reagent and sample.
- For the dispensing of Chromogenic solution and Stop solution, avoid pipettes with metal parts.
- High precision pipettes or automated pipetting equipment will improve the precision.
- Adhere the incubation times carefully. Tolerance limit of incubation time is $\pm 5\%$

LH ELISA KIT		IVD-REF: PS - LH
		Ver. No: 02

- Although, PSI LH EIA kit uses streptavidin coated plate technology and reaction does not proceed before adding conjugate solution, but to avoid drift, arising from low volume sample evaporation or incubation time difference between first and last well, it is strongly recommended not to run more than six strips when pipetting manually (Time delay).
- Perform all steps consecutively without any delay. Any delay in each step has deleterious effects on final results.
- Draw a calibration curve for each run, do not use stored curve from previous runs.
- Incubation with Chromogenic solution must be done in the dark.
- Only standard zero may be used to dilute any high serum samples. The use of any other reagent may lead to false results.
- Opened reagents are stable maximally for 8 weeks if they are not contaminated, and stored in resealed original containers and handled as prescribed. Return to 2-8°C immediately after use.

Measuring Procedure

1. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminum pouch containing a desiccant, reseal carefully and store at 2-8°C.
 2. Pipette 20 µL of each Calibrator (calibrators 0, 5, 10, 25, 50, 100 mIU/mL), controls and patient samples in duplicate into the appropriate wells.
 3. Pipette 100 µL of enzyme conjugate into each well. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid.
 4. Gently mix for 15 seconds
 5. Cover the plate with an adhesive seal and incubate for 45 minutes at room temperature (20-27°C).
 6. Wash the plate according to procedure below.
 - i. Carefully dispense 350 µL of working wash solution into each well.
 - ii. Aspirate the content of each well.
 - iii. Repeat 4 additional wash for a total of 5 wash cycles.
 - iv. After the last wash cycle invert the plate and tap it carefully against an absorbent paper.
- Note: Automatic strip washer: Follow the manufacturer's instructions for maintenance and wash the required number of wash cycles.
7. Add 100 µL of Chromogenic solution TMB to each well using the same pipetting technique as in procedural note 3. The Chromogenic solution TMB should be added to the wells as quickly as possible.
 8. Incubate for 15 min at room temperature preferentially in the dark.
 9. Add 50 µL of Stop Solution, mix and read the absorbance at 450 nm (with a reference filter 620, 630, or 650 nm) in a microplate spectrophotometer within 10 min after addition of Stop Solution.

LH ELISA KIT		IVD-REF: PS - LH
		Ver. No: 02

Important Note: always add stop solution in the same order as chromogenic solution TMB.


CALCULATION OF RESULTS

1. Construct a calibration curve by plotting the average O.D. (450 nm) for each calibrator on the vertical (Y) axis against the concentration on the horizontal (X) axis, and draw the best fit curve through the points on the graph.
2. To determine concentration of LH in each sample, first locate the O.D. value on the Y-axis and extend a horizontal line to the calibration curve. At the point of intersection, draw a vertical line to the X-axis and read the corresponding concentration.
3. If a microplate spectrophotometer with built-in data calculation program is used, refer to the manual for the spectrophotometer and create a program using the concentration stated on the label of each calibrator. For automatic calculation of LH results, it is recommended to use point-to-point method.

Typical Data

The following data are for illustration only and should never be used instead of the real-time calibration curve.

OD	Standards mIU/ml
0.02	0
0.28	5
0.53	10
1.11	25
2.01	50
2.74	100

Quality control

LH control serum should be used for validation of the assay series. Range of expected results are indicated on the vial label. If assigned range is not obtained, a complete check of reagents and reader performance should be made and the assay shall be repeated. Each laboratory can also prepare its own serum pools at different levels, which can be used as internal controls in order to ensure the precision of the assay.

Expected values

A study with PSI LH EIA kit was performed using samples from some apparently healthy people and following reference interval corresponding to 2.5th and 97.5th percentiles of results was extrapolated. Each laboratory should first investigate the transferability of the recommended expected value to its own patient population and if it does not match, determine its own reference interval:

Phase of menstrual cycle	2.5 th -97.5 percentile (mIU/mL)
Follicular phase	1.8 - 10
Mid-cycle	7.0-49
Luteal phase	0.8-12
Post-Menopause	7.5-42 mIU/ml

LH ELISA KIT		IVD-REF: PS - LH
		Ver. No: 02

Male	1.0-7.0mIU/ml
------	---------------

PERFORMANCE CRITERIA

1. Lower limit of measurement:

The Limit of Blank and Limit of Detection were determined according to CLSI (Clinical and Laboratory Standards Institute) EP17-A2 guideline.

In brief, Limit of Blank was determined as the 95th percentile value from n = 60 measurements on standard zero over several runs. The Limit of blank equals to the concentration that 95% of the repeated measurements on zero standard give value below it. In the next step, we determined standard deviation of the repeated measurements on three samples with LH contents less than 1.0 mIU/mL determined from a separate method study (ECLIA Method). By adding, the LOB calculated at first step to 1.645 times the greatest standards deviation calculated in the second step ($LOD = LOB + 1.645SD_L$) LOD calculated as 1.0 mIU/mL

Specificity:

Specificity evaluation was performed according to guideline CLSI EP07-A2 to assess the cross-reactivity of the assay with other glycoprotein hormones. Samples containing cross-reactants was prepared by spiking natural human sera with material containing high levels of cross-reactants. The modification of serum matrix was no more than 10%. The percentage cross-reactivity was calculated for each sample using the equation below and normalized to the LH content:

$$= \frac{\%cross - reactivity}{\frac{mean\ conc.\ of\ spiked\ sample - mean\ conc.\ of\ unspiked\ sample}{spiked\ concentration}} \times 100\%$$

The results are summarized in the following table:

Tested hormones	Cross Reactivity%	Concentration
FSH	0.05	500 mIU/ml
hCG	0.37	25000 IU/L
TSH	0.18	500 µIU/ml

Also interference evaluation from common sample abnormalities, hemolysis, Lipemia and icterus were performed according to aforementioned document and No interference were seen up to following listed concentration:

Interferent	source	Final concentration of interferent in test serum
Hemoglobin	Fresh haemolysed blood	50 mg/mL
Conjugated Bilirubin	Commercially crystalline Unconjugated Bilirubin	20 mg/dL
Triglycerides	Commercially Intralipid 10%	1000 mg/dL

Precision:

Precision was determined using PSI LH EIA kit reagents and three pooled human sera according to EP05-A3 Guideline.

LH ELISA KIT		IVD-REF: PS - LH
		Ver. No: 02

Three runs per day (morning, midday and afternoon) by 3 user in 10 working days (10 × 3 × 3). Results analyzed by fully nested ANOVA. Following results were obtained:

Sample Description	Mean (mIU/mL)	Repeatability		Within-Laboratory Precision	
		SD	%CV	SD	%CV
Patient Pool	3.1	0.26	8.39	0.33	10.65
Patient Pool	43	2.8	6.51	1.0	2.33
Patient Pool	72	3.7	5.14	3.1	4.31

Linearity:

Evaluation of Linearity was performed according to CLSI (Clinical and Laboratory Standards Institute) EP 06-A guideline

In brief, a sample with concentration of 93 mIU/mL LH as a high sample was diluted with another sample with 1.3 mIU/mL in different proportions and their LH was measured in duplicate on each dilution. Measured values were compared with expected values, which were estimated by dilution formula, then percentage bias and percentage recovery were calculated for each concentration. The following results are obtained:

No	Ratio	Expected (mIU/mL)	Rep 1 (mIU/mL)	Rep 2 (mIU/mL)	recovery%	%Bias
1	1	93	91	95	NA	NA
2	0.9	83.83	81.5	79	95.73%	-4.27%
3	0.8	74.66	72	79	101.13%	1.13%
4	0.7	65.49	63	60	93.91%	-6.09%
5	0.6	56.32	52	59	98.54%	-1.46%
6	0.5	47.15	51	45	101.80%	1.80%
7	0.4	37.98	36.2	39.7	99.92%	-0.08%
8	0.2	19.64	21	18.5	100.56%	0.56%
9	0.1	10.47	9.3	9.7	90.74%	-9.26%
10	0	1.3	1.2	1.4	NA	NA

NA: Not applicable

Method comparison:


A comparison of the PSI LH EIA kit (y axis) using patient samples (n=100 range:2.0-99.0 mIU/mL) with ECLIA-Cobas E411-Roche diagnostic (X axis) gave the following results by ordinary Liner regression:

The regression equation is

$$PSI = - 0.0208 + 0.9905 ECLIA$$

$$r = 0.9973$$

$$r^2 = 0.9945$$

High dose Hook effect:



High dose Hook effect has not seen up to 2000 mIU/mL.

References:

1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline. 2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline. 2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Desirable Specifications for Total Error, Imprecision, and Bias, derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
5. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests. 6th ed. Elsevier Inc.
6. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics. 6th ed. Elsevier Inc.
7. Scott M.G. et al. (1989). Hormonal evaluation of female infertility and reproductive disorders. Clin. Chem. 1989. 35:620-30

LH ELISA KIT		IVD-REF: PS - LH
		Ver. No: 02

SYMBOLS:

	BATCH NUMBER
	USE BY
	MANUFACTURER
	CONTAINS SUFFICIENT FOR <n> TESTS
	<i>IN VITRO</i> DIAGNOSTIC MEDICAL DEVICE
	TEMPERATURE LIMITATION
	CATALOGUE NUMBER
	CONSULT INSTRUCTIONS FOR USE

Procedure in-brief