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1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Desirable Specifications for Total Error, Imprecision, and Bias derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
5. Kasper D.L. et al. (2018). Harrison's principles of internal medicine.20th ed.(pp. 6608-12) Mc Graw Hill Education
6. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests.6th ed. Elsevier Inc.
7. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics. (pp. 1578-80). 6th ed. Elsevier Inc.

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INTENDED USE:

The Quantitative Determination total amount of Binding Sites Available for the Thyroid hormones in human serum or EDTA plasma

CLINICAL SIGNIFICANCE:

The thyroid hormone thyroxine (T4) is one of the thyroid hormones which plays role in general metabolism. The measurement of the T4 concentration is of importance in laboratory diagnostics for differential diagnosis between euthyroid, hyperthyroid, and hypothyroid conditions. As the major fraction of the total thyroxine is bound to transport proteins (TBG, transthyretin or prealbumin, and albumin), the measurement of thyroxine only provides correct information about total thyroxine. Any alteration in thyroid binding proteins specially TBG can affect total thyroxine, while the amount of free T4 which is physiologically active form of hormone, remains constant. A change in the TBG concentration can lead to elevated or lowered total T4 concentrations being measured although the free T4 concentration is in the euthyroid range. A T-uptake result provides useful information about serum thyroid binding capacity and therefore free thyroid index (FTI). In hyperthyroidism or hypothyroidism both T4 and T-up elevate or decrease, but in conditions in which there is a binding globulin abnormality, these two values diverge (e. g. in pregnancy or estrogen therapy total T4 increase while T-uptake decrease).

Measurement Principle:

First specified amount of thyroxine and HRP-T4 conjugate and patient serum as a source of binding proteins(P) are mixed with monoclonal antibody against T4(Ab) immobilized on microtiter plate. following mixing the enzyme-conjugate and thyroxine with the specimen, a binding reaction results between the patient's binding proteins and the added thyroxine but not with the enzyme conjugate.

this interaction is represented below:

The added T4 not consumed in reaction (1) compete with T4-HRP for limited site on immobilized solid phase antibody. The reaction can be explained by following equation:

In this equation T_4 is unreacted thyroxine in reaction(1)and Ab_{sp} is solid phase immobilized antibody. After equilibrium is attained, bound fraction separated from unbound $HRPT_4$ by washing. Whatever the binding capacity of patient serum is greater, more T4 consumed in reaction (1), less T_4Ab_{sp} formed in reaction (2) and therefore more $HRPT_4Ab_{sp}$ complexes produce more color (hypothyroidism). On the contrary, less binding capacity (Hyperthyroidism), less T4 consumed, and less $HRPT_4$ complexes produced and therefore less color is developed.

T-Uptake ELISA KIT		IVD-REF: PS – T.UP
		Ver. No: 02

Reagent and material supplied in the kit:

Reagent	Quantity/96 Test KIT	Reconstitution
Microplate (96 breakable wells) coated with anti T4	96 wells	Ready to Use
Calibrators A-D (Exact T-uptake percentage printed on the label) in buffer, containing preservative	6 Vial/0.5ml	Ready to Use
Controls N = 1 in buffer compatible with human serum matrix with preservatives	1 Vial/0.5 ml	Ready to Use
Enzyme Conjugate (red color)	1 Vial/12 ml	Ready to use
Concentrated Wash solution	1 Vial/30 ml	Dilute 20X with distilled water
Substrate-Chromogenic solution (Tetramethylbenzidine and hydrogen peroxide)	1 Vial/12ml	Ready to Use
Stop solution (H2So4 0.12M) Additional information about other additive is considered proprietary.	1 Vial/6.0 ml	Ready to Use

Materials required but not supplied with kit

The following materials are required but not provided in the kit:

- Distilled water with conductivity < 1 μ S/cm
- Precise micropipettes for delivery of 25, 50, 100, microliters. An 8-channel pipette or respenser pipette with disposable plastic tips for delivery of 25-100 μ L is useful but not essential.
- Pipettes for dispensing milliliter volumes.
- Vortex mixer
- Automatic washer or any other device which can deliver 350 microliter wash solution for microplates
- Microplate reader with a wavelength of 450 nm and 620-650 nm, and an absorbance range of 0-3.0

STORAGE AND EXPIRATION DATING OF REAGENTS

- If kept at 2-8 $^{\circ}$ C, before opening or reconstitution, all kit components are stable until the expiry date as indicated on their label.
- Calibrators concentration is displayed on the vial labels and may vary between lots.
- Do not freeze kit components.
- Freshly prepared Working Wash solution should be used on the preparation day.
- Each lot of reagents and calibrators has been standardized to produce the correct reaction. Do not interchange the reagents or calibrators between lots.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration. Chromogenic solution TMB should be colorless. Blue color indicates that the reagent has been contaminated and should be discarded.

SPECIMEN COLLECTION AND PREPARATION

- Human serum and EDTA plasma are samples of choice for PSI T-uptake EIA kit.
- Do not use grossly haemolysed or turbid samples.

T-Uptake ELISA KIT		IVD-REF: PS – T.UP
		Ver. No: 02

- If the measurement cannot be done within 24 hours after sample collection, samples should be kept at 2-8°C up to 48 hours. For longer storage serum or plasma must be stored at -20°C.
- Avoid repeated freeze and thaw cycles. Frozen specimens must be gently but thoroughly mixed after thawing.

PRECAUTIONS AND WARNINGS

Safety

- For Professional in vitro diagnostic use only.
- All the reagents within the kit are calibrated for the direct determination of T-Uptake in human serum. The kit is not calibrated for the determination of T-uptake in Non-EDTA plasma or other specimens of human or animal origin.
- Before starting the assay, read the working manual completely and carefully. Use the valid version of the instruction provided with the kit. Be sure that everything is understood.
- Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., USA) publication No (CDC) 88—8_95 on laboratory safety procedures or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Reagents contain sodium azide (NaN₃) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. Follow local guidelines for disposal of all waste material.
- Avoid any skin contact with all reagents. Stop solution contains strong acid in case of contact, wash thoroughly with soap and water.

Caution

Material used in the preparation of human source reagents has been tested and found to be non-reactive for HIV 1 & 2 antibodies, HCV antibody and hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they are potentially infectious.

Preparation of reagents

Wash Buffer:

According to the required volume, pour the desired volume of wash concentrate into a clean container and dilute 20 fold by adding adequate volume of distilled or deionized water to gain buffered Wash Solution. E.g., 10 mL wash concentrate to 190 mL distilled water. Wash solution should be prepared daily and discard unused working wash solution at the end of the day.

HANDELING NOTES

- Do not use the kit or components beyond the expiry date.
- Do not mix up vial caps.

T-Uptake ELISA KIT		IVD-REF: PS – T.UP
		Ver. No: 02

- Bring all the reagents to room temperature prior to use.
 - Thoroughly mix all reagents, calibrators/controls and samples by gentle agitation or swirling.
 - Assay calibrators, controls and samples in duplicate. Vertical alignment is recommended.
 - Use a clean container and high quality distilled or deionized water to prepare the working wash solution.
 - In order to avoid cross-contamination, use a clean disposable pipette tip for addition of each reagent and sample.
 - For the dispensing of Chromogenic solution and Stop solution, avoid pipettes with metal parts.
 - High precision pipettes or automated pipetting equipment will improve the precision.
 - Adhere the incubation times carefully. Tolerance limit of incubation time is $\pm 5\%$
1. To avoid drift, arising from low volume sample evaporation or incubation time difference between first and last well, it is strongly recommended not to run more than six strips when pipetting manually (Time delay).
 2. Perform all steps consecutively without any delay. Any delay in each step has deleterious effects on final results.
 3. Draw a calibration curve for each run, do not use stored curve from previous runs.
 4. Incubation with Chromogenic solution must be done in the dark.
 5. Opened reagents are stable maximally for 8 weeks if they are not contaminated, and stored in resealed original containers and handled as prescribed. Return to 2-8°C immediately after use.

Measuring Procedure

1. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminum pouch containing a desiccant reseal carefully and store at 2-8°C).
2. Pipette 25 μL of each Calibrator (calibrators A to D), controls and patient samples in duplicate into the appropriate wells.
3. Pipette 100 μL of enzyme conjugate into each well. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid.
4. Gently mix for 15 seconds
5. Cover the plate with an adhesive seal and incubate for 45 minutes at room temperature (20-27°C).
6. Wash the plate according to procedure below.
 - i. Carefully dispense 350 μL of working wash solution into each well.
 - ii. Aspirate the content of each well.
 - iii. Repeat four additional wash for five wash cycles.
 - iv. After the last wash cycle invert the plate and tap it carefully against an absorbent paper.

Note: Automatic strip washer: Follow the manufacturer's instructions for maintenance and wash the required number of wash cycles.

T-Uptake ELISA KIT		IVD-REF: PS – T.UP
		Ver. No: 02

7. Add 100 µL of Chromogenic solution TMB to each well using the same pipetting technique as in procedural note 3. The Chromogenic solution TMB should be added to the wells as quickly as possible.
8. Incubate for 15 min at room temperature preferentially in the dark.
9. Add 50 µL of Stop Solution, mix and read the absorbance at 450 nm (with a reference filter 620, 630, or 650 nm) in a microplate spectrophotometer within 10 min after addition of Stop Solution.

Important Note: always add stop solution in the same order as chromogenic solution TMB.


CALCULATION OF RESULTS

1. Construct a calibration curve by plotting the average O.D. (450 nm) for each calibrator on the vertical (Y) axis against the T-uptake percentage on the horizontal (X) axis, and draw the best-fit curve through the points on the graph.
2. To determine T-uptake percentage in each sample, first locate the O.D. value on the Y-axis and extend a horizontal line to the calibration curve. At the point of intersection, draw a vertical line to the X-axis and read the corresponding T-uptake percentage.
3. If a microplate spectrophotometer with built-in data calculation program is used, refer to the manual for the spectrophotometer and create a program using the concentration stated on the label of each calibrator. For automatic calculation of T-UPTAKE results, it is recommended to use point-to-point method.

Typical Data

The following data are for illustration only and should never be used instead of the real-time calibration curve.

ROW	OD	T-UPTAKE
1	2.5	14
2	2.2	24
3	0.9	34
4	0.5	44

Quality control

T-Uptake control serum should be used for validation of the assay series. Range of expected results are indicated on the vial label. If assigned range is not obtained, a complete check of reagents and reader performance should be made and the assay shall be repeated. Each laboratory can also prepare its own serum pools at different levels, which can be used as internal controls in order to ensure the precision of the assay.

T-Uptake ELISA KIT		IVD-REF: PS – T.UP
		Ver. No: 02

Reference interval:

A study with PSI T-uptake EIA kit was performed using samples from some apparently healthy people, following reference interval corresponding to 2.5th and 97.5th percentiles of results was extrapolated. Each laboratory should first investigate the transferability of the recommended expected value to its own patient population and if it does not match, determine its own reference interval.

Free thyroid index (FTI) as an indirect estimate of free thyroxine can also be calculated from following formula:

$$FTI = \frac{T - uptake \times Total T_4}{Mean T - uptake of Reference Lab}$$

You can calculate mean T-up by adding upper limit and lower limit of reference interval of T-uptake and divide the sum by two.

Analyte	Number	Mean (%)	2.5 th percentile (%)	97.5 percentile (%)
T-uptake	137	30	25	35
FTI	-		5.4	9.7

Precision:

Precision was determined using PSI T-uptake EIA kit reagents and three pooled human sera according to EP05-A3 Guideline.

Three runs per day (morning, midday and afternoon) by three user in 10 working days (10 × 3 × 3). Results analyzed by fully nested ANOVA. Following results were obtained:

Sample Description	Mean (%)	Repeatability		Within-Laboratory Precision	
		SD	%CV	SD	%CV
Patient Pool	22.4	1.9	8.48	1.8	8.04
Patient Pool	33.1	2.9	8.76	1.92	5.80
Patient Pool	43.2	2.11	4.88	2.73	6.32

specificity

Interference evaluation from common sample abnormalities, hemolysis, Lipemia and icterus were performed according to guideline CLSI EP07-A2 and percentage interference was calculated for each sample using the equation below and normalized to the sample ferritin content.

$$\%interference = \frac{measured\ value - true\ value}{true\ value} \times 100\%$$

No interference were seen up to following listed concentration:

T-Uptake ELISA KIT		IVD-REF: PS – T.UP
		Ver. No: 02


Interferent	source	Final concentration of interferent in test serum
Hemoglobin	Fresh haemolysed blood	50 mg/mL
Conjugated Bilirubin	Commercially crystalline Unconjugated Bilirubin	20 mg/dL
Triglycerides	Commercially Intralipid 10%	1000 mg/dL

Method comparison (Trueness)

A comparison of the PSI T-uptake EIA kit (y axis) using patient samples (n=52 range:14-43%) with Elecsis T-uptake Cobas E411-Roche diagnostic (X axis) gave the following results by Ordinary Linear regression:

$$PSI = 0.7625 + 0.9667ECLIA$$

$$R^2 = 0.9940 \text{ \& } R = 0.9970$$



T-Uptake ELISA KIT		IVD-REF: PS – T.UP
		Ver. No: 02


References:

1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Desirable Specifications for Total Error, Imprecision, and Bias derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
5. Kasper D.L. et al. (2015). Harrison’s principles of internal medicine.19th ed. Mc Graw Hill Education
6. Pagana K.D. et al. (2018). Mosby’s Manual of diagnostic and Laboratory tests.6th ed. Elsevier Inc.
7. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. Elsevier Inc.

SYMBOLS:

	BATCH NUMBER
	USE BY
	Date of Manufacture
	CONTAINS SUFFICIENT FOR <n> TESTS
	IN VITRO DIAGNOSTIC DEVICE
	TEMPERATURE LIMITATION
	CATALOGUE NUMBER
	CONSULT INSTRUCTIONS FOR USE

Procedure in Brief

Read at
450/630 nm