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Sample	Original concentration	10 ng/ml added		25 ng/ml added		50 ng/ml added	
		Recovery%	Bias%	Recovery%	Bias%	Recovery%	Bias%
1	34.8 ng/ml	109	1.84%	108	1.65%	109	1.84%
2	54 ng/ml	112	1.67%	115	2.04%	111	1.55%
3	76 ng/ml	109	0.96%	114	1.43%	112	1.25%

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
1. Adie Viljoen, et al, Analytical Quality Goals for 25-Vitamin D Based on Biological Variation, (2011). Journal of Clinical Laboratory Analysis 25 : 130–133
2. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
3. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
4. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
5. Marcus J.S. et al.(2006). Dynamics of Bone and Cartilage Metabolism. 2nd ed. Elsevier Inc.
6. Norman AW (August 2008). "From vitamin D to hormone D: fundamentals of the vitamin D endocrine system essential for good health". Am. J. Clin. Nutr. 88 (2): 491S–499S. PMID 18689389.
7. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory

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9. Royal college of pathologists of Australasia.

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INTENDED USE:

25-hydroxy (25-OH) Vitamin D ELISA KIT is intended for the quantitative determination of 25-hydroxyvitamin D in human serum and plasma.

INTRODUCTION:

Vitamin D is a fat-soluble vitamin responsible for maintaining calcium and phosphorus levels in blood by promoting their absorption from food in the intestines, and by promoting reabsorption of calcium in the kidneys. This enables normal mineralization of bone, needed for bone growth and bone remodeling.

In human, two major forms of vitamin D are the vitamin D3 (cholecalciferol) and vitamin D2 (ergocalciferol). The term vitamin D also refers to "hydroxy"-metabolites of these substances. Vitamin D2 is provided by dietary sources. Because only fish is naturally rich in vitamin D, most of the Vitamin D2 intake in industrialized world is from fortified products including milk and a few other food supplements.

Vitamin D3 is produced in skin exposed to sunlight, specifically ultraviolet B(UV-B) radiation. In this scenario 7-dehydrocholesterol reacts with UVB to produce vitamin D3. Melanin functions as a light filter in the skin so individuals with dark skin and high melanin content require more time in the sunlight to produce the same amount of vitamin D as individuals with lower melanin content. Each Vitamin D produced in the skin or consumed in food, is converted to 25-hydroxyvitamin D and 1,25-dihydroxyvitamin D in the liver and kidneys, respectively. 1,25-dihydroxyvitamin D is physiologically active form of vitamin D. Dihydroxy vitamin D has a short half-life and there are some problems with its measurement, but 25(OH) vitamin D has longer half-life and most clinicians have accepted its measurement as a good indicator of human vitamin D status.

Vitamin D deficiency can be a result of inadequate intake, inadequate sunlight exposure, malabsorption syndromes, liver or kidney disorders, or by a number of metabolic hereditary disorders. Its deficiency results in impaired bone mineralization and leads to bone softening diseases (rickets in children and osteomalacia in adults). Vitamin D deficiency may also contribute to development of osteoporosis. Vitamin D signaling to the cell is thought to be involved in cell proliferation/apoptosis and cell differentiation, so decrease in serum vitamin D level may be associated with colon, breast and pancreas cancer. Vitamin D deficiency is associated with an increase in blood pressure and cardiovascular risk. It also affects the immune system through cell signaling in monocytes and activated T and B cells.

PRINCIPLES OF THE METHOD:

The Pishgaman Sanjesh Isatis total 25-OH Vitamin D ELISA kit is a saturation competitive Enzyme Linked Immunosorbent assay performed on microplates coated with monoclonal antibodies against 25(OH) vitamin D. During the first 60-minutes incubation at 37°C temperature, total 25-OH Vitamin D (D2 and D3) present in samples or calibrators/controls are dissociated from binding proteins and bind to the solid phase coated antibody. At the same time, streptavidin conjugated with HRP is added to the wells. During the 2nd incubation phase, biotinylated 25(OH) vitamin D is added to the microtiter wells without washing step. The greater the amount of 25-hydroxyvitamin D in the sample, the more places of antibody were covered on the solid phase,

and the biotinylated 25-hydroxyvitamin D had lower opportunities to bind to the solid phase, so the less biotinylated 25(OH) vitamin D bound to streptavidin-HRP conjugate can bind to the solid-phase monoclonal antibody, resulting in the production of few signal. After washing, buffered Substrate/Chromogen reagent (hydrogen peroxide and tetra-methylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction, a blue color will develop and the intensity of the color is reversely proportional to the amount of 25(OH) vitamin D present in the samples.

REAGENTS PROVIDED:

Reagent	Quantity/48 Test KIT	Quantity/96 Test KIT	Quantity/192 Test KIT	Reconstitution
Microplate (96 breakable wells) coated with anti 25-OH Vitamin D	48 wells	96 wells	192 wells	Ready to Use
Calibrators 1-6 (see exact values on vial labels) in Fetal Bovine Serum with preservatives	6 Vial/0.5 ml	6 Vial/0.5 ml	6 Vial/1.0 ml	Ready to Use
2 Serum controls, in Hormone free Fetal Bovine Serum with preservatives	2 Vial/0.5 ml	2 Vial/0.5 ml	2 Vial/1.0 ml	Ready to Use
Conjugate diluent (blue color)	1 Vial/3.0 ml	1 Vial/7 ml	1 Vial/13 ml	Ready to use
enzyme Conjugate Concentrate 11X (streptavidin-horseradish peroxidase conjugate in preservatives): pale yellow color	1 Vial/0.5 mL	1 Vial/0.75 mL	1 Vial/ 1.5 mL	Dilute 11X with conjugate diluent
Biotin solution(25-OH- vitamin D conjugated to Biotin in preservative): Deep yellow color	1 Vial/3.0 ml	1 Vial/7 ml	1 Vial/13 ml	Ready to Use
Concentrated Wash solution	1 Vial/30ml	1 Vial/30 ml	1 Vial/50 ml	Dilute 1:20 with distilled water
Substrate-Chromogenic solution (Tetra-methylbenzidine and hydrogen peroxide)	1 Vial/6.0ml	1 Vial/12ml	2 Vial/12 ml	Ready to Use
Stop solution (H ₂ So ₄ 0.12M) Additional information about other additive is considered proprietary.	1 Vial/6.0 ml	1 Vial/6.0 ml	1 Vial/12.0 ml	Ready to Use

Materials required but not supplied with kit

The following materials are required but not provided in the kit:

- Distilled water with conductivity < 1µS/cm
- Precision micropipettes for delivery of 25, 50, 100 microliters. An 8-channel pipette or resenser pipette with disposable plastic tips for delivery of 50-100 µL is useful but not essential. Pipettes for dispensing milliliter volumes.
- Vortex mixer
- Automatic washer or any other device which can deliver 350 microliter wash solution for microplates
- Microplate reader with a wavelength of 620-650 nm and 450 nm, and an absorbance range of 0-3

STORAGE AND EXPIRATION DATE OF REAGENTS

- If kept at 2-8 °C, before opening or reconstitution, all kit components are stable until the expiry date as it is indicated on their label.
- **Freshly prepared Working conjugate solution should be used within 30 minutes after preparation. Do not prepare more working conjugate solution than the amount which will be used within this period and unused remaining solution must be discarded.**
- Freshly prepared Working Wash solution should be used on the preparation day.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration. Chromogenic solution TMB should be colorless. Blue color indicates that the reagent has been contaminated and should be discarded.

SPECIMEN COLLECTION AND PREPARATION

- Human serum and EDTA plasma are samples of choice for PSI 25(OH) vitamin D kit.
- If the measurement cannot be done within 24 hours after sample collection, samples should be kept at 2-8 °C up to 48 hours. For longer storage serum or plasma must be stored at -20°C.
- Avoid repeated freeze and thaw cycles. Frozen specimens must be gently but thoroughly mixed after thawing.

PRECAUTIONS AND WARNINGS

Safety

- For Professional in vitro diagnostic use only
- Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., USA) publication No (CDC) 88—8_95 on laboratory safety procedures or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Reagents contain sodium azide (NaN₃) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. Follow local guidelines for disposal of all waste material.
- All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal products should be treated as potentially infectious.
- Avoid any skin contact with all reagents, stop solution contains strong acid. In case of contact, wash thoroughly with soap and water.

Caution

Materials used in the preparation of human source reagents have been tested and found to be non-reactive for HIV 1 & 2 antibodies, HCV antibody and hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

Preparation of reagents:

Wash Buffer:

According to volume required, pour the desired volume of wash concentrate into a clean container and dilute 20 fold by adding adequate volume of distilled or deionized water to gain a buffered

Wash Solution. E.g. 10 mL wash concentrate to 190 mL distilled water. Wash solution should be prepared daily and discard unused working wash solution at the end of the day.

Working Conjugate solution:

Prepare the required quantity of working conjugate solution by adding 100 μ L of Enzyme Conjugate Concentrate 11X to 1 mL of conjugate diluent. Mix gently but thoroughly. Working conjugate solution must be used within 30 minutes after preparation. Following table can be used as a guide to prepare the working conjugate solution:

Conjugate diluent (mL)	Concentrated conj.(11X) (μ L)	No of strips
1	100	2
2	200	4
3	300	6
4	400	8
5	500	10
6	600	12

Caution: Conjugate diluent contains some ingredients that may precipitate in cold, so it is strongly recommended that you **bring conjugate diluent to the room temperature, before preparing the working conjugate solution.**

Note: Working conjugate solution is highly sensitive to metal ions, **so it is strongly recommended to use disposable plastic tubes to prepare the working conjugate solution.**

HANDELING NOTES:

- Do not use the kit or components beyond expiry date.
- Do not mix up vial caps.
- Do not mix reagents from different kit lots.
- Bring all the reagents to the room temperature prior to use. Conjugate diluent contains some ingredients that may precipitate in cold, so it is strongly recommended that you **bring conjugate diluent to the room temperature, before preparing the working conjugate solution.**
- Thoroughly mix all reagents, calibrators/controls and samples by gentle agitation or swirling.
- Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
- Use a clean container and high-quality distilled or deionized water to prepare the Working Wash solution.
- In order to avoid cross-contamination, use a clean disposable pipette tip for addition of each reagent and sample.
- For dispensing of the Chromogenic solution and Stop solution avoid pipettes with metal parts.
- High precision pipettes or automated pipetting equipment will improve the precision.
- Adhere the incubation times carefully. Tolerance limit of incubation time is $\pm 5\%$.

- To avoid drift, the time between pipetting of the first calibrator and the last sample must not be more than 10 minutes (Time delay). It is strongly recommended not to run more than 5 strips when pipetting manually.
- Perform all steps consecutively without any delay. Any delay in each step has deleterious effects on the final results.
- Draw a calibration curve for each run, do not use stored curve from previous runs.
- Incubation with Chromogenic solution, must be done in the dark.
- Opened reagents are stable maximally for 8 weeks if they are not contaminated, stored in resealed original containers, and handled as prescribed. Return to 2-8°C immediately after use

Measurement Procedure:

1. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminum pouch containing a desiccant, reseal carefully and store at 2-8°C.
2. Pipette 25 µL of each Calibrator (calibrators 0, 5, 10, 30, 60, 120 ng/mL), controls and patient samples in duplicate into the appropriate wells.
3. Pipette 50 µL of working conjugate solution into each well. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid.
4. Gently mix for 15 seconds
5. Cover the plate with an adhesive seal and incubate for 60-minutes at 37°C.
6. After 1st incubation, add 50 µL biotin solution to each well ***without washing***.
7. Gently mix for 15 seconds.
8. Cover the plate with an adhesive seal and incubate it for 30 minutes at 37°C ± 1°C.
9. Wash the plate according to the procedure below.
 - i. Carefully dispense 350 µL of working wash solution into each well.
 - ii. Aspirate the content of each well.
 - iii. Repeat four additional wash for a total of 5 wash cycles.
 - iv. At the end of the last wash cycle, invert the plate and tap it carefully against an absorbent paper.

Note: Automatic strip washer: Follow the manufacturer's instructions for maintenance and wash the required number of wash cycles after second incubation step.

10. Add 100 µL of Chromogenic solution TMB to each well using the same pipetting technique as in item 3. The Chromogenic solution TMB should be added to the wells as quickly as possible and the time between addition to the first and last well should not exceed 5 min.
11. Incubate for 15 min at room temperature preferably in the dark. Avoid exposure to direct sunlight.
12. Add 50 µL of Stop Solution, mix and read the absorbance at 450 nm in a microplate spectrophotometer within 10 min after addition of Stop Solution.

Important Note: always add stop solution in the same order as chromogenic solution TMB.

CALCULATION OF RESULTS


1. Construct a calibration curve by plotting the average O.D. (450 nm) for each calibrator on the vertical (Y) axis against the concentration on the horizontal (X) axis, and draw a best fit curve through the points on the graph.
2. To determine concentration of 25(OH) Vit D in each sample, first locate the O.D. value on the Y-axis and extend a horizontal line to the calibration curve. At the point of intersection, draw a vertical line to the X-axis and read the corresponding concentration.
3. If a microplate spectrophotometer with built-in data calculation program is used, refer to the manual of the spectrophotometer and create a program using the concentration stated on the label of each Calibrators. For automatic calculation of 25(OH) vitamin D results, it is recommended to use point-to-point method.

Typical Data:

The following data are for illustration only and should never be used instead of the real-time calibration curve.

Total 25-OH Vitamin D	OD Units
0.0 ng/ml	2.9
5.0 ng/ml	2.1
10.0 ng/ml	1.4
30.0 ng/ml	0.6
60.0 ng/ml	0.3
120.0 ng/ml	0.15

Note: 1 ng/ml=2.5 pmol/L

Quality control:

25(OH) vitamin D low and high control serums should be used for validation of the assay series. Ranges of expected results are indicated on the vial labels. If assigned range is not obtained, a complete check of reagents and reader performance should be made and the assay shall be repeated. Each laboratory can also prepare its own serum pools at different levels, which can be used as internal controls in order to ensure the precision of the assay

Health based reference values

Currently there is no absolute definition of the optimal vitamin D status. It should be taken into account that differences in 25(OH) D levels may exist with respect to gender, age, season, geographical latitude and ethnic groups.

Each laboratory should first investigate the transferability of the recommended expected value to its own patient population and if it doesn't match, determine its own reference interval. Population based reference interval should not be taken as clinical guide to recommend vitamin D supplementation. The values given here is based on the published data, gathered from healthy population from IR Iran-Tehran and must be considered as information only.

LEVEL	REFERENCE VALUE
DEFICIENT	<10 ng/ml
INSUFFICIENT	10-30 ng/ml
SUFFICIENT	31-100 ng/ml
TOXICITY	>100 ng/ml

PERFORMANCE CRITERIA:

1. Lower limit of measurement:

The Limit of Blank and Limit of Detection were determined according to CLSI (Clinical and Laboratory Standards Institute) EP17-A2 guideline.

In brief, Limit of Blank was determined as the 95th percentile value from $n = 60$ measurements on standard 0 over several runs. The Limit of Blank equals to the concentration that 95% of repeated measurements on zero standard give value below it. In the next step we determined standard deviation of repeated measurements on 3 samples with Vitamin D contents less than 4 ng/mL determined from a separate method study (HPLC Method). By adding the LOB calculated at 1st step to 1.645 times the greatest standards deviation calculated in the 2nd step ($LOD = LOB + 1.645SD_L$) LOD calculated as 4.0 ng/mL.

2. Precision:

Precision was determined using PSI 25(OH) Vitamin D reagents and 3 pooled human sera ranging in different points of measuring range (4-120 ng/mL) according to EP05-A3 Guideline.

In brief, 2 human pooled serum prepared without any modification and another one prepared by

spiking an extra human serum pooled. 3 runs per day (morning, midday and afternoon) by 3 user in 10 working days (10 × 3 × 3). Results analyzed by fully nested ANOVA using a spreadsheet. The results are summarized in table below:

Sample Description	Mean (ng/mL)	Repeatability		Within-Laboratory Precision	
		SD	%CV	SD	%CV
Patient Pool	9.38817	1.694	18.04	1.857	19.78
Patient Pool	32.8041	3.128	9.53	4.849	14.78
Spiked patient Pool	91.3913	5.098	5.58	8.939	9.78

3. Specificity:

Specificity evaluation was performed according to guideline CLSI EP07-A2 to assess the cross-reactivity of the assay with other vitamin D metabolites. Samples containing cross-reactants were prepared by spiking natural human sera with material containing high levels of cross-reactants. Modification of serum matrix was no more than 10%. The % cross-reactivity was calculated for each sample using the equation below and normalized to the 25(OH) vitamin D3:

$$\% \text{cross-reactivity} = \frac{\text{mean conc. of spiked sample} - \text{mean conc. of unspiked sample}}{\text{spiked concentration}} \times 100\%$$

The results are summarized in the following table:

Compound	% Cross reaction
25OH Vitamin D3	100
25OH Vitamin D2	83
Vitamin D3	<0.1
Vitamin D2	<0.1

No interference, from common serum abnormalities, were seen up to following listed concentration:

Interferent	source	Final concentration of interferent in test serum
Hemoglobin	Fresh haemolysed blood	50 mg/mL
Conjugated Bilirubin	Commercially crystalline Unconjugated Bilirubin	20 mg/dL
Triglycerides	Commercially Intralipid 10%	1000 mg/dL

4. Trueness:

4.1 Recovery Evaluation

Recovery was assessed by adding different levels of 25-OH vitamin D to 3 different samples in 3 points of measuring range. The percentage recovery and % bias was calculated for each sample using the equations below:

$$\%Recovery = \frac{\text{Concentration Recovered}}{\text{Concentration Added}} \times 100$$

$$\%Bias = \frac{\bar{X} \text{ measured} - \text{expected}}{\text{Expected}} \times 100$$

The results are summarized in the following table:

Sample	Original concentration	10 ng/ml added		25 ng/ml added		50 ng/ml added	
		Recovery %	Bias%	Recovery %	Bias%	Recovery %	Bias %
1	34.8 ng/ml	109	1.84%	108	1.65%	109	1.84%
2	54 ng/ml	112	1.67%	115	2.04%	111	1.55%
3	76 ng/ml	109	0.96%	114	1.43%	112	1.25%


4.2 Method comparison:

A comparison of PSI 25(OH) vitamin D assay (y axis) using 48 patient samples with concentrations assigned by the HPLC Vitamin D assay (X axis) gave the following correlations:

$$Y = 0.9771x - 0.4658$$


$$r = 0.977$$

$$r^2 = 0.9505$$

REFERENCES:

1. Adie Viljoen, et al, Analytical Quality Goals for 25-Vitamin D Based on Biological Variation,(2011). Journal of Clinical Laboratory Analysis 25 : 130–133
2. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
3. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
4. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
5. Norman AW (August 2008). "From vitamin D to hormone D: fundamentals of the vitamin D endocrine system essential for good health". Am. J. Clin. Nutr. 88 (2): 491S–499S. PMID 18689389.
6. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests.6th ed. Elsevier Inc.
7. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. Elsevier Inc.
8. Royal college of pathologists of Australasia. <https://www.datainnovations.com/allowable-total-error-table>
9. Wolf G (June 2004). "The discovery of vitamin D: the contribution of Adolf Windaus". J. Nutr. 134 (6): 1299–302. PMID 15173387.

SUMMARY OF THE PROTOCOL