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$$\%Bias = \frac{\bar{X}\ measured - expected}{Expected} \times 100$$

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Sample	Original concentration	50 IU/L added		100 IU/L added	
		Recovery%	Bias%	Recovery%	Bias%
1	45.8 IU/L	109.00%	4.70%	109.00%	6.17%
2	345 IU/L	95.00%	-0.63%	95.00%	-1.12%
3	725 IU/L	112.00%	0.77%	98.00%	-1.21%

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
1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Kasper D.L. et al. (2015). Harrison's principles of internal medicine.19th ed. Mc Graw Hill Education
5. Marcilac I. Troalen F. et al.(1992) Free human chorionic gonadotropin beta subunit in gonadal and nongonadal neoplasms. Canc Res. 1992;52:3901-3907.
6. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests.6th ed. Elsevier Inc.
7. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. Elsevier Inc

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INTENDED USE:

For the quantitative determination of human chorionic gonadotropin (hCG) concentration in human serum.

CLINICAL SIGNIFICANCE:

Similarly, to LH, FSH and TSH, human chorionic gonadotropin (hCG) is a member of the glycoprotein family and consists of 2 subunits (α - and β -chains) which are associated to make the intact hormone. The α -chains in all four of these glycoprotein hormones are virtually identical, whereas the β -chains have greatly differing structures and are responsible for the respective specific hormonal functions.

Human chorionic gonadotropin consists of a number of isohormones with differing molecular size. hCG serves to maintain The biological action of the corpus luteum during pregnancy. It also influences steroid production. The serum of pregnant women contains mainly intact hCG.

Measurement of the hCG concentration permits the diagnosis of pregnancy just one week after conception. The determination of hCG in the 1st trimester of pregnancy is of particular importance. Elevated values here serve as an indication of hydatid form mole or multiple pregnancy. Depressed values indicate threatening or missed abortion, ectopic pregnancy or intra-uterine death. Elevated values in the absence of pregnancy are indicative of a tumor. The specific monoclonal antibodies, which are used, recognize only the holo-hormone, therefore, the PSI hCG Rapid ELISA kit should be used in particular in the diagnosis and monitoring of pregnancy.

PRINCIPLE:

The PSI hCG rapid EIA kit is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators, controls and patient samples are first incubated with a biotinylated Anti-hCG monoclonal antibody in a streptavidin-coated microstrips and a bound/free separation is performed by a solid-phase washing. Then during the 2nd incubation, HRP-conjugated Anti-hCG monoclonal antibody which recognizes hCG from different sites is added to microstrips. Existing hCG in calibrators/controls or samples reacts with both antibodies to form a sandwich complex. This complex is adsorbed to the streptavidin coated microstrips via the biotinylated Anti-hCG Mab. The strips are then washed and buffered substrate/chromogen reagent (hydrogen peroxide and tetramethylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue color will develop and the intensity of the color is directly proportional to the amount of hCG present in the samples.

hCG Rapid ELISA KIT		IVD-REF: PS – HCG.R
		Ver. No: 02

Reagent and material supplied in the kit:

Reagent	Quantity/48 Test KIT	Quantity/96 Test KIT	Quantity/192 Test KIT	Reconstitution
Microplate (96 breakable wells) coated with streptavidin	48 wells	96 wells	192 wells	Ready to Use
Standard 0 (Buffer without hCG compatible to human serum matrix containing preservative)	1 ×5.0 mL	1 ×5.0 mL	1 ×5.0 mL	Ready to Use
Calibrators 2-6 (20, 100, 250, 500, 1000 IU/L) in buffer, containing preservative with traceability to reference material “3 rd ISO WHO 75/539”	5 Vial/0.5 ml	5 Vial/1.0 ml	5 Vial/2.0 ml	Ready to Use
Controls N = 1 in buffer compatible with human serum matrix containing preservatives	1 Vial/0.5 ml	1 Vial/1.0 ml	1 Vial/2.0 ml	Ready to Use
Assay Buffer-Biotinylated Anti-hCG	1 Vial/3.0 ml	1 Vial/6 ml	1 Vial/12 ml	Ready to use
Enzyme Conjugate-HRP conjugated Anti-hCG (red color)	1 Vial/6.0 ml	1 Vial/12 ml	2 Vial/12 ml	Ready to use
Concentrated Wash solution	1 Vial/30ml	1 Vial/30 ml	1 Vial/50 ml	Dilute 1:20 with distilled water
Substrate-Chromogenic solution (Tetramethylbenzidine and hydrogen peroxide)	1 Vial/6.0ml	1 Vial/12ml	2 Vial/12 ml	Ready to Use
Stop solution (H ₂ So ₄ 0.12M). Additional information about the other additive is considered proprietary.	1 Vial/6.0 ml	1 Vial/6.0 ml	1 Vial/12.0 ml	Ready to Use

Materials required but not supplied with kit

The following materials are required but not provided in the kit:

- Distilled water with conductivity < 1µS/cm
- Precise micropipettes for delivery of, 50, 100, microliters. An 8-channel pipette or respenser pipette with disposable plastic tips for delivery of 50-200 µL is useful but not essential.
- Pipettes for dispensing milliliter volumes.
- Vortex mixer
- Automatic washer or any other device which can deliver 350 microliter wash solution for microplates
- Microplate reader with a wavelength of 450 nm and 620-650 nm, and an absorbance range of 0-3.0

STORAGE AND EXPIRATION DATING OF REAGENTS

- If kept at 2-8 °C, before opening or reconstitution, all kit components are stable until the expiry date as indicated on their label.
- Calibrators’ concentration is displayed on the vial labels and may vary between lots.
- Do not freeze kits.
- Freshly prepared Working Wash solution should be used on the preparation day.
- Each lot of reagents and calibrators has been standardized to produce the correct reaction. Do not interchange the reagents or calibrators between lots.

hCG Rapid ELISA KIT		IVD-REF: PS – HCG.R
		Ver. No: 02

- Alterations in physical appearance of kit reagents may indicate instability or deterioration. Chromogenic solution TMB should be colorless. Blue color indicates that the reagent has been contaminated and should be discarded.

SPECIMEN COLLECTION AND PREPARATION

- Human serum and EDTA plasma are samples of choice for PSI hCG Rapid kit.
- Do not use grossly haemolysed or turbid samples.
- If the measurement cannot be done within 24 hours after sample collection, samples should be kept at 2-8°C up to 48 hours. For longer storage, serum or plasma must be stored at -20°C.
- Avoid repeated freeze and thaw cycles. Frozen specimens must be gently but thoroughly mixed after thawing.

PRECAUTIONS AND WARNINGS

Safety

- For Professional in vitro diagnostic use only.
- All the reagents within the kit are calibrated for the direct determination of intact hCG in human serum. The kit is not calibrated for the determination of hCG in urine, Non-EDTA plasma or other specimens of human or animal origin.
- Before starting the assay, read the working manual completely and carefully. Use the valid version of the instruction provided with the kit. Be sure that everything is understood.
- Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., USA) publication No (CDC) 88—8_95 on laboratory safety procedures or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Reagents contain sodium azide (NaN₃) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. Follow local guidelines for disposal of all waste material.
- Avoid any skin contact with all reagents. Stop solution contains strong acid in case of contact, wash thoroughly with soap and water.

Caution

Material used in the preparation of human source reagents has been tested and found to be non-reactive for HIV 1 & 2 antibodies, HCV antibody and hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they are potentially infectious.

Preparation of reagents

Wash Buffer:

According to the required volume, pour the desired volume of wash concentrate into a clean container and dilute 20 fold by adding adequate volume of distilled or deionized water to gain buffered Wash Solution. E.g. 10 mL wash concentrate to 190 mL distilled water. Wash solution should be prepared daily and discard unused working wash solution at the end of the day.

HANDELING NOTES

1. Do not use the kit or components beyond the expiry date.

hCG Rapid ELISA KIT		IVD-REF: PS – HCG.R
		Ver. No: 02

2. Do not mix up vial caps.
3. Bring all the reagents to room temperature prior to use.
4. Thoroughly mix all reagents, calibrators/controls and samples by gentle agitation or swirling.
5. Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
6. Use a clean container and high-quality distilled or deionized water to prepare the working wash solution.
7. In order to avoid cross-contamination, use a clean disposable pipette tip for addition of each reagent and sample.
8. For the dispensing of Chromogenic solution and Stop solution, avoid pipettes with metal parts.
9. High precision pipettes or automated pipetting equipment will improve the precision.
10. Adhere the incubation times carefully. Tolerance limit of incubation time is $\pm 5\%$
11. Although, PSI hCG EIA kit uses streptavidin coated plate technology and reaction does not proceed before adding conjugate solution, but to avoid drift, arising from short incubation times or low volume sample evaporation, it is strongly recommended not to test more than 3 strips, each run (Time delay).
12. Perform all steps consecutively without any delay. Any delay in each step has deleterious effects on results.
13. Draw a calibration curve for each run, do not use stored curve from previous runs.
14. Incubation with Chromogenic solution must be done in the dark.
15. Only standard 0 may be used to dilute any high serum samples. The use of any other reagent may lead to false results.
16. Opened reagents are stable maximally for 8 weeks if they are not contaminated, and stored in resealed original containers and handled as prescribed. Return to 2-8°C immediately after use.

PROCEDURE:

1. Take out the microplate with the required number of strips. The left strips are placed in the plastic pouch along with the silica gel bag, and seal it.
2. Pipette 50 μ l of standards, control and sera into the wells.
3. Pipette 50 μ l of Assay Buffer into each well.
4. Mix gently and cover the strips with a plate sealer. Incubate at room temperature for 5 minutes.
5. Wash the plate according to procedure below.
 - i. Carefully dispense 350 μ L of working wash solution into each well.
 - ii. Aspirate the content of each well.
 - iii. Repeat 4 additional wash for a total of 5 wash cycles.
 - iv. After the last wash cycle invert the plate and tap it carefully against an absorbent paper.

Note: Automatic strip washer: Follow the manufacturer’s instructions for maintenance and wash the required number of wash cycles

6. Pipette 100 μ l of Enzyme conjugate into each well and Incubate at room temperature for 5 minutes.
7. Wash as described in procedural note 5.


8. Pipette 100 µl of Substrate-Chromogenic solution into each well and Incubate in the dark at room temperature for 5 minutes
9. Add 50 µL of Stop Solution, mix and read the absorbance at 450 nm (with 620, 630, or 650 nm as reference filter) in a microplate spectrophotometer within 10 min after addition of Stop Solution.

Important Note: always add stop solution in the same order as chromogenic solution TMB.

Calculation of Results:

1. Construct a calibration curve by plotting the average O.D. (450 nm) for each calibrator on the vertical (Y) axis against the concentration on the horizontal (X) axis, and draw a best-fit curve through the points on the graph.
2. To determine concentration of hCG in each sample, first locate the O.D. value on the Y-axis and extend a horizontal line to the calibration curve. At the point of intersection, draw a vertical line to the X-axis and read the corresponding concentration.
3. If a microplate spectrophotometer with built-in data calculation program is used, refer to the manual for the spectrophotometer and create a program using the concentration stated on the label of each calibrator. For automatic calculation of hCG results, it is recommended to use point-to-point method.
4. The following data are for illustration never should be used instead of the real-calibration curve.

Standards IU/L	OD 450nm
0	0.025
20	0.110
100	0.38
250	0.8
500	1.45
1000	2.6

only and time

Quality control:

hCG control serum should be used for validation of the assay series. Ranges of expected results are indicated on the vial labels. If assigned range are not obtained, a complete check of reagents and reader performance should be made and the assay shall be repeated. Each laboratory can also prepare its own serum pools at different levels, which can be used as internal controls in order to ensure the precision of the assay.

Expected values:

A study with PSI hCG Rapid ELISA kit was performed using samples from some healthy non-menopausal non-pregnant and pregnant women and also referring to some published scientific literature following reference interval corresponding to 2.5th and 97.5th percentiles of results was extrapolated.

hCG Rapid ELISA KIT		IVD-REF: PS – HCG.R
		Ver. No: 02

Each laboratory should first investigate the transferability of the recommended expected value to its own patient population and if it does not match, determine its own reference interval:

Pre-menopausal Non-pregnant women	N=158	<10 IU/L
Suspicious	-	10-25 IU/L
Pregnancy	N=120	>25 IU/L

Pregnant women in different pregnancy weeks

Weeks of pregnancy	2.5th percentile(IU/L)	97.5th percentile(IU/L)
4 weeks	5.0	100
5 weeks	200	3000
6 weeks	10000	80000
7-14 weeks	90000	500000
15-26 weeks	5000	80000
27-40 weeks	3000	15000

PERFORMANCE AND CHARACTERISTICS:

Lower limit of measurement:

The Limit of Blank and Limit of Detection were determined according to CLSI (Clinical and Laboratory Standards Institute) EP17-A2 guideline.

In brief, Limit of Blank was determined as the 95th percentile value from n = 60 measurements on standard 0 over several runs. The Limit of blank equals to the concentration that 95% of repeated measurements on zero standard give value below it. In the next step, we determined standard deviation of repeated measurements on three samples with hCG contents less than 2.0 IU/L determined from a separate method study (ECLIA Method). By adding, the LOB calculated at first step to 1.645 times the greatest standards deviation calculated in the second step ($LOD = LOB + 1.645SD_L$) LOD calculated as 1.0 IU/L.

Specificity:

Specificity evaluation was performed according to guideline CLSI EP07-A2 to assess the cross-reactivity of the assay with common hCG cross reactants. Samples containing cross-reactants was prepared by spiking natural human sera with material containing high levels of cross-reactants. Modification of serum matrix was no more than 10%. The % cross-reactivity was calculated for each sample using the equation below and normalized to the hCG content:

$$\%cross - reactivity = \frac{mean\ conc.\ of\ spiked\ sample - mean\ conc.\ of\ unspiked\ sample}{spiked\ concentration} \times 100\%$$

The results are summarized in the following table:

Cross-reactants	concentration	% Cross-reaction
------------------------	----------------------	-------------------------

hCG Rapid ELISA KIT		IVD-REF: PS – HCG.R
		Ver. No: 02

FSH	500 m IU/ml	0.05
LH	500 m IU/ml	0.37
TSH	500 µ IU/ml	0.18

No interference, from common serum abnormalities, were seen up to following listed concentration:

Interferent	source	Final concentration of interferent in test serum
Hemoglobin	Fresh haemolysed blood	50 mg/mL
Conjugated Bilirubin	Commercially crystalline Unconjugated Bilirubin	20 mg/dL
Triglycerides	Commercially Intralipid 10%	1000 mg/dL

Precision and reproducibility:

Precision was determined using PSI hCG Rapid ELISA kit reagents and three pooled human sera according to EP05-A3 Guideline.

Three runs per day (morning, midday and afternoon) by 3 users in 10 working days (10 × 3 × 3). Results analyzed by fully nested ANOVA. Following results were obtained:

Sample Description	Mean (IU/L)	Repeatability		Within-Laboratory Precision	
		SD	%CV	SD	%CV
Patient Pool	78.2	5.9	7.54	4.2	5.37
Patient Pool	353.8	20.8	5.88	8.1	2.29
Patient Pool	850.6	39.1	4.6	38.3	4.5

Recovery Evaluation (Proportional error)

Recovery was assessed by adding different amounts of hCG to 3 samples with different concentrations. Then recovery and bias was calculated for each sample using the equations below:

$$\%Recovery = \frac{\text{Concentration Recovered}}{\text{Concentration Added}} \times 100$$

$$\%Bias = \frac{\bar{X} \text{ measured} - \text{expected}}{\text{Expected}} \times 100$$

The results are summarized in the following table:

Sample	Original concentration	50 IU/L added		100 IU/L added	
		Recovery%	Bias%	Recovery%	Bias%
1	52	108%	3.92%	111%	7.24%
2	337	95%	-0.65%	93%	-1.60%
3	845	110%	0.56%	106%	0.63%

hCG Rapid ELISA KIT

IVD-REF: PS – HCG.R

Ver. No: 02


Method comparison:

Comparison of the PSI hCG Rapid kit (y axis) using patient samples (n=88 941 IU/L) with ECLIA-Cobas E411-diagnostic(HCG stat) (X axis) gave the following results by Linear regression

$$PSI\ ELISA = 0.9803ECLIA - 1.4885$$

$$r = 0.9988$$

$$r^2 = 0.9967$$

ELISA
range:7-
Roche


Hook effect:

There is no high dose “hook effect” up to the hCG concentration of 1000000 IU/L.

References:

1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Desirable Specifications for Total Error, Imprecision, and Bias, derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
5. Kasper D.L. et al. (2015). Harrison’s principles of internal medicine.19th ed.(page. 124e-5) Mc Graw Hill Education
6. Marcilac I. Troalen F. et al.(1992) Free human chorionic gonadotropin beta subunit in gonadal and nongonadal neoplasms. Canc Res. 1992;52:3901-3907.

- Pagana K.D. et al. (2018). Mosby’s Manual of diagnostic and Laboratory tests.6th ed. Elsevier Inc.
- Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. Elsevier Inc.
- Royal college of pathologists of Australasia. <https://www.datainnovations.com/allowable->

[total-error-table](#)

Procedure in brief



hCG Rapid ELISA KIT

IVD-REF: PS – HCG.R

Ver. No: 02

SYMBOLS:

	BATCH NUMBER
	USE BY
	MANUFACTURER
	CONTAINS SUFFICIENT FOR <n> TESTS
	<i>IN VITRO</i> DIAGNOSTIC MEDICAL DEVICE
	TEMPERATURE LIMITATION
	CATALOGUE NUMBER
	CONSULT INSTRUCTIONS FOR USE