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1. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
2. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
3. Kao C.Y, Sheu B, Wu J-J. Helicobacter pylori infection an overview of bacterial virulence factors and pathogenesis. Biomed J(2016); 39(1) 14-23
4. Kasper D.L. et al. (2015). Harrison’s principles of internal medicine.19th ed. (pp. 1038-1041) Mc Graw Hill Education
5. Pagana K.D. et al. (2018). Mosby’s Manual of diagnostic and Laboratory tests.6th ed. (pp. 1049-1051). Elsevier Inc.
6. Stefan Reidel. et al. (2019). Jawetz, Melnick, & Adalbert’ s, Medical Microbiology.28th ed. (pp. 268-272). Mc Graw Hill (LANGE Medical Book)

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**INTENDED USE:**

The Quantitative Determination of specific Antibody against *H.pylori* from IgA class Concentration in Human Serum.

**CLINICAL SIGNIFICANCE:**

Helicobacter pylori is a spiral Gram-negative bacterium (2 - 6.5 µm in size, flagellated) formerly classified as campylobacter pylori. H.pylori which colonizes the human gastric mucosa, has been recognized as the principal cause of duodenitis and duodenal ulcers as well as being strongly associated with type B chronic antral gastritis, gastric ulcers, non-ulcer dyspepsia, gastric carcinoma and mucosa associated lymphoid tissue(MALT) lymphomas. The use of nonsteroidal anti-inflammatory drugs (NSAIDS) causes or aggravates peptic and gastric inflammation and ulceration. Hypersecretory states are a much rarer cause of acid peptic disease. Data gathered by history and physical examination may initially suggest peptic ulcer disease.

Infection rate is high and as many as about 50% worldwide and is more prevalent in countries with low socioeconomic status. Approximately 10% of healthy persons younger than 30 years of age have H. pylori without disease or symptoms. Gastric “colonization” increases with age, with people older than age 60 years having rates at a percentage similar to their age. Testing should only be performed on symptomatic patients because a large percentage of Helicobacter pylori-colonized individuals would have positive results. All positive tests should be interpreted with regard to patient clinical findings. At present, H. pylori serology is generally utilized to screen for H. pylori infection and breath test are used to confirm eradication after treatment unless endoscopy allows collection of tissues for rapid urease testing or histologic review.

**PRINCIPLE:**

The quantitative immunoenzymatic determination of specific antibodies is based on the Indirect ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound components HRP-labeled Anti-Human IgA is added to each well. Anti-Human IgA binds to the captured antibodies. In a second washing step, unbound conjugate is removed. The strips are then washed and, buffered substrate/chromogen reagent (hydrogen peroxide and Tetramethylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction, a blue color will develop and then enzymatic reaction will be stopped by adding stop solution. The intensity of the color is directly proportional to the amount of specific antibodies against H. pylori antigens from IgA class present in the sample.

<b>H.Pylori IgA ELISA KIT</b>		IVD-REF: PS - HPA
		Ver. No: 02

### Reagent and material supplied in the kit:

Reagent	Quantity/96 Test KIT	Reconstitution
Microplate (96 breakable wells) coated with H. pylori specific antigens	96 wells	Ready to Use
Standard set 1-5 (0, 10, 50, 100, 200 AU/mL) in buffer, containing preservative	5 Vial/2.0 ml	Ready to Use
Controls N = 2 (cut-off and high control- range indicated on the label) in buffer compatible with human serum matrix containing preservatives	2 Vial/2.0 ml	Ready to Use
Sample diluent (Red bottle)	1 Vial/ 30 ml	Ready to Use
Enzyme Conjugate- Anti-Human IgA conjugated to HRP (red color)	1 Vial/12 ml	Ready to use
Concentrated Wash solution	1 Vial/30 ml	Dilute 1:20 with distilled water
Substrate-Chromogenic solution (Tetramethylbenzidine and hydrogen peroxide)	1 Vial/12ml	Ready to Use
Stop solution (HCl 1.0 M). Additional information about the other additive is considered proprietary.	1 Vial/6.0 ml	Ready to Use

### Materials supplied

- 1 Adhesive cover
- 1 Instruction for use (IFU)

### Materials required but not supplied with kit

- The following materials are required but not provided in the kit:
- Distilled water with conductivity < 1µS
- Precise micropipettes for delivery of 20, 50, 100 and 200 microliters. An 8-channel pipette or resenser pipette with disposable plastic tips for delivery of 50-200 µL is useful but not essential. Pipettes for dispensing milliliter volumes.
- Vortex mixer
- Automatic washer or any other device which can deliver 350 microliter wash solution for microplates
- Microplate reader with the wavelength of 620-650 nm and 450 nm, and an absorbance range of 0-3

### Preparation of the Sample:

- Human serum and EDTA plasma are samples of choice for PSI H. pylori-IgA kit.
- Do not use grossly haemolysed or turbid samples.
- If the measurement cannot be done within 24 hours after sample collection, samples should be kept at 2-8°C up to 48 hours. For longer storage, serum or plasma must be stored at -20°C.
- Avoid repeated freeze and thaw cycles. Frozen specimens must be gently but thoroughly mixed after thawing

### PRECAUTION:

1. For Professional in vitro diagnostic use only.

<b>H.Pylori IgA ELISA KIT</b>		IVD-REF: PS - HPA
		Ver. No: 02

2. All the reagents within the kit are calibrated for the direct determination of human specific IgA against H. pylori antigens in human serum. The kit is not calibrated for the determination of specific IgA in Non-EDTA plasma or other specimens of human or animal origin.
3. Before starting the assay, read the working manual completely and carefully. Use the valid version of the instruction provided with the kit. Be sure that everything is understood.
4. Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., USA) publication No (CDC) 88—8\_95 on laboratory safety procedures or any other local or national regulation.
5. Handle all patient specimens as potentially infectious.
6. Reagents contain sodium azide (NaN<sub>3</sub>) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up. Follow local guidelines for disposal of all waste material.
7. Avoid any skin contact with all reagents. Stop solution contains strong acid in case of contact, wash thoroughly with soap and water.

### **Preparation of reagents:**

#### **Wash Buffer:**

According to the required volume, pour the desired volume of wash concentrate into a clean container and dilute 20 fold by adding adequate volume of distilled or deionized water to gain a buffered Wash Solution. E.g. 10 mL wash concentrate to 190 mL distilled water. Wash solution should be prepared daily. Please discard the unused working wash solution at the end of the day.

#### **PROCEDURE:**

1. Take out the microplate with the required number of strips. The left strips are placed in the plastic pouch along with the silica gel bag, and sealed.
2. Pipette 200 µl of Sample Diluent into each well except Standards and Controls well.
3. Pipette 200 µl of standards and controls (ready to use) into the appropriated wells.
4. Pipette 20 µl of Patients Serum into appropriated well.
5. It is strongly recommended that mix serum by pipetting up and down 5 times in sample diluent. As an alternative procedure, you can also dilute patient sample 1:11 out-of well in a small clean tube, mix thoroughly and gently, then transfer 200 µl, diluted serum to appropriate wells. Cover the strips with a plate sealer. Incubate at room temperature for 20 minutes.
6. Wash the plate according to procedure below.
  - a. Carefully dispense 350 µL of working wash solution into each well.
  - b. Aspirate the content of each well.
  - c. Repeat 4 additional wash for a total of 5 wash cycles.
7. Pipette 100 µl of Enzyme Conjugate into each well. Incubate at room temperature for 20 minutes.
8. Wash the wells 5 times as discussed in procedural note 6.
9. Pipette 100 µl of Substrate-chromogen into each well and Incubate in the dark at room temperature for 10 minutes.
10. Add 50 µL of Stop Solution, mix and read the absorbance at 450 nm (with a reference filter 620, 630, or 650) in a microplate spectrophotometer within 10 min after addition of Stop Solution.

<b>H.Pylori IgA ELISA KIT</b>		IVD-REF: PS - HPA
		Ver. No: 02

**Important Note:** always add stop solution in the same order as chromogenic solution TMB.

**Quantitative calculation of results:**


1. Construct a calibration curve by plotting the average O.D. (450 nm) for each calibrator on the vertical (Y) axis against the concentration on the horizontal (X) axis, and draw a best-fit curve through the points on the graph.
2. To determine concentration of H. pylori specific IgA in each sample, first locate the O.D. value on the Y-axis and extend a horizontal line to the calibration curve. At the point of intersection, draw a vertical line to the X-axis and read the corresponding concentration.

If a microplate spectrophotometer with built-in data calculation program is used, refer to the manual for the spectrophotometer and create a program using the concentration stated on the label of each calibrator. For automatic calculation of H.pylori specific IgA results, it is recommended to use point-to-point method.

**Typical Data:**

The following data are for illustration only and never should be used instead of the real-time calibration curve.

Row	OD	AU/mL
1	0.03	0
2	0.28	10
3	1.33	50
4	2.10	100
5	2.78	200

**Alternative qualitative evaluation**

The Cut-off value is given by the optical density (OD) of the cut-off control. The Cut-off index(COI) is calculated from the mean optical densities of the sample wells and Cut-off value. If the optical density of the sample is within a range of 20% around the Cut-off value (grey zone), the sample has to be considered as borderline. Samples with higher ODs are positive, samples with lower ODs are negative. For a semi quantification, the Cut-off index (COI) of the samples can be determined as follows:

$$COI = \frac{OD\ sample}{OD\ cut - off\ control}$$

**Antibody Isotype and state of infection**

serology	significance
<b>IgG</b>	Characteristic of secondary Immune response. High titers may indicate an active infection if clinically indicated, but titers may persist for several years. Low Positive results may indicate past infection.

<b>H.Pylori IgA ELISA KIT</b>		IVD-REF: PS - HPA
		Ver. No: 02

<b>IgA</b>	Produced in mucosal surfaces throughout the body and acts as a protective barrier. Usually produced early in infection course. High titers may indicate an active infection if clinically indicated.
------------	--

### Reference interval

In an in-house study apparently healthy subject showed the following results.

Qualitative interpretation	Quantitative results (AU/ml)	Semi quantitative results (COI)
Negative	<10	<0.8
Undetermined	10-15 AU/ml	08-1.2
Positive	>15 AU/ml	>1.2

The H.pylori serology results, should not be used as a sole basis for any therapeutic invention. They have to be interpreted in context of other clinical observations and diagnostic tests.

### Acceptable quality Criteria:

standard	OD
Standard Zero	< 0.1
Cut-off control	> 0.1
Standard 200	> 1.2

### Quality control

The test results are only valid if the test has been performed following the exact instruction. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable standards/laws. All validation criteria must be met as stated on the QC Certificate for standard 0 and standard 200. If the criteria are not met, the run is not valid and should be repeated. Each laboratory should use known samples as further controls. It is recommended to participate at appropriate quality assessment trials. In case of any deviation the following technical issues should be proven: Expiration dates of (prepared) reagents, storage conditions, pipettes, devices, incubation conditions and washing methods.

### PERFORMANCE CHARACTERISTICS

#### Diagnostic specificity

The diagnostic specificity is defined as the percentage (number fraction multiplied by 100) of subjects scoring negative in the absence of specific analyte. It is 97.0% (95% CI 89.25%-99.15%)

#### Diagnostic sensitivity

The diagnostic sensitivity is defined as the percentage (number fraction multiplied by 100) of subjects scoring positive in the presence of specific analyte. It is 98.0% (95% CI 87.11%-99.14%)

<b>H.Pylori IgA ELISA KIT</b>		IVD-REF: PS - HPA
		Ver. No: 02

**Precision:**

Precision was determined using PSI H.pylori IgA kit reagents and three pooled human sera according to EP05-A3 Guideline.

Three runs per day (morning, midday and afternoon) by three users in 10 working days (10 × 3 × 3). Results analyzed by fully nested ANOVA. Following results were obtained:

Sample Description	Mean (AU/mL)	Repeatability		Within-Laboratory Precision	
		SD(AU/mL)	%CV	SD(AU/mL)	%CV
Patient Pool	13.59	0.8	5.89	1.05	<b>7.73</b>
Patient Pool	18.42	1.0	5.43	1.4	<b>7.60</b>
Patient Pool	56.8	2.88	5.07	3.56	<b>6.27</b>

**1. Specificity:**

Interference evaluation from common sample abnormalities, hemolysis, Lipemia and icterus were performed according to guideline CLSI EP07-A2 and percentage interference was calculated for each sample using the equation below and normalized to the anti-H.pylori IgA content.

$$\%interference = \frac{measured\ value - true\ value}{true\ value} \times 100\%$$


No interference were seen up to following listed concentration:

Interferent	source	Final concentration of interferent in test serum
Hemoglobin	Fresh haemolysed blood	<b>50 mg/mL</b>
Conjugated Bilirubin	Commercially crystalline Unconjugated Bilirubin	<b>20 mg/dL</b>
Triglycerides	Commercially Intralipid 10%	<b>1000 mg/dL</b>

**References:**




1. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
2. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
3. Kao C.Y, Sheu B, Wu J-J. Helicobacter pylori infection An overview of bacterial virulence factors and pathogenesis . Biomed J(2016); 39(1) 14-23
4. Kasper D.L. et al. (2015). Harrison’s principles of internal medicine.19th ed. (pp. 1038-1041) Mc Graw Hill Education
5. Pagana K.D. et al. (2018). Mosby’s Manual of diagnostic and Laboratory tests.6th ed. (pp. 1049-1051). Elsevier Inc.
6. Stefan Reidel. et al. (2019). Jawetz, Melnick, & Adelberg’ s, Medical Microbiology.28th ed. (pp. 268-272). Mc Graw Hill (LANGE Medical Book)

### Brief Assay Procedure

<b>H.Pylori IgA ELISA KIT</b>		IVD-REF: PS - HPA
		Ver. No: 02

**SYMBOLS:**

	<b>BATCH NUMBER</b>
	<b>USE BY</b>
	<b>MANUFACTURER</b>
	<b>CONTAINS SUFFICIENT FOR &lt;n&gt; TESTS</b>
	<b><i>IN VITRO</i> DIAGNOSTIC MEDICAL DEVICE</b>
	<b>TEMPERATURE LIMITATION</b>
	<b>CATALOGUE NUMBER</b>
	<b>CONSULT INSTRUCTIONS FOR USE</b>