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$$\%Bias = \frac{\bar{X}\ measured - expected}{Expected} \times 100$$

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Sample	Original Concentration n (ng/mL)	2 ng/mL added		4 ng/mL added	
		Recovery%	Bias	Recovery %	Bias
1	0.33	97.5	-5.58%	99.2	2.77%
2	1.4	97.75	-1.33%	97.6	-1.76%
3	4.4	95	-1.56%	93.7	-2.98%
4	11.25	97.5	-0.38%	91.2	-0.69%

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
1. CLSI. Evaluation of detection capability for clinical laboratory measurement procedures; approved guideline.2nd ed. CLSI document EP 17-A2. Pierson-Perry J.F. Clinical and Laboratory Standard Institute, 2012.
2. CLSI. Evaluation of precision of quantitative measurement procedures; approved guideline. CLSI document EP 05-A3. McEnroe R.J. Clinical and Laboratory Standard Institute, 2014.
3. CLSI. Interference Testing in clinical chemistry; approved guideline.2nd ed. CLSI document EP07-A23. McEnroe R.J. Clinical and Laboratory Standard Institute, 2005.
4. Desirable Specifications for Total Error, Imprecision, and Bias, derived from intra- and inter-individual biologic variation, <https://www.westgard.com/biodatabase1.htm>
5. Kasper D.L. et al. (2018). Harrison's principles of internal medicine.20th ed.(pp. 3/17-4/17) Mc Graw Hill Education
6. Pagana K.D. et al. (2018). Mosby's Manual of diagnostic and Laboratory tests.6th ed. (pp. 378-82). Elsevier Inc.
7. Rifai Nader. et al. (2018). Clinical chemistry and molecular diagnostics.6th ed. (pp. 472-5). Elsevier Inc.

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PSA ELISA KIT

Direct immunoenzymatic determination of PSA in serum

Σ = 96 tests

REF PSPS-01

INTENDED USE

The Quantitative Determination of PSA Concentration in Human Serum.

CLINICAL SIGNIFICANCE

Prostate Specific antigen (PSA) is a serine protease with chymotrypsin-like activity. The protein is a single chain glycoprotein with a molecular weight of 28.4 kDA. PSA derives its name from the observation that it is a normal antigen of the prostate but is not found in any other normal or malignant tissue. PSA is found in benign, malignant and metastatic prostate cancer. Since prostate cancer is the second most prevalent form of male malignancy, the detection of elevated PSA levels plays an important role in the early diagnosis. Serum PSA levels have been found to be more useful than prostatic acid phosphatase (PAP) in the diagnosis and management of patients due to increased sensitivity. In this method, PSA calibrator, patient specimen or control is first added to a streptavidin coated well. Biotinylated monoclonal and enzyme labeled antibodies (directed against distinct and different epitopes of PSA) are added and the reactants mixed. Reaction between the various PSA antibodies and native PSA forms a sandwich complex that binds with the streptavidin coated to the well. After the completion of the required incubation period, the enzyme- PSA antibody bound conjugate is separated from the unbound enzyme-PSA conjugate by aspiration or decantation. The activity of the enzyme present on the surface of the well is quantitated by reaction with a suitable substrate to produce color. The employment of several serum references of known prostate specific antigen (PSA) levels permits the construction of a dose response curve of activity and concentration. From comparison to the dose response curve, an unknown specimen's activity can be correlated with PSA concentration.

PRINCIPLE

The technology uses streptavidin-Biotin with two high affinity monoclonal antibodies in an immunometric assay system. One of these antibodies is conjugated to biotin and the other is conjugated to horseradish peroxidase enzyme. The two antibodies react simultaneously with the antigen present in standards or samples and at the same time biotinylated antibody is captured to the solid phase by biotin-streptavidin interaction. This reaction leads to the formation of a streptavidin – biotinylated *capture antibody - antigen - signal antibody* complex, also referred to as a “sandwich”. For replacement of unreacted reagents, microtiter wells are washed repeatedly. After the addition of a ready-to-use tetramethylbenzidine substrate, TMB is catalyzed by the HRP to produce a blue color product that changes to yellow after adding stop solution. The signal is measured using an ELISA reader at 450nm and 630nm reference filter.

Reagent and material supplied in the kit

1.Microwell 1 plate: 12 x 8 wells coated with streptavidin.

2.Standards Set : 6 vials, 0.5 ml per vial, containing 0, 1.0 ,5.0 ,10 ,20 and 40 ng/ml PSA, calibrated against NIBSC/WHO, 668/96 serum with preservative.

3.Enzyme Conjugate: 1 vial, 12 ml, monoclonal anti-PSA Conjugated to biotin plus monoclonal anti-PSA Conjugated to HRP. **ready to use**

4.Substrate - Chromogen: 1 vial, 11 ml. Hydrogen peroxide and (TMB) solution, **ready to use**

5.Stop Solution: 1 vial, 6 ml, 1N HCL, **ready to use**

6.Wash Buffer: 1 vial, 30 ml, concentrated buffer with Tween-20, **(20x)**.

7.Control Serum: 1 vial, 0.5 ml, containing PSA, in serum with preservative. the exact concentration is printed on the label.

Preparation of the Sample

Serum samples can be prepared according to common procedures used routinely in clinical laboratory practice. Samples can be stored at 2-8 °C if the assay is carried out within 24 hours, otherwise aliquots should be prepared and stored deep frozen (-20 °C). Frozen samples should be thawed and thoroughly mixed before assaying. Repeated freezing and thawing should be avoided.

PRECAUTION

1. In one test-run do not combine strips, conjugate and standards from kits which have different lot numbers.
2. Solutions containing TMB and/or peroxide should not come into contact with metals or metal- ions, since this may give rise to unwanted color formation.
3. Allow that samples and all reagents come to ambient temperature before use.

Preparation of Wash Solution


1. Dilute the concentrated wash buffer 1:20 with distilled water. The diluted wash buffer must be at 20 to 25°C when used.
2. Incomplete washing will adversely affect the test results. Operating instructions for the washing equipment should be followed. Wash five times each plate. Blot dry the microplate by inverting the plate onto absorbent tissue, and striking a hard surface several times. If no automatic washer is available, washing can be performed manually as follows: Aspirate thoroughly the liquid in the wells. Then fill each well (0.3ml) with diluted wash solution. Repeat these steps 5 times. Complete the wash cycle by blotting the plate dry with absorbent tissue.

PROCEDURE

1. take out the microplate with the required number of strips. The left strips are placed in the plastic pouch along with the silica gel bag, and sealed
2. Pipette **20 µl of standards, control and sera** into the wells.
3. Pipette **100 µl of enzyme conjugate** into each well.
4. Cover the strips with a plate sealer. Incubate **at room temperature for 45 minutes**.
5. Wash each well with the solution **above five times** and then blot dry by pressing plate onto absorbent tissue. (See Directions for washing).
6. Pipette **100 µl of Substrate-chromogen** into each well and Incubate **in dark and room temperature for 15 minutes**
7. Stop the reaction by adding **50 µl of stop solution** to each well and mixing completely.
8. Put the plate in the microplate reader and read (within 10 minutes after step 7) the absorbance of the solution in the wells at **450nm and 630 nm**.

Calculation of Results

1. Calculate the mean absorbance value for each set of standards and samples at 450 nm.
2. Construct a standard curve by plotting the mean absorbance obtained from each reference standard against its concentration in ng/ml on graph paper, with absorbance values on the vertical or Y axis, and concentrations on the horizontal or X axis.
3. Use the mean absorbance values for each sample to determine the corresponding concentration of PSA in ng/ml from the standard curve.

Standards ng/ml	OD 450nm
0	0.022
1.0	0.342
5.0	1.19
10	1.94
20	2.50
40	2.90

REFERENCE VALUE

Healthy males are expected to have values below 4 ng/ml.

PERFORMANCE AND CHARACTERISTICS

Sensitivity

For the analytical sensitivity **0.15 ng/ml** has been obtained by assaying 15 replicates of the zero standard. It is defined as the concentration corresponding to the sum of the mean OD and its double standard deviation.

. Specificity

No interference was detected with the performance of this kit upon other proteins.

Precision and reproducibility

Intra Assay Variation

Within run variation was determined by replicate determination of 3 different sera in one assay.

CV %	SD	Mean (ng/ml)	Number of determinations	Sample
8.46	0.11	1.25	15	1
7.61	0.28	3.65	15	2
6.22	0.58	9.33	15	3

Inter Assay Variation

Between run variations was determined by replicate measurements of three different sera.

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Linearity

linearity was investigated by serial dilution of three samples of known PSA concentrations. All samples were diluted with zero standard of FSH ELISA Kit.

Recovery %			First Conc. (ng/ml)	Samples
1:8	1:4	1:2		
105	100	97	4.2	1
104	106	99	10.5	2
104	101	103	34.8	3

Hook Effect

In this kit up to 200 ng/ml there isn't hook effect.

References

1. Wild D, *The Immunoassay Handbook.*, Stockton Press (1994) p452.
2. Junker R, Brandt B, Zechel C, Assmann G, *Clin Chem*, **43**, 1588- 94 (1997).
3. Prestigiacomo AF, Stamey TA, ' *Physiological variations of serum prostate antigen in the (4-10 ng/ml) range in male volunteers.*' J. Urol 1996;**155**:1977-80.
4. Stamey TA, McNeal JE, Yemoto CM, Sigal BM, Johnstone IM, ' *Biological determinants of cancer progression in men with prostate cancer.*' JAMA 1999;**281**:1395-1400.