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$$\frac{\text{abs Sample } (\Delta\text{abs/min})}{\text{abs Calibrator } (\Delta\text{abs/min})} \times \text{Calibrator activity (IU/L)} = \text{Sample activity (IU/L)}$$

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Measuring Range: 10-450 IU/L
Limit Of Blank (LOB): 0.0 IU/L
Limit Of Detection (LOD): 2.5 IU/L
Limit Of Quantification (LOQ): 10 IU/L

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Level	n	Mean (IU/L)	Within Run CV%	Between Run CV%
Low	80	34	1.7	3.4
Medium	80	68	0.8	1.9
High	80	535	0.4	2.0

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Correlation Coefficient: (r)= 0.999
 Linear regression: Y= 0.927 (x) – 0.3 IU/L

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- Burtis CA, Ashwood ER, Bruns DE. Tietz Fundamentals of Clinical Chemistry 6th Edition, Saunders Elsevier, St. Louis, MO. 2008.317
- Tietz, N.W., Clinical guide to laboratory tests, 4th Ed. (W.B. Saunders eds. Philadelphia USA), (2006), 154.
- Itani, O., Tsang, R.C., The Liver: Function and clinical pathology, Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 586 and appendix.
- Anderson SC, Cockayne S. Clinical chemistry: concepts and applications. McGraw-Hill/Appleton & Lange; 1993, 238.
- Schumann G, et al. International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. International Federation of Clinical Chemistry and Laboratory Medicine: Part 4. Reference procedure for the measurement of catalytic concentration of alanine aminotransferase. Clin Chem Lab Med. 2002;40(7):725-33.
- World Health Organization. Use of Anticoagulants in Diagnostic Laboratory Investigations & Stability of blood, plasma and serum samples (WHO/DIL/LAB/99.1 Rev. 2), 2002.

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	Temperature limitation		Catalogue number
	Manufacture address		Expiration date
	Batch code		Date of manufacture
	In vitro diagnostic medical device		Reagent 1
	Consult instruction for use		Reagent 2