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- Rifai, N., et al. Lipids, Lipoproteins, and Apolipoproteins. Titez fundamentals of clinical Chemistry, 5th Ed., Burtis, C.A & Ashwood, E.R (W.B. Saunders eds. Philadelphia USA), (2001), 463.
- Naito, H.K., Coronary Artery Disease and Disorder of Lipid Metabolism. Clinical Chemistry: Theory, Analysis, Correlation, 4th Ed. Kaplan, L.A., Pesce, A.J, Kazmierczak, S. C. (Mosby, Inc.eds. St Louis USA), (2003), 603.
- Allain, C. C., et al, Clin Chem., (1974), 20, 470.
- Tietz, N.W., Clinical guide to laboratory tests, 4th Ed. (W.B. Saunders eds. Philadelphia USA), (2006), 244
- Expert Panel on Detection, JAMA, (2001), 285, 2486.
- Vassault A., et al., Ann. Biol. Clin., (1986), 44, 686.
- Vassault A., et al., Ann. Biol. Clin., (1999), 57, 685.
- Young, D. S., Effects of preanalytical variables on clinical laboratory tests, 2nd. AACC Press, (1997).
- Young, D. S., Effects of drugs on clinical laboratory tests, 4th Ed. AACC Press, (1995).
- <https://www.mayocliniclabs.com>
- Berth, M. & Delanghe, J., Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Act Clin Belg., (2004), 59, 263.
- Evans, K., Mitcheson, J., & Laker, M. F. (1997). Effect of storage at—70° C on lipid, lipoprotein and apolipoprotein concentrations, Clinica Chimica Acta, 258(2), 219-229.
- Nauck, M., Warnick, G. R., & Rifai, N. (2002). Methods for measurement of LDL-cholesterol: a critical assessment of direct measurement by homogeneous assays versus calculation. Clinical chemistry, 48(2), 236-254.

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Temperature limitation

Catalogue number

Manufacture address

Expiration date

Batch code

Date of manufacture

In vitro diagnostic medical device

Reagent 1

Consult instruction for use

Reagent 2