



مشخصات کالا

نام کالا : تجهیزات آزمایشگاهی/Immunology/Immunology Reagents/ معرف های تشخیص نشانگرهای
تومور/Other Tumour Markers/Alphafetoprotein

نام انگلیسی وسیله پزشکی: Alphafetoprotein

Immunoassay reagents intended to perform qualitative and/or quantitative analyses of a sample of body fluids (e.g., maternal and fetal serum, amniotic fluid) to determine alpha-fetoprotein (AFP), a fetal glycoprotein secreted during pregnancy, initially by the fetal yolk sac and later by the fetal liver. Abnormally high levels of AFP in maternal serum are associated with many conditions, including fetal open neural tube defects such as anencephaly and spina bifida (meningomyelocele); low levels are found in conditions such as trisomy-21 (Down's syndrome) and trisomy-18 pregnancies. AFP is secreted after birth by neoplastic tumor cells (i.e., a tumor marker). High levels (typically greater than 1,000 microgram/liter) are indicative of hepatocellular carcinoma; the level of AFP is used for prognosis and monitoring of hepatocellular carcinoma.

شرح وسیله:

کد (های) UMDNS: 17259

کمپانی سازنده قانونی
(Legal License Holder)
پیش‌تاز طب زمان :

کشور کمپانی سازنده قانونی
(Legal License Holder): Iran

مدل : 96 تستی

نوع شرکت به عنوان : تولید کننده

وضعیت نمایندگی : نماینده قانونی

واحد شمارش کالا : False

حیطه کاربرد : تعیین میزان آلفا فیتو پروتئین در سرم یا پلاسما به روش الیزا

نام های رایج :

نام کالا در برچسب : AFP

برند : پیش‌تاز طب زمان

شماره کاتالوگ : PT-AFP-96

مشخصه
های
فنی

RIA ☐	ECLIA ☐	CLIA ☐	EIA ☒
Nephelometric immunoassay ☐	IFA ☐	Immunochromatographic assay ☐	LIA ☐
Immunoblot ☐	Immunochemistry ☐	Immunoprecipitation ☐	Turbidimetric immunoassay ☐
Immunodiffusion ☐	Microlymphocytotoxicity assay ☐	Complement fixation ☐	Immunofluorometric assay (IFMA) ☐
		Others ☐	Molecular-based assay ☐
Urine ☐	Whole Blood ☐	Plasma ☒	Serum ☒
	Other Body Fluids ☐	CSF ☐	Stool ☐

روش
سنجش
:

نوع
نمونه :

تعداد
تست : 96

نوع
سیستم
بر اساس
سازگاری
کیت : Open

اجزاء
کیت : پلیت-کالیبراتور-محلول شستشو و متوقف کننده و رنگزا

شماره
کاتالوگ
: PT-AFP-96

نوع کیت
: Quantitative

کلیه حقوق معنوی این نرم افزار متعلق به اداره کل تجهیزات پزشکی وزارت بهداشت، درمان و آموزش پزشکی می باشد.

طراحی و تولید نرم افزار : نوآوران فناوری اطلاعات امروز (<http://noict.com>)